TPHA-DAT ADVANCED

DIESSE

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For *In Vitro* Diagnostic Use Only

1. INTENDED USE

TPHA-DAT ADVANCED (REF 26041) is a hemagglutination test for the semi-quantitative detection of specific antibodies directed against *Treponema pallidum*.

The test is performed on human serum by passive hemagglutination using the AUTO-DAT instrument (REF 26001). The kit is intended to determine exposure to *Treponema pallidum* infection and to be used as an aid in the diagnosis of Syphilis.

The test is not intended to be used to detect exposure to *Treponema pallidum* infection in blood, blood components, cells, tissues, organs, or any derivative thereof for the purpose of assessing their suitability for transfusion, transplantation, or cell administration.

It should be used only by professional laboratory personnel.

2. INTRODUCTION

Syphilis, caused by *Treponema pallidum* (subsp. *pallidum*) spirochetes, is a chronic infection with many different clinical manifestations occurring in distinct stages. Generally, this systemic infectious disease is contracted through direct sexual contact and has lesions containing treponemes. The only known hosts are humans.

Many tests are available for direct and indirect diagnosis of syphilis, but there is still no single optimal test.

Direct diagnostic methods include detection of *T. pallidum* by microscopic examination of fluid or smears from lesions, histological examination of tissues, or nucleic acid amplification methods such as polymerase chain reaction (PCR).

Indirect diagnosis is based on serologic tests for the detection of antibodies. Serological tests fall into two categories:

1. nontreponemal screening tests
2. confirmatory treponemal tests.

Treponemal tests are specific tests for syphilis in that they detect specific antibodies to *Treponema pallidum* antigens. Generally, these tests are not used for assessment of treatment response because of persistent reactivity throughout the patient's life.

Diagnosis should be based on correlation of test results with other clinical findings.

3. PRINCIPLE OF THE METHOD

The test is based on the principle of hemagglutination.

Specific antibodies in the patient's serum bind to treponemal antigens (Nichols strain) attached to avian erythrocytes.

The reaction is positive in case of agglutination of cells. The reaction is negative in case the aggregate has the shape of a compact button.

The test is applied to the AUTO-DAT instrument (REF 26001), which captures the image of the reaction that has occurred, analyzes it, provides and prints a result calculated on the basis of a specific algorithm.

The results are expressed in titer, and confirmation of the visually obtained titer is possible.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry.

All materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use, bring the reagents to room temperature (18-30°C).

1. The test should be used together with the AUTO-DAT instrument (REF 26001), strictly following the Instructions for Use and User's Manual of the instrument.
2. Do not change the Test Procedure. Do not substitute reagents with those from other suppliers or other batches unless it is specifically stated that the reagent is interchangeable between batches.
3. Check that the instrument is set up correctly (see User's Manual).
4. Avoid the use of self-defrosting freezers for sample storage.
5. Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the directions in the chapter "Analytical Specificity").
6. Do not use the device after the expiry date.
7. After use, store reagents at 2-8°C.
8. Avoid microbial pollution of reagents as this reduces the validity of the product and may result in erroneous results.

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 50 tests.

TEST CELLS

3 x 4 mL

Content: avian erythrocytes sensitized with *Treponema pallidum* antigen.

CONTROL CELLS

1 x 4 mL

Content: non-sensitized avian erythrocytes.

Warning: Store vials upright when not in use. Bring to room temperature before use. Do not shake vigorously.

SAMPLE DILUENT

1 x 13 mL

Contents: saline solution containing substances inhibiting nonspecific reactions.

CONTROL +

1 x 0.3 mL

Contents: diluted human serum containing anti-*Treponema pallidum* antibodies. Liquid, ready-to-use.

CONTROL -

1 x 0.3 mL

Contents: diluted animal (rabbit) serum not containing anti-*Treponema pallidum* antibodies. Liquid, ready-to-use.

MATERIALS REQUIRED BUT NOT PROVIDED

- AUTO-DAT INSTRUMENT (REF 26001)
- AUTO-DAT MICROPLATES TPHA (REF 26004)
- Micropipettes capable of accurately taking volumes of 5-100 µl.
- Disposable gloves
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

After first use, microplates should be stored at room temperature.

The expiry date is printed on each component and on the kit label.

The stability of the reagents does not change after opening the bottle, as long as the user takes care to keep the product away from possible microbial contamination.

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be stored for 6 days at 2/8°C, or frozen for longer periods at temperature ≤ -20°C for a least 11 months

The sample can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

Evaluation is done by observing sample agglutination at dilutions 1/80, 1/640 and 1/5120. To perform pre-dilution of samples, use the wells in rows A and E, following the steps below.

1. Add 125 µl of Sample Diluent to well A1 and 35 µl of Sample Diluent to wells C1 and D1.
2. Add 5 µl of sample to well A1. Be careful not to transfer any cellular elements. Homogenize the solution while avoiding the formation of bubbles.
3. Take 35 µl of solution in A1 and transfer to well B1. Avoid the formation of bubbles.
4. Take 5 µl of solution in A1 and transfer to well C1. Homogenize thoroughly.
5. Take 5 µl of solution in C1 and transfer to well D1. Homogenize thoroughly and discard 5 µl of solution.
6. Discard 55 µl from well A1.
7. Repeat the same steps for the second sample using wells E1 to H1 and columns 2 to 12 for the other samples.

NOTE: Positive and negative controls are supplied ready-to-use and do not require pre-dilution. They should be used at each run as described below.

- a) Add 35 µl of Sample Diluent in C and D.
- b) Add 35 µl of Control (positive or negative) in A and B.
- c) Add 5 µl of Control (positive or negative) in C. **Homogenize without forming bubbles.**
- d) Take 5 µl from C and transfer to well D. **Homogenize thoroughly and discard 5 µl of solution.**
- e) Perform the same steps for the other Control using wells E to H.

Once the dilutions for samples and controls have been made, proceed to add the reagents:

8. Add 75 µl of CONTROL CELLS in A and E.
9. Add 75 µl of TEST CELLS to the other wells.

Sample dispensing scheme

		Step 1	Step 2	Step 3
		1	1	1
SAMPLE	A	125µl DILUENT +5µl SAMPLE	DISCARD 55µl	+75µl CONTROL CELLS
	B	35µl A1		+75µl TEST CELLS
	C	35µl DILUENT +5µl A1		+75µl TEST CELLS
	D	35µl DILUENT +5µl C1	DISCARD 5µl	+75µl TEST CELLS

Control dispensing scheme

		Step 1	Step 2	Step 3
		1	1	1
CONTROL +	A	+35µl CONTROL +		+75µl CONTROL CELLS
	B	+35µl CONTROL +		+75µl TEST CELLS
	C	35µl DILUENT +5µl CONTROL +		+75µl TEST CELLS
	D	35µl DILUENT +5µl C1	DISCARD 5µl	+75µl TEST CELLS
CONTROL -	E	+35µl CONTROL -		+75µl CONTROL CELLS
	F	+35µl CONTROL -		+75µl TEST CELLS
	G	35µl DILUENT +5µl CONTROL -		+75µl TEST CELLS
	H	35µl DILUENT +5µl G1	DISCARD 5µl	+75µl TEST CELLS

Once dilutions have been made and reagents added place the plate in the AUTO-DAT instrument (REF 26001) and follow the instrument's instructions.

In case the instrument is occupied, the test incubation can be carried out at room temperature on a vibration-free bench. In such a case, after 60 minutes of incubation, insert the microplate into the AUTO-DAT instrument with utmost caution, avoiding accidental shocks.

At the end of the analysis, remove the plates from the instrument and cover the wells where the reaction took place with the plate cover.

WARNING: If antigen is not added to the sample, the instrument will register an error.

9. TEST VALIDATION

Use the controls provided at each session.

If the result is different than expected, contact Customer Care.

Tel: 0039 0577 319554
email: scientificsupport@diesse.it;
customercare@diesse.it

10. INTERPRETATION OF THE RESULTS

The instrument gives the result in titer.

The test on the test sample can be interpreted as follows:


POSITIVE: when the result is $\geq 1/80$

NEGATIVE: when the result is $< 1/80$

It is recommended to visually check the results obtained to confirm the report (especially in the case of NEG result $< 1/80$).

In the absence of agglutination the erythrocytes sediment, forming a compact button in the center of the well.

In the case of agglutination, the erythrocytes sediment, forming a homogeneously diffuse layer at the bottom of the well and/or forming a characteristic ring structure.

Negative NEG $< 1/80$

Positive: $\geq 1/80$

Agglutination of Test Cells and simultaneous absence of agglutination of Control Cells indicates the presence of anti-Treponema pallidum antibodies.

Control Cells sediment forming a particularly compact button; in case their agglutination occurs, the test is to be considered invalid.

The presence of a compact button of erythrocytes sedimented at the bottom is interpreted as a negative result.

11. LIMITATIONS

The product should be used only by professional laboratory personnel.

The test is not suitable for samples different from human serum. The test is only applicable to the AUTO-DAT instrument (REF 26001).

Treponemal tests become positive 5-15 days after the onset of symptoms. No single test or reference standard is available for each stage of the disease.

No hemagglutination test on serum can distinguish antibodies arising from infections caused by Treponema pallidum from those arising from infections caused by other treponemal pathogens.

All values obtained need careful interpretation that does not prescind from other indicators related to the same patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. REFERENCE RANGE

Measurement Range: NEG - 1/20480 (Titer)

13. ANALYTICAL SPECIFICITY

3 samples (1 Negative, 1 Cut-Off and 1 Positive) were spiked with the following potentially interfering factors and then tested :

Rheumatoid Factor (44 – 220 IU/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (250 mg/dl – 1500 mg/dl)
Hemoglobin (2.5 mg/ml – 10 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

14. CROSS-REACTIONS

18 samples, positive to Borrelia (3), Herpes Virus (3), Toxoplasma (3), Epstein-Barr (3), Cytomegalovirus (3), Rubella (3) were tested.

No significant cross-reactions were found except for Epstein-Barr virus.

15. METHOD COMPARISON

In an experimentation 220 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Total
Diesse	+	127	0	127
	-	0	93	93
	Total	127	93	220

Overall Percent Agreement = 100.0% CI95%: 93.3-100.0
 Percent Positive Agreement (~Sensibilità Diagnostica)= 100.0%
 CI95%: 97.1-100.0
 Percent Negative Agreement (~Specificità Diagnostica)
 =100.0% CI95%:96.0-99.9
 Positive Predictive Value (PPV) = 100.0% CI95%: 100.0-100.0
 Negative Predictive Value (NPV) = 100.0% CI95%: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.00.

16. PRECISION AND REPEATABILITY

Precision and repeatability were evaluated on replicates of negative and positive samples, in three different batches, with different operators performing the test, on different days, and with different instruments.

All results obtained fell within the acceptance criteria (Agreement % $\geq$ 95% between expected and obtained results) demonstrating the accuracy and repeatability of the method.

17. REFERENCES

- I. Ratnam S. The laboratory diagnosis of syphilis. Can J Infect Dis Med Microbiol. 2005 Jan;16(1):45-51. doi: 10.1155/2005/597580. PMID: 18159528; PMCID: PMC2095002.
- II. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
- III. WHO guidelines for the treatment of *Treponema pallidum* (Syphilis) (201).
- IV. Janier M, Unemo M, Dupin N, Tiplica GS, Potočník M, Patel R. 2020 European guideline on the management of syphilis. J Eur Acad Dermatol Venereol. 2021 Mar;35(3):574-588. doi: 10.1111/jdv.16946. Epub 2020 Oct 22. PMID: 33094521.
- V. Morshed MG, Singh AE. Recent trends in the serologic diagnosis of syphilis. Clin Vaccine Immunol. 2015 Feb;22(2):137-47. doi: 10.1128/CVI.00681-14. Epub 2014 Nov 26. PMID: 25428245; PMCID: PMC4308867.
- VI. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

18. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

19. SUMMARY OF SAFETY AND PERFORMANCE

This document, which will be made available on the EUDAMED database (when this is fully implemented and functioning), is part of the Technical Documentation and can be requested from the manufacturer.

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