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Overall Percent Agreement = 97.3% % CI_{95%}: 93.8-98.8
 Percent Positive Agreement (~Sensibilità Diagnostica) = 96.1%
 CI_{95%}: 89.1 -98.6
 Percent Negative Agreement (~Specificità Diagnostica) =98.1%
 CI_{95%}: 93.5 -99.5
 Positive Predictive Value (PPV) = 97.4% CI_{95%}: 95.1 -99.7
 Negative Predictive Value (NPV) = 97.2% CI_{95%}: 94.8 -99.6

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- Jukic A, Bakiri L, Wagner EF, Tilg H, Adolph TE. Calprotectin: from biomarker to biological function. Gut. 2021;70(10):1978-1988.
- Franzosa EA, Sirota-Madi A, Avila-Pacheco J, et al. Gut microbiome structure and metabolic activity in inflammatory bowel disease [published correction appears in Nat Microbiol. 2019 May;4(5):898]. Nat Microbiol. 2019;4(2):293-305.
- Carroccio A, Iacono G, Cottone M, et al. Diagnostic accuracy of fecal calprotectin assay in distinguishing organic causes of chronic diarrhea from irritable bowel syndrome: a prospective study in adults and children. Clin Chem. 2003;49(6 Pt 1):861-867.
- Alibrahim B, Aljasser MI, Salh B. Fecal calprotectin use in inflammatory bowel disease and beyond: A mini-review. Can J Gastroenterol Hepatol. 2015;29(3):157-163.
- Lundgren D, Eklöf V, Palmqvist R, Hultdin J, Karling P. Proton pump inhibitor use is associated with elevated faecal calprotectin levels. A cross-sectional study on subjects referred for colonoscopy. Scand J Gastroenterol. 2019;54(2):152-157.

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INSTRUCTIONS FOR USE

CHORUS CALPROTECTIN 3RD GENERATION

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

CHORUS CALPROTECTIN 3RD GENERATION is an immunoassay kit for the automated quantitative determination of Calprotectin in stools.

The test is performed on feces using a disposable device attached to the CHORUS and CHORUS TRIO instruments.

High concentration of Calprotectin is observed in the feces of patients with chronic inflammatory diseases (IBD). The kit, intended for the determination of Calprotectin concentration is used as an aid to diagnosis and management of related therapeutic measures.

It is to be used only by professional laboratory personnel.

2. INTRODUCTION

Calprotectin is a 36-kDa heterodimeric protein found predominantly in neutrophils. Its release can occur in the presence of numerous cytotoxic effects that make the cell membrane permeable. Its concentration thus tends to increase in plasma, CSF, synovial fluid, urine, and feces. Its elevated concentration in feces may indicate intestinal inflammation, and its correct determination may aid in a diagnosis of chronic inflammatory disease and guide and control drug choices in support of clinical evidence and other relevant diagnostic investigations. In addition, increased values can be observed in all diseases that involve the activation of an acute or chronic inflammatory process of the gastrointestinal system. Administration of anti-inflammatory drugs, celiac disease, and bacterial infections can change Calprotectin levels.

3. PRINCIPLE OF THE METHOD

The CHORUS Calprotectin 3rd generation device is ready-to-use for the determination of Calprotectin in feces in CHORUS instruments.

The test is based on the ELISA (Enzyme Linked ImmunoSorbent Assay) principle. Calprotectin-specific antibodies are bound to the solid phase. Calprotectin binds to the solid phase following incubation with the stool extract. After washes to remove unreacted proteins, incubation with the conjugate consisting of monoclonal anti-calprotectin antibodies conjugated with horseradish peroxidase is performed. The conjugate that has not bound is removed and the substrate for the peroxidase is added. The blue color that develops is proportional to the concentration of Calprotectin presents in the test sample.

The disposable devices contain all the reagents to perform the test when applied to the Chorus/Chorus TRIO instruments.

Results are expressed in mg/Kg.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use, bring the package containing the devices to room temperature (18-30°C) for at least 30 minutes for at least 30 minutes.

Use the device within 60 minutes.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed, and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a

Release (Rel.) subsequent to the one reported in the table published on Diesse website (<https://www.diesse.it/it/downloads/downloads/strumeto:39/>)


5. Check that the instrument is set up correctly (see CHORUS Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606)**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space available for barcode label

Position 7: Well with STOP SOLUTION

Position 6: MICROPLATE WELL.

Sensitized with monoclonal anti-calprotectin antibodies (maximum concentration 20µg/Well)

Position 5: MICROPLATE WELL

Not sensitized.

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer.

Position 3: SAMPLE DILUENT

Content: Protein buffer solution containing Proclin (0.1%) and Tween (0.1%).

Position 2: CONJUGATE

Contents: anti-Calprotectin monoclonal antibodies (maximum concentration 0.35µg/ml) labeled with peroxidase, in protein buffer solution containing Proclin (0.1%) and Tween (0.05%).

Position 1: EMPTY WELL.

Where the sample is transferred to.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.425 ml

Contents: Calprotectin solution (maximum concentration 20µg/ml) containing Proclin (0.1%) and Tween (0.1%). Liquid, ready-to-use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Calprotectin solution (maximum concentration 20µg/ml) containing Proclin (0.1%) and Tween (0.1%). Liquid, ready-to-use.

SAMPLE DILUENT/NEGATIVE CONTROL 1 x 1.8 ml

Content: Contents: Protein buffer solution containing Proclin (0.1%) and Tween (0.1%). Liquid, ready-to-use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

The Calibrator and Positive Control are produced from a solution containing Calprotectin diluted to achieve a specific concentration, the range of which is batch-dependent and assigned during the quality control release phase using a series of secondary calibrators and are prepared and characterized in accordance with a panel of human reference fecal samples with different antigen levels.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS/CHORUS TRIO/CHORUS EVO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
SAMPLE DILUENT/NEGATIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample type is stool extract, prepared as reported in the Instructions for Use of the product Chorus Coprocollect (REF 86602 - 86602/40) and handled as required in standard laboratory procedures.

The sample can be kept for 4 days at 2/8°C; for longer storage periods, freeze at -20°C for up to 3 months. The sample may undergo up to a maximum of 3 thawings. Avoid the use of self-defrosting freezers for sample storage. After thawing, shake the sample carefully before assaying.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	50 µl/device
POSITIVE CONTROL	50 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Customer Care.

Tel: 0039 0577 319554

email: scientificsupport@diesse.it;
customercare@diesse.it

10. INTERPRETATION OF THE RESULTS

The instrument provides the result in mg/Kg, calculated on the basis of a batch-dependent curve stored in the instrument. The test on the sample under examination can be interpreted as follows:

POSITIVE: when the result is > 119 mg/Kg

NEGATIVE: when the result is < 51 mg/Kg

DOUBT/EQUIVOCAL: when the result is between 51 and 119 mg/Kg

In case of doubtful/equivocal result, repeat the test.

Levels found in infants' first months of life may be significantly higher.

11. LIMITATIONS

The product should be used only by professional laboratory personnel.

The test is not suitable for samples different from faecal sample. All values obtained need careful interpretation that does not prescind from other indicators related to the same patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

The test is applicable only to the CHORUS and CHORUS TRIO instruments; it is not applicable to the CHORUS EVO instrument. Patients who are taking nonsteroidal anti-inflammatory drugs (NSAIDs) or proton pump inhibitors (PPIs) may have increased fecal calprotectin levels.

Patients with the following diseases may have elevated calprotectin levels: microscopic colitis, diverticular disease, celiac disease, gastrointestinal infections, and colorectal cancer.

12. CALIBRATION RANGE

Calibration range 24.-1000 mg/Kg.

For samples > 1000 mg/Kg Repeat the test on the sample diluted 1:5 in sample dil/neg control.

The result obtained from the instrument should be multiplied by 5.

13. REFERENCE RANGE

The expected values in the normal population, determined by examining 120 healthy donor samples, were <24mg/Kg.

14. ANALYTICAL SPECIFICITY

3 samples (1 Negative and 2 Positive) to which the following interferers were added:

-	Vancomycin	(0.67 mg/50 mg stool)
-	Ciprofloxacin	(0.50 mg/50 mg stool)
-	Vitamin E	(0.10 mg/50 mg stool)
-	Lansoprazolo	(0.02 mg/50 mg stool)
-	Pentasa	(1.33 mg/50 mg stool)
-	Asacol	(1.33 mg/50 mg stool)
-	Prednisone	(0.01 mg/50 mg stool)
-	Vitamin A	(8.0 UI/50 mg stool)
-	Vitamin C	(0.05 mg/50 mg stool)
-	Vitamin D	(1.1 UI/50 mg stool)
-	Hemoglobin	(5.83 mg/50 mg stool)

The presence of the above interfering substances in the faecal sample under examination does not alter the test result.

15. CROSS-REACTIONS

Three samples (1 Negative and 2 Positives) were tested to which *Escherichia coli*, *Salmonella enteritidis*, *Shigella flexneri*, *Yersinia enterocolitica*, *Klebsiella pneumoniae*, *Citrobacter freundii*, *Campylobacter jejuni* were added.

No significant cross-reactions were detected.

16. METHOD COMPARISON

In an experimentation 185 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Total
Diesse	+	74	2	76
	-	3	106	109
	Total	77	108	185

Overall Percent Agreement = 97.3% CI_{95%}: 93.8-98.8

Percent Positive Agreement = 96.1% CI_{95%}: 89.1 -98.6

Percent Negative Agreement =98.1% CI_{95%}: 93.5 -99.5

Positive Predictive Value (PPV) = 97.4% CI_{95%}: 95.1 -99.7

Negative Predictive Value (NPV) = 97.2% CI_{95%}: 94.8 -99.6

The agreement between the two methods is excellent with a Cohen's Kappa of 0.94.

Of the 185 samples, correlation between the Diesse kit and the trade kit was tested out of 100 quantifiable samples.

The data are summarized in the following table:

Correlazione	r	IC 95%
Pearson	0.98	0.97-0.99
Spearman	0.97	0.96-0.98

The correlation between the two methods is very high. This correlation is confirmed by the results obtained from the Passing-Bablok and Bland-Altman tests.

17. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean mg/Kg	CV%	Mean mg/Kg	CV%
1	24.0*	-	24.0*	-
2	44.1	3.3	41.6	13.6
3	189.9	7.3	199.1	4.4
4	919.7	9.6	878.2	4.1

Sample	Precision between batches		Precision between instruments	
	Mean mg/Kg	CV%	Mean mg/Kg	CV%
1	24.0*	-	24.0*	-
2	43.4	4.8	67.8	13.0
3	187.9	12.5	142.2	13.4
4	890.4	9.4	521.9	14.4

*The obtained experimental values <24.0 mg/Kg, were treated as 24.0 mg/Kg to allow the calculation of the results.

18. REFERENCES

- Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
- Fagerhol MK. Nomenclature for proteins: is calprotectin proper name for the elusive myelomonocytic protein? Clin Mol Pathol 1996;49:M74-9.
- Barassi A, Melzi d'Eril GV. La calprotectin nella malattia infiammatoria e nel cancro del colon. RIME/LIJLaM 2008;4:104-8
- Hoskin TS, Crowther JM, Cheung J, et al. Oxidative cross-linking of calprotectin occurs in vivo, altering its structure and susceptibility to proteolysis. Redox Biol. 2019;24:101202.
- Jukic A, Bakiri L, Wagner EF, Tilg H, Adolph TE. Calprotectin: from biomarker to biological function. Gut. 2021;70(10):1978-1988.
- Franzosa EA, Sirota-Madi A, Avila-Pacheco J, et al. Gut microbiome structure and metabolic activity in inflammatory bowel disease [published correction appears in Nat Microbiol. 2019 May;4(5):898]. Nat Microbiol. 2019;4(2):293-305.
- Carroccio A, Iacono G, Cottone M, et al. Diagnostic accuracy of fecal calprotectin assay in distinguishing organic causes of chronic diarrhea from irritable bowel syndrome: a prospective study in adults and children. Clin Chem. 2003;49(6 Pt 1):861-867.
- Alibrahim B, Aljasser MI, Salh B. Fecal calprotectin use in inflammatory bowel disease and beyond: A mini-review. Can J Gastroenterol Hepatol. 2015;29(3):157-163.
- Lundgren D, Eklöf V, Palmqvist R, Hultdin J, Karling P. Proton pump inhibitor use is associated with elevated faecal calprotectin levels. A cross-sectional study on subjects referred for colonoscopy. Scand J Gastroenterol. 2019;54(2):152-157.

19. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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CHORUS CALPROTECTIN 3RD GENERATION

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- Hoskin TS, Crowther JM, Cheung J, et al. Oxidative cross-linking of calprotectin occurs in vivo, altering its structure and susceptibility to proteolysis. Redox Biol. 2019;24:101202.
- Jukic A, Bakiri L, Wagner EF, Tilg H, Adolph TE. Calprotectin: from biomarker to biological function. Gut. 2021;70(10):1978-1988.
- Franzosa EA, Sirota-Madi A, Avila-Pacheco J, et al. Gut microbiome structure and metabolic activity in inflammatory bowel disease [published correction appears in Nat Microbiol. 2019 May;4(5):898]. Nat Microbiol. 2019;4(2):293-305.
- Carroccio A, Iacono G, Cottone M, et al. Diagnostic accuracy of fecal calprotectin assay in distinguishing organic causes of chronic diarrhea from irritable bowel syndrome: a prospective study in adults and children. Clin Chem. 2003;49(6 Pt 1):861-867.
- Alibrahim B, Aljasser MI, Salh B. Fecal calprotectin use in inflammatory bowel disease and beyond: A mini-review. Can J Gastroenterol Hepatol. 2015;29(3):157-163.
- Lundgren D, Eklöf V, Palmqvist R, Hultdin J, Karling P. Proton pump inhibitor use is associated with elevated faecal calprotectin levels. A cross-sectional study on subjects referred for colonoscopy. Scand J Gastroenterol. 2019;54(2):152-157.

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