CHORUS

ENA-6 S

REF 86012

REF 86012/12

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For In Vitro Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS ENA-6 S is an immunoassay kit for the automated qualitative determination of IgG class antibodies against 6 cellular and nuclear antigens (SS-A, SS-B, Scl-70, Sm, snRNP/Sm and Jo-1).

As ENA antibodies are widely used as a serologic marker of systemic autoimmune rheumatic disease, the kit is used as an aid to related diagnosis.

The test performed in human serum, using a disposable device applied on the chorus and CHORUS TRIO instruments, must be used by professional laboratory users only.

2. INTRODUCTION

Anti-nuclear antibodies (ANA) are an important tool for the differential diagnosis of systemic rheumatic diseases. Indirect immunofluorescence test (IFT) on eukaryotic cells like HeLa has been the established method for the detection of ANAs. Single antibody specificities are distinguished by fluorescence patterns but more specific testing by ELISAs employing the target antigens are available too for a simple and reliable differentiation of ANAs.

ANAs are especially found in active and inactive systemic lupus erythematosus (SLE), mixed connective tissue diseases (MCTD), scleroderma, Sjögren's syndrome, polymyositis.

Antibodies against:

- Sm (Smith antigen) are directed against core proteins (B,B', D1-D3, E,F,G) of small nuclear ribonucleoproteins (snRNPs). Anti-Sm as well as antibodies against double stranded DNA (dsDNA) are highly specific for SLE and thus are included in diagnostic and classification criteria for SLE.
- snRNP/Sm are directed against Sm and U1-snRNP proteins (70 kDa, A and C). They occur in SLE, Sjögren's syndrome, scleroderma and polymyositis.
- SS-A (Ro; soluble cytoplasmic and/or nuclear ribonucleoproteins of 52 kDa and 60 kDa) and antibodies against SS-B (La; 48 kDa protein associated with RNA polymerase III) are mainly found in high titers for primary and secondary Sjögren's syndrome but also in SLE, congenital heartblock and neonatal lupus.
- Scl-70 are directed against DNA-topoisomerase I. They are highly specific for systemic scleroderma and give a hint for a severe course.

- Jo-1 are directed against histidyl-tRNA synthetase (cytoplasmic protein involved in protein biosynthesis) and are found in 20-40% of patients with polymyositis and dermatomyositis.

3. PRINCIPLE OF THE METHOD

The CHORUS ENA-6 S device is ready to use for the detection of IgG antibodies against 6 cellular and nuclear antigens, in the CHORUS/CHORUS TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). The antigens are bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum. After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human immunoglobulins antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate is added. The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the CHORUS/CHORUS TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off) calculated in reference to the CDC Atlanta.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS/CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet available on DIESSE website: www.diesse.it for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.

- Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.

- Discard devices which show the substrate (well 4) blue colored.**
- Adding the sample into the well verify that it is perfectly distributed on the bottom.
- Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
- The devices are for use with the CHORUS/CHORUS TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
-
- Check that the CHORUS/CHORUS TRIO instrument is set up correctly (see Operating Manual).
- Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
- Avoid using self-defrosting freezers for the storage of the samples.
- Defective barcodes can be inserted manually in the instrument (see Operating Manual).
- Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
- Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the guidance in the chapter "Analytical Specificity").
- Do not use the device after the expiry date.
- Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86012)


The kit is sufficient for 12 tests. (REF 86012/12)

DD DEVICES

6 packages each containing 6 devices (REF 86012)

2 packages each containing 6 devices (REF 86012/12)

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: saline proteic solution with Proclin (0.1%)

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies 18-13 (0.6 µg/ml-0.075 µg/ml concentrated) labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%

Position 1: EMPTY WELL

in which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-SS-A, SS-B, Scl-70, Sm, snRNP/Sm and Jo-1 and preservative (Proclin = 0.0%, Tween-20 = 0.2%, methylorange = 7.5 µg/ml) . Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti-SS-A, SS-B, Scl-70, Sm, snRNP/Sm and Jo-1 and preservative Proclin = 0.1%, Tween-20 = 0.2%, methylorange = 7.5 µg/ml) .. Liquid, ready for use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human antibodies against 6 cellular and nuclear antigen in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different antibodies against 6 cellular and nuclear antigens levels.

MATERIALS REQUIRED BUT NOT PROVIDED:

- WASHING BUFFER AUTOIMUNITY REF 86004
- CLEANING SOLUTION 2000 REF 83609
- SANITIZING SOLUTION REF 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT REF 83607
- CHORUS/CHORUS TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for longer periods at temperature ≤-20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.

3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the CHORUS/CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Customer care

Tel: 0039 0577 319554
email: scientificsupport@diesse.it
customercare@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS/CHORUS TRIO instrument expresses the result in

Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALITICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)
Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (10 mg/dl - 250 mg/dl)
Hemoglobin (5 mg/ml - 30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

33 samples, positive to dsDNA, SS-A, SS-B, Sm, Scl-70, Jo-1, U1-RNP, U3-RNP, Centromere Pattern FANA, PM1, rRNP/Ribosomal P, MPO-ANCA, PR3-ANCA, Gliadin A, Gliadin G, Transglutaminase IgA, Cardiolipin IgG, Cardiolipin IgM, β-2-

Glycoprotein IgG, β -2-Glycoprotein IgM, Anti-Tg and Anti-TPO were tested.

No significant cross-reactions were found.

14. METHOD COMPARISON

In an experimentation 482 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Tot.
Diesse	+	155	7	162
	-	8	312	320
	Tot.	163	319	482

Percent Positive Agreement (~Diagnostic Sensitivity):

95.1% $CI_{95\%}$: 90.6- 97.5

Percent Negative Agreement: (~Diagnostic Specificity):

97.8% $CI_{95\%}$: 95.5- 98.9

Positive Predictive Value (PPV) =95,7% $CI_{95\%}$ =93,9-97,5

Negative Predictive Value (NPV) = 97,5% $CI_{95\%}$ =96,1-99,2

The agreement between the two methods is excellent with a Cohen's Kappa of 0.93.

15. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean Index	CV%	Mean Index	CV%
1	0.6	8.3	0.6	8.3
2	0.7	7.1	0.6	8.3
3	0.6	6.7	0.5	10.0
4	0.4	10.0	0.5	10.0
5	0.8	8.8	0.9	7.8
6	1.0	8.0	1.0	9.0
7	1.2	9.2	1.1	6.4
8	1.2	5.0	1.2	7.5
9	3.1	10.0	3.5	7.1
10	4.7	8.9	4.5	9.1
11	5.7	5.1	5.5	4.5
12	7.9	9.5	7.7	5.5

Sample	Between lots		Between Instruments	
	Mean Index	CV%	Mean Index	CV%
1	0.6	-	0.6	10.0
2	0.6	-	0.6	-
3	0.5	-	0.5	-
4	0.5	12.0	0.4	15.0
5	0.9	-	0.9	-
6	1.0	6.0	1.0	6.0
7	1.2	-	1.2	-
8	1.1	-	1.1	5.5
9	3.2	4.7	3.1	1.9
10	4.7	3.2	4.7	-
11	5.7	6.1	5.7	2.6
12	7.5	6.7	7.5	1.3

16. REFERENCES

1. Antinuclear antibody. The Lancet 1984: Sept. 15, 611-13.
2. Froelich C.H., Wallmann H., Skosey J.I., Teodorescu M. (1990) J. Rheumatology 17: 192.

3. Mierau R., Genth E. (1998). In: Thomas L. Labor und Diagnose. TH-Books, Frankfurt, 15. Auflage: 843-851.
4. Schmolke M., Oppermann M., Helmke K., Guder WG (2000) Poster P59, 5th Dresden Symposium on Autoantibodies.
5. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38
6. Gabriela Riemekasten, Jens Y. Humrich, Falk Hiepe, 28 - Antibodies Against ENA (Sm, RNP, SSA, SSB), Editor(s): Daniel J. Wallace, Bevra Hannahs Hahn, Dubois' Lupus Erythematosus and Related Syndromes (Ninth Edition), Elsevier, 2019.
7. Karen J. Tietze, Chapter 5 - Review of Laboratory and Diagnostic Tests, Editor(s): Karen J. Tietze, Clinical Skills for Pharmacists (Third Edition), Mosby, 2012.

17. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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