CHORUS

**MYCOPLASMA
PNEUMONIAE**

IgA

REF 81033

DIESSE

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16. BIBLIOGRAFIA

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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For In Vitro Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS Mycoplasma pneumoniae IgA (REF 81033) is an immunoassay kit for the automated semi-quantitative determination of IgA antibodies against Mycoplasma pneumoniae.

Mycoplasma Pneumoniae is the most common etiological agent causing pneumonia acquired in community environments. The IgA represent the highest sensitivity marker in case of an infection in progress, a reinfection, or a recent infection; therefore, the kit is used as an aid in the diagnosis of pneumonia infection.

The test, performed in human serum using a disposable device applied to the CHORUS and CHORUS TRIO instruments, must be used by professional laboratory personnel only.

2. INTRODUCTION

Mycoplasma pneumoniae is the most common etiological agent causing pneumonia acquired in community environments, particularly in the age range of 5 to 30 years; it can cause epidemics which develop slowly as the incubation period varies from 10 to 14 days, and the infection involves subjects living in close contact or in segregated groups (schools, military barracks, family groups). Pneumonia caused by Mycoplasma is also known as atypical primary pneumonia or Eaton's pneumonia.

Mycoplasma pneumoniae attacks and destroys ciliated epithelial cells of the respiratory tract mucosa. Microscopic examination evidences interstitial pneumonia, bronchitis and bronchiolitis.

Due to the recurrence of symptoms caused by different etiological agents, additional diagnostic measures such as serological tests are necessary for the diagnosis of an acute pneumonia infection.

The immune response, following a Mycoplasma pneumoniae infection, depends on the type of infection: the IgA-class antibodies are more frequently found in case of primary infections as demonstrated by their presence in younger patients. The IgM are low or not found in older patients with higher probability of reinfection. On the other hand the IgA represent the highest sensitivity marker in case of an infection in progress, a reinfection or a recent infection.

Specific IgA-class antibodies appear early during the infection and reach peak levels during the first 4 weeks, then declining

quite rapidly, earlier than the IgG; diagnosis of an acute infection can often be made on the basis of a single sample test. The IgG appear after the IgM and IgA and they generally reach the peak level during the 5° week. A high titer together with a significant increase between two samples taken approximately 2 weeks apart, confirms that an infection is in progress.

The antigen used in the test for the serological determination of the specific IgA derives from an extract of M. pneumoniae membrane. This membrane contains significant quantities of P1 adhesin that represents the immunodominant antigen and corresponds to a transmembrane protein which is the main responsible of the adherence of M. pneumoniae to the host respiratory epithelium.

3. PRINCIPLE OF THE METHOD

The test is based on the indirect ELISA (Enzyme Linked Immunosorbent Assay) principle.

The antigen, composed of an extract of inactivated Mycoplasma pneumoniae is bound to the solid phase and is incubated with the human serum sample so that the specific IgA are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate constituted of anti-human IgA monoclonal antibodies conjugated to peroxidase.

The conjugate which has not reacted is eliminated and the peroxidase substrate is added. The blue colour which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents necessary to perform the test when applied on the Chorus instruments.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instruments.

4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) for at least 30 minutes before the use; use within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
The Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website (<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument and do not use devices with defective labels.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
10. Do not use hemolyzed, lipemic or, icteric samples with a concentration of interferents higher than that tested (according to the indications given in the "Analytical specificity" chapter).
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer REF 83606.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 determinations.

DD DEVICES 6 packages each containing 6 devices

Description of device:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Mycoplasma pneumoniae antigen: maximum concentration 2.51 µg/ml

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgA (maximum concentration 2 µg/m) labelled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid, ready for use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human antigens in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different antigens levels.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER REF 83606
- CLEANING SOLUTION 2000 REF 83609
- SANITIZING SOLUTION REF 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT REF 83607
- CHORUS/CHORUS TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated, and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice.

According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at temperature ≤ -20°C, a least 32 months, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used. The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.

2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
e-mail: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument gives a semiquantitative result in AU/mL.

The test serum can be interpreted as follows:

POSITIVE: when the result is >18 AU/mL
NEGATIVE: when the result is < 12 AU/mL
DOUBTFUL: when the result is in the range 12-18 AU/mL

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample after 1-2 weeks.

11. LIMITATIONS

The test results should be used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

In order to obtain a better evaluation of the antibody profile of the sample it is advisable to detect the IgM and IgG concentration as well.

Samples collected during the early phase of the infection may not have developed sufficient antibodies to be detectable; in the case of doubtful or negative results, collect a new sample after 2-3 weeks and repeat the test if clinical suspects persist.

It is advisable to centrifuge or filter lipemic or cloudy samples before testing.

Interpretation of the results based on the combined assay of IgG, IgM and IgA:

M. pneumoniae antibody level			
IgG	IgM	IgA	Interpretation
negative	negative	negative	No indication of infection
negative or positive	positive	negative or positive	Indication of infection in progress
positive	negative	negative	Indication of

			past infection
negative or positive	negative	positive	Indication of infection in progress

12. CALIBRATION RANGE

Calibration range 0-100.0 AU/ml.

For samples > 100.0 AU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – not supplied with the kit).

13. CROSS-REACTIONS

2 samples (1cut-off and 1negative) containing potentially interfering substances were tested:

Bilirubin (15 mg/dl – 60mg/dl)

Triglycerides (7.5 mg/ml – 30 mg/ml)

Hemoglobin (2.5 mg/ml – 10 mg/ml)

RF (30 IU/ml – 150 IU/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

14. METHOD COMPARISON

In an experimentation, 102 samples were analyzed with the Diesse kit as well as with another commercial method.

Below are the schematized results of the trial:

		Reference		
		+	-	Total
Diesse	+	20	8	28
	-	0	74	74
	Total	20	82	102

Diagnostic Sensitivity: 100.0 % CI_{95%}: 83.9 – 100.0

Diagnostic Specificity: 90.2 % CI_{95%}: 81.9 – 95.0

Positive Predicted Value (PPV) = 71.4 CI_{95%}: 62.7 – 80.2

Negative Predicted Value (NPV) = 100.0 CI_{95%}: 100.0 – 100.0

15. PRECISION AND REPEATABILITY

WITHIN-RUN PRECISION

	No. of replicates	Mean AU/ml	S.D.	CV%
Sample 1	5	10.0	0.00	0.0
Sample 2	5	20.5	2.47	12.0
Sample 3	4	69.4	3.46	5.0
Sample 4	4	35.7	4.81	13.5

BETWEEN-RUNS PRECISION

		Mean AU/ml			Mean	S. D.	CV%
		Day 1	Day 2	Day 3			
Sample 1		20.3	20.5	24.1	21.6	2.14	9.9

Sample 2	48.9	43.9	44.7	45.8	2.69	5.9
Sample 3	11.9	10.3	10.0	10.7	1.02	9.5
Sample 4	10.0	10.0	10.0	10.0	0.00	0.0
Sample 5	55.0	55.4	57.8	56.1	1.51	2.7
Sample 6	18.2	16.4	17.5	17.4	0.91	5.2

16. REFERENCES

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

17. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

18. SUMMARY OF SAFETY AND PERFORMANCE

This document, which will be made available on the EUDAMED database (when this is fully implemented and functioning), is part of the Technical Documentation and can be requested from the manufacturer.

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16. BIBLIOGRAPHIE

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).

2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of *Mycoplasma pneumoniae*. *Microbiology* (2000),146,741-747
3. S. Rotten: Interaction of *Mycoplasma* With Host Cells. *Physiol.Rev.*83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of *Mycoplasma pneumoniae* Infections. *J. Clin. Microbiol.* 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of *Mycoplasma pneumoniae* by Fractionated Solubilization and size exclusion Chromatography. *Biol.Chem* Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIE

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev. 83:417-432, 2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAPHIE

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol. Rev. 83: 417-432, 2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40: 165-171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol. Chem. Vol. 369, pp. 1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFÍA

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of *Mycoplasma pneumoniae*. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of *Mycoplasma* With Host Cells. Physiol. Rev. 83: 417-432, 2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of *Mycoplasma pneumoniae* Infections. J. Clin. Microbiol. 2002 40: 165-171
5. E. Jacobs: Isolation of the Adherence Protein of *Mycoplasma pneumoniae* by Fractionated Solubilization and size exclusion Chromatography. Biol. Chem. Vol. 369, pp. 1295-1299, December 1988.

6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIA

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of *Mycoplasma pneumoniae*. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of *Mycoplasma* With Host Cells. Physiol.Rev.83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of *Mycoplasma pneumoniae* Infections. J. Clin. Microbiol. 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of *Mycoplasma pneumoniae* by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIJA

- G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
- S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
- J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40:165-171
- E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
- CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIA

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40:165-171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

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16. VIITED

1. G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741–747
3. S. Rotten: Interaction of Mycoplasma With Host Cells. Physiol. Rev. 83: 417–432, 2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J. Clin. Microbiol. 2002 40: 165–171
5. E. Jacobs: Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol. Chem. Vol. 369, pp. 1295–1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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