**RPR-DAT
ADVANCED**

DIESSE

REF 26031

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Percent Positive Agreement (~Sensibilità Diagnostica): 100.0%
CI95%: 96.2-100.0

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17. BIBLIOGRAFIA

1. Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1-21
2. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
3. WHO guidelines for the treatment of Treponema pallidum (Syphilis) (2016)
4. Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2
5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

RPR-DAT ADVANCED (REF 26031) is a non-treponemal agglutination assay for the qualitative and semi-quantitative detection of reagines in human serum to be applied to the AUTO-DAT (REF 26001) instrument.

Because non-treponemal tests are based on the detection of reagines, a class of antibodies found in Syphilis (a disease caused by *Treponema pallidum* infection), the kit is intended to determine exposure to *Treponema pallidum* infection and to be used as an aid in its diagnosis.

The test is not intended to be used to detect exposure to *Treponema pallidum* infection in blood, blood components, cells, tissues, organs, or any derivatives thereof in order to assess their suitability for transfusion, transplantation, or cell administration.

It should only be used by professional laboratory personnel.

2. INTRODUCTION

Syphilis, caused by *Treponema pallidum* (subsp. *pallidum*) spirochetes, is a chronic infection with many different clinical manifestations occurring in distinct stages. Generally, this systemic infectious disease is contracted through direct sexual contact and has lesions containing treponemes. The only known hosts are humans.

Many tests are available for direct and indirect diagnosis of syphilis, but there is still no single optimal test.

Direct diagnostic methods include detection of *T. pallidum* by microscopic examination of fluid or smears from lesions, histological examination of tissues, or nucleic acid amplification methods such as polymerase chain reaction (PCR).

Indirect diagnosis is based on serologic tests for the detection of antibodies. Serological tests fall into two categories:

1. nontreponemal screening tests
2. confirmatory treponemal tests.

Nontreponemal tests are based on the detection of reagins, a class of antibodies found in syphilis and occasionally in other acute and chronic conditions. Reagins can be detected in serum about 4-6 weeks after infection, or 1-3 weeks after the appearance of primary syphiloma.

Purified beef heart extract (cardiolipin) fortified with lecithin and cholesterol, is used as an antigen for reagin tests. Nontreponemal tests are routinely used in syphilis serology because of their practical advantages and reproducibility, although they are not always 100% specific.

Since the product RPR-DAT ADVANCED is a non-treponemal test, based on the flocculation principle, it is intended to determine exposure to *Treponema pallidum* infection as an aid in the related diagnosis.

3. PRINCIPLE OF THE METHOD

The test is based on the principle of flocculation.

The serum sample is reacted in the microplate well with a colloidal suspension of cardiolipin, lecithin and cholesterol mixed with carbon microparticles.

In the presence of antibodies (reagins), which agglutinate the antigen, an aggregate is formed (positive reaction). The reaction is negative if a circle appears inside the microplate well.

The test is applied to the AUTO-DAT instrument (REF 26001), which captures the image of the reaction that has occurred, analyzes it, provides and prints a result calculated on the basis of a specific algorithm.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry.

All materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use, bring the reagents to room temperature (18-30°C).

1. The test should be used together with the AUTO-DAT instrument (REF 26001), strictly following the Instructions for Use and User's Manual of the instrument.
2. Do not change the Test Procedure. Do not substitute reagents with those from other suppliers or other batches unless it is specifically stated that the reagent is interchangeable between batches.
3. Check that the instrument is set up correctly (see User's Manual).
4. Avoid the use of self-defrosting freezers for sample storage.

5. Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the directions in the chapter "Analytical Specificity").
6. Do not use the device after the expiry date.
7. After use, store reagents at 2-8°C.
8. Avoid microbial pollution of reagents as this reduces the validity of the product and may result in erroneous results.

5. KIT COMPOSITION AND REAGENT PREPARATION

The reagents are sufficient for 500 determinations.

ANTIGEN

ANTIGEN

2 x 4 mL

Content: The antigen is composed of a stabilized mixture of cardiolipin 0.003%, cholesterol 0.09%, lecithin 0.02% and carbon microparticles 0.2 g/L

Caution: Store vials in an upright position when not in use. Bring to room temperature before use. Do not shake vigorously.

CONTROL +

POSITIVE CONTROL

1 x 0.5 mL

Contents: Human serum reactive for syphilis, diluted up to a maximum of 15% in a protein solution with 0.09% sodium azide.

Liquid, ready to use.

CONTROL -

NEGATIVE CONTROL/SAMPLE DILUENT

1 x 4 mL

Content: 100% non-reactive human serum containing Proclin 300 0.029% and gentamicin 0.05% as preservatives.

Liquid, ready-to-use

MATERIALS REQUIRED BUT NOT PROVIDED

- AUTO-DAT INSTRUMENT (REF 26001)
- AUTO-DAT MICROPLATES RPR (REF 26006)
- Micropipettes capable of accurately taking volumes of 15-25 µl.
- Disposable gloves
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

After first use, microplates should be stored at room temperature.

The expiry date is printed on each component and on the kit label.

The stability of the reagents does not change after opening the bottle, as long as the user takes care to keep the product away from possible microbial contamination.

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and

complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum can be kept for 7 days at 2/8°C; for longer storage periods, freeze at -20°C for at least 10 months.

The sample can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The operator can choose how to carry out the analysis:

- Qualitative test
- Semi-quantitative test

The evaluation of the result is done by observing the agglutination of the sample in the single well, in the case of a qualitative test, or in the pattern of wells, in the case of a semi-quantitative test.

In any case, bring reagents and samples to room temperature. Resuspend the reagent.

QUALITATIVE TEST

1. Add 25 µl of test serum to the bottom of the microplate well. Be careful not to transfer any cellular items. **Avoid the formation of bubbles.**
2. Repeat the step 1 for each sample to be tested.
3. Add 15 µL of the antigen in each well. **Avoid the formation of bubbles.**
4. Verify that the sample-antigen mixture is **evenly distributed at the bottom of the well**. Otherwise, tap the plate on the side and/or tap gently on the work surface.
5. Place the plates in the AUTO-DAT (REF 26001) tool and follow the instructions.
6. At the end of the analysis, remove the plates from the AUTO-DAT instrument and discard the wells where the reaction occurred. Previously unused wells can be used for later analysis.

NOTE: the positive and negative controls are supplied ready to use and do not require pre-dilution. They must be tested at each session like the samples under examination.

In case of a positive sample, it is advisable to carry out the semi-quantitative analysis.

SEMI-QUANTITATIVE TEST

Perform serial dilutions of the test sample from indilute to 1/128:

1. Add 25 µL of Negative Control/Diluent from B1 to H1.
2. Add 25 µl of test serum to well A1 and well B1, stirring thoroughly. Be careful not to transfer any cellular items. **Avoid the formation of bubbles.**
3. Take 25 µl from B1 and transfer it to the C1 well. **Mix the solution before transferring. Avoid the formation of bubbles.**
4. Perform serial dilutions: Repeat step 3 for subsequent wells up to H1. **Discard 25 µl from the last well.**
5. Follow the procedure described for the qualitative test in step 3 to step 6.

Repeat the operation for all samples to be analyzed using the other columns on the plate.

	step 1	step 2	step 3	
	1	1	1	DIL
A		25 µl sample	15 µl ANTIGEN	1/1
B	25 µl CONTROL -	25 µl sample	15 µl ANTIGEN	1/2
C	25 µl CONTROL -	+25 µl B1	15 µl ANTIGEN	1/4
D	25 µl CONTROL -	+25 µl C1	15 µl ANTIGEN	1/8
E	25 µl CONTROL -	+25 µl D1	15 µl ANTIGEN	1/16
F	25 µl CONTROL -	+25 µl E1	15 µl ANTIGEN	1/32
G	25 µl CONTROL -	+25 µl F1	15 µl ANTIGEN	1/64
H	25 µl CONTROL -	+25 µl G1	15 µl ANTIGEN	1/128

DISCARD 25 µl

NOTE: the positive and negative controls are supplied ready to use and do not require pre-dilution. They must be tested at each session like the samples under examination.

The titer will be given by the last dilution found reactive with the instrument.

9. TEST VALIDATION

Use the controls at every session.

Proceed as described in the "PROCEDURE" section.

If the result is different than expected, contact Customer Care.

Tel: 0039 0577 319554

email: scientificsupport@diesse.it;
customercare@diesse.it

10. INTERPRETATION OF THE RESULTS

For semi-quantitative test, the instrument gives the result in titer.

The test on the test sample can be interpreted as follows:

POSITIVE: when the result is $\geq 1/1$

NEGATIVE: when the result is NEG

The instrument can detect the presence of pro-zone effect with a limit of 1/8; refer to the user manual.

For qualitative test, the instrument gives the result as follows:

POSITIVE (POS): Reactive sample

NEGATIVE (NEG): Non-reactive sample

It is advisable to visually verify the results obtained to confirm the report.

POSITIVE

NEGATIVE

11. LIMITATIONS

The product should be used only by professional laboratory personnel.

The test is not suitable for samples different from human serum. The test is only applicable to the AUTO-DAT instrument (REF 26001).

Treponemal tests become positive 5-15 days after the onset of symptoms. No single test or reference standard is available for each stage of the disease.

In qualitative test, the presence of false negative results due to pro-zone effect of highly positive samples cannot be excluded, to exclude false negativity for pro-zone samples should be tested by the semi-quantitative procedure.

All values obtained need careful interpretation that does not prescind from other indicators related to the same patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. REFERENCE RANGE

SEMI-QUANTITATIVE TEST : 1/1 – 1/128 (Titer)

13. ANALYTICAL SPECIFICITY

3 samples (1 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (4.5 mg/dl – 45 mg/dl)

Triglycerides (250 mg/dl – 1500 mg/dl)

Hemoglobin (2.5 mg/ml – 10 mg/ml)

RF (44 – 220 IU/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

14. CROSS-REACTIONS

18 samples positive to Borrelia (3), Herpes virus (3), Toxoplasma (3), Cytomegalovirus (3), Epstein-Barr Virus (3) and Rubella (3) were tested.

No significant cross-reactions were found.

15. METHOD COMPARISON

In an experimentation 185 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	98	3	101
	-	0	84	84
	Total	98	87	185

Percent Positive Agreement (~Diagnostic Sensitivity): 100.0%
CI95%: 96.2-100.0

Percent Negative Agreement: (~Diagnostic Specificity): 96.6%
CI95%: 90.3-98.8

Positive Predictive Value (PPV): 97.0% CI95%: 94.5-99.5

Negative Predictive Value (NPV): 100.0% CI95%: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.96.

16. PRECISION AND REPEATABILITY

Precision and repeatability has been evaluated on replicates of negative and positive samples, in three different batches, performed in different days with different instruments.

All the obtained results were within the acceptance criteria (Agreement % $\geq$ 95% between the expected and the obtained results) demonstrating the precision and repeatability of the method.

17. REFERENCES

1. Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1–21
2. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
3. WHO guidelines for the treatment of Treponema pallidum (Syphilis) (2016)
4. Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2
5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

18. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

19. SUMMARY OF SAFETY AND PERFORMANCE

This document, which will be made available on the EUDAMED database (when this is fully implemented and functioning), is part of the Technical Documentation and can be requested from the manufacturer.

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17. BIBLIOGRAPHIE

1. Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1–21

2. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
3. WHO guidelines for the treatment of Treponema pallidum (Syphilis) (2016)
4. Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2
5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

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17. BIBLIOGRAPHIE

1. Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1–21
2. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
3. WHO guidelines for the treatment of *Treponema pallidum* (Syphilis) (2016)
4. Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2
5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

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