CHORUS

Calprotectin

REF 86600

DIESSE

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Chorus/Chorus TRIO: necessaria Software Release 3.2.1 e successive

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16. BIBLIOGRAFIA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, IO-09/352-C IFU 86600 – Ed. 28.10.2016

- simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol.* 2014 Apr;49(4):418-23.
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol.* 1996 Apr;31(4):339-44.

INSTRUCTIONS FOR USE

CHORUS Calprotectin

**For the quantitative determination of
Calprotectin**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of calprotectin in stools samples, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Calprotectin is a protein present in neutrophils and its high concentration in the stools is a sign of inflammation. High levels are observed typically in patients with UC (ulcerative colitis), Crohn's disease (CD), colic, gastric and oesophageal neoplasms. Further, an increase of values can be observed in all the pathologies that involve the activation of an acute or chronic inflammatory process of the gastrointestinal system. Anti-inflammatory medications may vary the Calprotectin levels.

3. PRINCIPLE OF THE METHOD

The Chorus Calprotectin device is ready to use for the detection of calprotectin in stools samples in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). The polyclonal antibodies are bound to the solid phase. Calprotectin is bound to the solid phase through incubation with the extract of faeces. After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-calprotectin antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated and the peroxidase substrate is added. The colour which develops is proportional to the concentration of specific antibodies present in the stools sample.

The disposable devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in mg/Kg.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory

practice should be followed when handling reagents and samples.

Waste disposal: samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument (**Software Release 3.2.1 and subsequent ones**), the Instructions for Use and the Instrument Operating Manual must be carefully followed.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).


9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

coated with polyclonal anti-calprotectin antibodies

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: buffer solution at pH 7.9

Position 2: CONJUGATE

Contents: anti-calprotectin IgG antibodies labeled with horseradish peroxidase, in buffer solution.

Position 1: EMPTY WELL

in which sample must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2/8°C.

CALIBRATOR CALIBRATOR 1 x 0.3 ml

Contents: calprotectin. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.6 ml

Contents: calprotectin. Liquid, ready for use.

SAMPLE DILUENT SAMPLE DILUENT 1 x 1.8 ml

Contents: diluent ready to use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Chorus COPROCOLLECT **REF** 86602
- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)

- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
SAMPLE DILUENT	3 months at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of extract of faeces, prepared as reported in the Instructions for Use of the product Chorus Coprocollect (**REF** 86602) and handled with all precautions dictated by good laboratory practice.

The sample may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in mg/Kg, calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 60 mg/Kg

NEGATIVE: when the result is < 40 mg/Kg

EQUIVOCAL: for all values between 40 and 60 mg/Kg

If the result is equivocal, repeat the test.

Levels detected in newborns might be significantly higher.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

Chorus/Chorus TRIO: necessary Software Release 3.2.1 and subsequent ones.

12. CALIBRATION RANGE

Calibration range 12.5 - 250 mg/Kg.

For results > 250 mg/Kg repeat the test on the sample diluted 1:5 with the ready to use sample diluent. Multiply by 5 the result obtained by the instrument.

13. REFERENCE RANGE

As reported by literature, normal values of calprotectin in the faeces correspond to 25 mg/Kg.

14. METHOD COMPARISON

In an experimentation 42 samples have been tested with DIESSE kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Tot.
DIESSE	+	23	0	23
	-	1	18	19
	Tot.	24	18	42

Percent Positive Agreement (~Diagnostic Sensitivity):

95.8% CI_{95%}: 79.7-99.1

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI_{95%}: 82.4-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 0.96.

15. PRECISION AND REPEATABILITY

Sample	Within run	
	Mean (mg/Kg)	CV%
1	<12.5	-
2	52.2	4.3
3	105.0	5.4
4	47.1	4.2
5	63.8	8.0

6	55.9	5.6
7	95.1	12.2
8	67.0	4.2

Sample	Between run	
	Mean (mg/Kg)	CV%
1	<12.5	-
2	49.7	5.8
3	101.2	9.5
4	43.9	13.7
5	74.5	14.0

Sample	Between lots		Between Instruments	
	Mean (mg/Kg)	CV%	Mean (mg/Kg)	CV%
1	12.6	0.8	13.6	14.9
2	54.9	9.8	51.5	8.0
3	119.5	10.6	109.4	12.1
4	55.4	14.7	50.5	14.5
5	87.9	7.3	84.6	14.4

16. REFERENCES

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamara Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. Scand J Gastroenterol. 2014 Apr;49(4):418-23

11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. Scand J Gastroenterol. 1996 Apr;31(4):339-44.

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1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23–36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct; 47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. Scand J Gastroenterol. 2014 Apr; 49(4):418-23
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. Scand J Gastroenterol. 1996 Apr; 31(4):339-44.

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16. LITERATUR

- Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
- Tamara Pang, Steven T. Leach, Tamara Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
- Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
- J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
- Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.

6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol* 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr* 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol*. 2014 Apr;49(4):418-23
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol*. 1996 Apr;31(4):339-44.

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3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease *Ther Adv Gastroenterol* 2015, Vol. 8(1) 23–36.
4. J Tibble, K Teahon, B Thjodleifsson, A Røseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. *Gut*. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. *Scand J Gastroenterol* 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol* 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr* 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol*. 2014 Apr;49(4):418-23
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol*. 1996 Apr;31(4):339-44.

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16. BIBLIOGRAFÍA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.

9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol.* 2014 Apr;49(4):418-23
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol.* 1996 Apr;31(4):339-44.

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16. BIBLIOGRAPHIE

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.

5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. *Scand J Gastroenterol* 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol* 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr* 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol*. 2014 Apr;49(4):418-23
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol*. 1996 Apr;31(4):339-44.

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16. BIBLIOGRAFIA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamara Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic

- diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
 10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol*. 2014 Apr;49(4):418-23
 11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol*. 1996 Apr;31(4):339-44.

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Percent Positive Agreement (~diagnostiline tundlikkus):

95.8% CI_{95%}: 79.7-99.1

Percent Negative Agreement: (~diagnostiline spetsiifilisus):

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16. BIBLIOGRAAFIA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 Jan 14; 20(2): 363-375
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day ja Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri ja Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798
6. Røseth A.G. et al: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.
9. Sohnle P.G. et al: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. Scand J Gastroenterol. 2014 Apr;49(4):418-23.
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. Scand J Gastroenterol. 1996 Apr;31(4):339-44.

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Percent Positive Agreement (~ diagnostinis jautrumas):

95,8 % CI₉₅ %: 79,7–99,1

Percent Negative Agreement: (~ diagnostinis specifiškumas):

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16. BIBLIOGRAFIJA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23–36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers

- in ulcerative colitis during high-dose steroid treatment. Scand J Gastroenterol. 2014 Apr;49(4):418-23.
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. Scand J Gastroenterol. 1996 Apr;31(4):339-44.

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16. BIBLIOGRAFIA

1. Jan Däbritz, Jason Musci, Dirk Foell. World J Gastroenterol 2014 January 14; 20(2): 363-375.
2. Tamara Pang, Steven T. Leach, Tamarah Katz, Andrew S. Day and Chee Y. Ooi. Fecal biomarkers of intestinal health and disease in children January 2014.
3. Frank S. Lehmann, Emanuel Burri and Christoph Beglinger. The role and utility of faecal markers in inflammatory bowel disease Ther Adv Gastroenterol 2015, Vol. 8(1) 23-36.
4. J Tibble, K Teahon, B Thjodleifsson, A Roseth, G Sigthorsson, S Bridger, R Foster, R Sherwood, M Fagerhol, and I Bjarnason. A simple method for assessing intestinal inflammation in Crohn's disease. Gut. 2000 Oct;47(4):506-13.
5. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992; 27: 793-798.
6. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.

7. Bunn S.K. et al.: Fecal calprotectin: Validation as a non invasive measure of bowel inflammation in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr* 2001; 33: 14-22.
8. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95: 2831-2837.
9. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. *Journal of Infectious Diseases* 1991; 164: 137-142.
10. Theede K1, Kiszka-Kanowitz M, Nielsen AM, Nordgaard-Lassen I.: The correlation between fecal calprotectin, simple clinical colitis activity index and biochemical markers in ulcerative colitis during high-dose steroid treatment. *Scand J Gastroenterol*. 2014 Apr;49(4):418-23.
11. Meling TR1, Aabakken L, Røseth A, Osnes M.: Faecal calprotectin shedding after short-term treatment with non-steroidal anti-inflammatory drugs. *Scand J Gastroenterol*. 1996 Apr;31(4):339-44.

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