# ENZY-WELL

## HERPES SIMPLEX 2 IgG RECOMBINANT

Prodotto da/  
Manufactured by:

**REF 91029**  
**(96 tests)**

**REF 91183**  
**(960 tests)**

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### 19. BIBLIOGRAFIA

1. Arvaja M., Lehtinen M., Koskela P., Lappalainen M., Paavonen J., Vesikari T.; "Serological evaluation of Herpes Simplex Virus type 1 and type 2 infections in pregnancy"; Sexual Transmission Infections , 1999; 75: 168-171
2. Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
3. Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L.; "Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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## INSTRUCTIONS FOR USE

### ENZY-WELL HERPES SIMPLEX 2 IgG RECOMBINANT

**For the qualitative determination of IgG-class antibodies to Herpes Simplex Virus (type 2)**

**For In Vitro Diagnostic Use Only**

#### 1. INTENDED USE

Immunoenzymatic kit for the qualitative determination of IgG-class antibodies to Herpes Simplex Virus (type 2) in human serum.

#### 2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cherratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgG is important to establish the serological state of the patient.

#### 3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme Linked ImmunoSorbent Assay).

The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum. After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue color which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, the yellow color which develops can be easily read using an ELISA microplate reader.

#### 4. WARNINGS AND PRECAUTIONS

##### FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.**

##### Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. The following reagents contain low concentrations of harmful or irritant substances:
  - a) The Wash Buffer contains detergents
  - b) The Conjugate contains phenol
  - c) The Substrate is acid
  - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
5. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
6. Sulphuric acid required for the Stop Solution and 2 M hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
7. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
8. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

### Analytical Precautions

1. The OD of calibrator, controls and samples can be slightly different among different plates. For such reason if during the same run strips from different plates are used, even if the lot is the same, it is necessary to repeat the determination of the calibrator.
2. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.**
3. Open the package containing the strips after 30 minutes at room temperature at least.
4. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
5. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
6. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
7. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
8. Do not allow wells to become dry during the assay procedure.
9. Do not keep the samples in auto-defrosting freezers.
10. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the substrate and conjugate.
11. Care should be taken to avoid touching or splashing the rim of the well with conjugate.
12. Do not "blow-out" from microplates.
13. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO2 incubators must not be used.
14. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
15. Use of highly hemolyzed samples or samples containing bilirubin, incompletely clotted sera, samples with microbial contamination may give rise to erroneous results.
16. Care should be taken to avoid contaminating the microplate wells with disposable gloves powder.
17. The use of the kit with automated equipment has to be validated by the user.
18. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:

- installation and particular requisites
- operating principles, instructions, precautions and risks
- manufacturer's specifications and instrument performance
- servicing and maintenance

### 5. KIT COMPOSITION AND REAGENT PREPARATION

Reagents are sufficient for 96 determinations.

Reagents are sufficient for 960 determinations.

**Bring to room temperature before use.**

#### **MT PLATE** MICROPLATE (PF93029)

12x8 (REF 91029)

10x96 (REF 91183)

Content: microplate (96 wells) coated with HSV 2 IgG recombinant antigen, composed of recombinant G2 glycoprotein specific for Herpes Simplex Virus Type 2.

Use: open the package at the opposite end from the code (H2R followed by the lot number) which is useful for identification purposes, remove the support and strips required for the assay from the foil package and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

#### **CONJ** CONJUGATE (PF93529)

1x16mL (REF 91029)

8x16 mL (REF 91183)

Content: anti-IgG monoclonal antibodies, labeled with peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%. Liquid, ready for use.

#### **CONTROL +** POSITIVE CONTROL (PF93829)

1x1.6 mL (REF 91029)

4x1.6 mL (REF 91183)

Content: Human serum, containing a known concentration of anti-HSV 2 IgG antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use.

#### **CONTROL IgG -** IgG NEGATIVE CONTROL (PF93910)

1x 1.6 mL (REF 91029)

4x1.6 mL (REF 91183)

#### **INTERCHANGEABLE BETWEEN LOTS**

Content: Human serum, not containing anti-HSV 2 IgG antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use.

#### **CONTROL CUT-OFF** CUT OFF CONTROL (PF91829)

1x2 mL (REF 91029)

5x2mL (REF 91183)

Content: Human serum, containing a known concentration of anti-HSV 2 IgG antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use.

**SAMP DIL 3** DILUENT 3 (PF30029)

1x100 mL (REF 91029)

6x100 mL (REF 91183)

**INTERCHANGEABLE BETWEEN LOTS**

Content: Proteic solution in phosphate buffer with sodium azide 0.09% containing methyl orange as dye and cells extracts.

Use: To be used to dilute samples. Liquid, ready for use.

**WASH BUF 10x** WASH BUFFER 10X (PF93603)

1x100 mL (REF 91029)

5x100 mL (REF 91183)

**INTERCHANGEABLE BETWEEN LOTS**

Content: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

**SUBS TMB** SUBSTRATE (PF93619)

1x12 mL (REF 91029)

1x120 mL (REF 91183)

**INTERCHANGEABLE BETWEEN LOTS**

Content: Tetramethylbenzidine 0.26 mg/mL and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8). Liquid, ready for use.

**H<sub>2</sub>SO<sub>4</sub> 0.3 M** STOP SOLUTION (PF93602)

1x16 mL (REF 91029)

1x120 mL (REF 91183)

**INTERCHANGEABLE BETWEEN LOTS**

Content: H<sub>2</sub>SO<sub>4</sub> 0.3 mol/L. Liquid, ready for use.

ADHESIVE FILMS (2).

POLYETHYLENE BAG (1).

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Incubator at 37-40 °C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD ≥ 2.000)
- Microplate washer (optional) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl of solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue

**6. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2-8°C.

The expiry date is printed on each component and on the box label.

**Reagents have a limited stability after opening and/or preparation:**

|                      |                                                                      |
|----------------------|----------------------------------------------------------------------|
| MICROPLATE           | 8 weeks at 2-8°C in polyethylene bag                                 |
| CUT-OFF CONTROL      | 8 weeks at 2-8°C                                                     |
| POSITIVE CONTROL     | 8 weeks at 2-8°C                                                     |
| IgG NEGATIVE CONTROL | 8 weeks at 2-8°C                                                     |
| CONJUGATE            | 8 weeks at 2-8°C                                                     |
| SUBSTRATE            | up to the expiry date at 2-8°C, 1 week at 15-30°C; store in the dark |
| DILUENT 3            | up to the expiry date at 2-8°C                                       |
| WASH BUFFER 10x      | 2 weeks 2-8°C; 5 days 15-30°C                                        |
| STOP SOLUTION        | up to the expiry date at 2-8°C                                       |

**7. SPECIMEN COLLECTION AND STORAGE**

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Avoid the use self-defrosting freezers for the storage of the samples. Defrosted samples must be carefully shaken before the test. Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, or contaminated samples cannot be used.

The test cannot be applied to plasma.

**8. ASSAY PROCEDURE****PREPARATION**

Bring all the reagents and samples to room temperature (18-30°C) before use.

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H<sub>2</sub>O).
- Dilute samples 1:51 distributing 10 µL of serum into 500 µL of diluent.

1. Distribution of the samples: Dispense 100 µL of diluted sample per well (duplicate testing is recommended) or 100 µL of control (positive, negative or cut-off). The minimum requisite is 1 negative control, 1 positive control and 2 cut-off controls. Leave one well for the blank (100 µL of substrate).

2. Incubation: Incubate the plate covered with an adhesive protective film for 45 minutes at 37°C.

3. Washing: Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.

4. Distribution of the conjugate: Dispense 100 µL of conjugate in each well.

5. Conjugate Incubation: Incubate the plate covered with an adhesive protective film for 45 minutes at 37°C.

6. Washing: Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.
7. Distribution of the substrate: Dispense 100 µL of the HRP substrate in each well.
8. Substrate incubation: Incubate the plate for 15 minutes at room temperature.
9. Interruption of the reaction: Dispense 100 µL of blocking solution, in the same order as that followed for point 4.
10. Reading: Read the O.D. at 450 nm or 450/620 within 30 minutes. Repeat the reading at 405 nm in case of O.D. > 2.000.

## 9. SCHEME OF TEST PROCEDURE FOR HERPES SIMPLEX 1 IgG RECOMBINANT

- STEP 1 Place 100 µL of diluted sample/controls (positive, negative and cut-off) in the wells of the strips. Mix well.
- 
- Incubate for 45 min. at 37°C
- 
- Wash 4 times (300 µL)
- 
- STEP 2 Add 100 µL of conjugate to each well
- 
- Incubate for 45 min. at 37°C
- 
- Wash 4 times (300 µL)
- 
- STEP 3 Add 100 µL of Substrate to each well
- 
- Incubate for 15 min. at R.T.
- 
- STEP 4 Add 100 µL of Stop Solution
- 
- Read the O.D. at 450 nm or 450/620 nm within 30 min

## 10. TEST VALIDATION

Negative control: the ratio between the negative control OD and cut-off OD must be < 0.6.

Positive control: the positive control must have an OD at least 1.5 times higher than the cut-off

Cut-off control: the OD of the Cut-off must be > 0.2

Repeat the test, if one of the results of the control sera is not within the acceptability range. If the problem persists contact the Scientific Support.

Tel: 0039 0577 319554  
Fax: 0039 0577 366605  
email: [scientificsupport@diesse.it](mailto:scientificsupport@diesse.it)

## 11. INTERPRETATION OF THE RESULTS

Calculate the ratio between the OD value of the sample and that of the Cut-off (Index).

The sample is considered:

POSITIVE: when the Index is > 1.2

DOUBTFUL: for all the values between 0.8 and 1.2

NEGATIVE: when the Index is < 0.8

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

## 12. LIMITATIONS

A serum sample obtained during the acute phase of infection, when only IgM-class antibodies are present, may be negative by this procedure. The antibodies anti- Herpes Simplex Virus IgM level should be determined using Enzy-Well Herpes Simplex Virus IgM kits. Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

All the positive results require a careful interpretation.

The test cannot be used as only method for a clinical diagnosis. Negative results may not exclude an eventual infection.

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures.

## 13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 0.2 and 0.8 Index.

## 14. ANALYTICAL SPECIFICITY

5 samples (2 negative, 1 Cut-Off and 2 Positive) containing the following interfering substances were tested:

Rheumatoid Factor (44 – 220 UI/ml)  
Bilirubin (4.5 mg/dl – 45 mg/dl)  
Triglycerides (10 mg/dl – 250 mg/dl)  
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test, except for bilirubin that may affect the result.

## 15. CROSS-REACTIONS

205 samples positive to Rubella, Epstein-Barr Virus, Varicella, Cytomegalovirus, HSV 1 and Toxoplasma were tested.

No significant cross-reactions were found.

## 16. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation, 290 samples were analyzed with the Diesse kit as well as with another commercial method. Below are the schematized results of the trial:

|        |       | Reference |     |       |
|--------|-------|-----------|-----|-------|
|        |       | +         | -   | Total |
| Diesse | +     | 45        | 6   | 51    |
|        | -     | 4         | 235 | 239   |
|        | Total | 49        | 241 | 290   |

Percent Positive Agreement (~Diagnostic Sensitivity):

91.8% CI<sub>95%</sub>: 80.8- 96.7

Percent Negative Agreement: (~Diagnostic Specificity):

97.5% CI<sub>95%</sub>: 94.7- 98.8

## 17. PRECISION

| Sample | Within-run Precision |      | Between-run precision |      | Precision between batches |      |
|--------|----------------------|------|-----------------------|------|---------------------------|------|
|        | Mean (Index)         | CV%  | Mean (Index)          | CV%  | Mean (Index)              | CV%  |
| 1      | 0.4                  | 8.6  | 0.5                   | 10.0 | 0.4                       | 15.0 |
| 2      | 0.5                  | 7.5  | 0.4                   | 10.0 | 0.3                       | -    |
| 3      | 0.7                  | 7.5  | 0.4                   | 12.5 | 0.4                       | 15.0 |
| 4      | 0.9                  | 9.4  | 0.9                   | 14.4 | 1.1                       | 5.5  |
| 5      | 1.3                  | 14.4 | 1.4                   | 8.6  | 1.3                       | 9.2  |
| 6      | 2.9                  | 9.5  | 2.1                   | 4.8  | 2.1                       | 4.8  |
| 7      | 5.3                  | 7.0  | 5.0                   | 10.8 | 5.1                       | 7.8  |
| 8      | 7.8                  | 8.0  | 6.8                   | 10.0 | 7.0                       | 5.4  |
| 9      | 11.3                 | 5.4  | 9.7                   | 13.4 | 9.5                       | 4.4  |
| 10     | 13.3                 | 3.3  | 12.5                  | 14.5 | 12.3                      | 2.9  |

## 18. TROUBLESHOOTING

| PROBLEM                      | POSSIBLE ERROR                                            | TEST OR ACTION                                                                               |
|------------------------------|-----------------------------------------------------------|----------------------------------------------------------------------------------------------|
| Invalid run (all negative)   | One or more reagents not added or added in wrong sequence | Recheck procedure. Check for unused solutions. Repeat test.                                  |
|                              | Unreactive plate                                          | Check the code on the package containing the plate (see package insert for correct code).    |
|                              |                                                           | Check for moisture in unused plate. (Silica gel desiccant must be pale yellow). Repeat test. |
| Invalid run (all positive)   | Contamination of substrate                                | Take new aliquot of substrate.                                                               |
|                              | Inadequate washing                                        | Ensure that plate washer works well                                                          |
| Poor precision               | Inadequate aspiration of wells                            | Ensure that plate washer works well                                                          |
|                              | Pipetting error                                           | Check pipette function                                                                       |
|                              | Addition of reagents too slow                             | Avoid drying of the plate after washing step. Add reagents immediately                       |
|                              | Presence of air-bubbles                                   | Avoid air-bubbles formation during pipetting.                                                |
|                              | Optical pathway not clean                                 | Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.  |
| Inadequate color development | Incorrect incubation times or temperature                 | Check for temperature control and time monitoring                                            |
|                              | Inadequate volume of substrate added                      | Check pipette function                                                                       |

## 19. REFERENCES

1. Arvaja M., Lehtinen M., Koskela P., Lappalainen M., Paavonen J., Vesikari T.; "Serological evaluation of Herpes Simplex Virus type 1 and type 2 infections in pregnancy"; Sexual Transmission Infections , 1999; 75: 168-171
2. Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
3. Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L; " Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

## 20. GLOSSARY OF LABELING SYMBOLS

|                                                                                   |                                         |
|-----------------------------------------------------------------------------------|-----------------------------------------|
|  | Date of manufacture                     |
|  | Use By                                  |
|  | Caution, consult accompanying documents |
|  | Manufacturer                            |
|  | Contains sufficient for <n> tests       |
|  | Temperature limitation                  |
|  | Consult instructions for use            |
|  | Biological risks                        |
|  | Catalogue number                        |
|  | In vitro diagnostic medical device      |
|  | Batch code                              |

Manufactured by  
**DIESSE Diagnostica Senese S.p.A.**  
 Strada dei Laghi 39  
 53035 Monteriggioni (Siena)  
 Italy