

Chorus SARS-CoV-2

IgA, IgG and IgM

For the serological diagnosis of Covid-19 infection




The new CHORUS serological test kits for diagnosing the SARS-CoV-2 infection were developed in collaboration with the National Institute for Infectious Disease Lazzaro Spallanzani in Rome, where the virus circulating in Italy was first isolated.

The antigen used in the solid phase of the kit is prepared, according to a patented protocol, in our BSL3 high biological containment laboratories. The native virus is first deactivated, to ensure complete safety.

Our kits detect the various serum antibody isotypes: IgA, IgG and IgM produced in the immune response to the virus. In general, the detection of IgM - antibodies produced in the early post-infection phases - is a marker of the acute phase, while IgG are produced in a later phase. If a patient has been in contact with other Coronaviruses in the past, however, the IgG titre may rise more quickly than that of the IgM and IgA.

Detecting IgA is particularly important in SARS-CoV-2 infections, as these antibodies are produced in the respiratory mucosa. The purpose of these antibodies is to bind to the virus and prevent its entry into the body, and are therefore an effective marker of the acute phase.


Indirect ELISA



SARS-CoV - 2 is a "new" virus, the use of the "whole" virion allows us to determine antibodies directed against all viral proteins, even in their conformational state.

Molecular tests are most sensitive in the early stages of SARS-CoV-2 infection. The positivity rate, higher than 90% in the first three days of the onset of symptoms, declines quickly; less than 80% on the sixth day and well under 50% after 14 days.

The sensitivity of molecular tests is higher than serological tests if performed within 5.5 days of onset of symptoms, while the determination of IgM and IgA is more sensitive in later stages. In asymptomatic patients, the molecular test has a sensitivity of 52%, if performed alone but determining IgM and IgA simultaneously sensitivity levels rise to more than 96%



Performing serological tests and molecular tests contemporarily is particularly useful in the case of asymptomatic patients, allowing the immune response against the virus to be highlighted, and thus estimation of the real incidence rate of the infection in the population.

Graphic elaboration of the evolution of antibody titers in relation to positivity to molecular tests (see bibliography).

Chorus SARS-CoV-2 serology kits highlights

Single-test ready-to-use diagnostic devices

Convenient packaging: 6 devices per 6 resealable aluminum triplex pouches

Reduced analysis times: results in just 37 minutes

Flexibility and simplicity in the execution of the tests

Complete traceability of results thanks to the Identisystem technology

The new CHORUS kits allow the detection of the various antibody classes produced as a result of the immune response to the SARS-CoV-2 infection:

IgA especially important in respiratory viruses, because they are produced in respiratory mucosa, the entrance door of the virus;

IgG produced in the later stages of infection, a good marker of prior infection;

IgM, acute phase marker, the first antibody produced by the immune system in response to an infection.

Available kits:

Chorus SARS-CoV-2 IgG

81400

36 tests

For the determination of IgG antibodies against SARS-CoV-2 in human serum.
Contains calibrator and positive control.

Chorus SARS-CoV-2 IgM

81401

36 tests

For the determination of IgM antibodies against SARS-CoV-2 in human serum.
Contains calibrator and positive control.

Chorus SARS-CoV-2 IgA

81402

36 tests

For the determination of IgA antibodies against SARS-CoV-2 in human serum.
Contains calibrator and positive control.

Bibliography:

Zhang Wei et Al. Molecular and serological investigation of 2019-nCoV infected patients: implication of multiple shedding routes. Emerging Microbes & Infection. 2020; 9 (1). Ahead of printing.

Ng Lisa & Hiscox JA. Coronaviruses in animals and humans. The British Medical Journal 2020. Ahead of printing.

Guo L. et al. Profiling Early Humoral Response to Diagnose Novel Coronavirus Disease (COVID-19). Clin Infect Dis. 2020, 21 mar. Ahead of printing

BèNè MC et Al. Good IgA bad IgG in SARS-CoV-2 infection? Clin Infect Dis. 2020, 11 Apr. Ahead of printing.



www.diesse.it