

**DIESSE DIAGNOSTICA SENESE SpA**  
 53035 Monteriggioni (SI) – Italia – Via delle Rose 10  
 Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione  
 20144 Milano (Italy) Via A. Solari 19 sc. 6  
 Tel. +39-02-4859121 r.a. Fax +39-02-48008530

**CERTIFICATO DI CONTROLLO QUALITÀ/QUALITY CONTROL CERTIFICATE**

**PRODOTTO/PRODUCT** ESR CONTROL CUBE

**REF:** 10435 (4x9 ml)  
 10436 (2x9 ml)

**REAGENTI/REAGENTS**

	Kit	Level 1	Level 2
Lot:	443	443	443
	2018-10	2018-10	2018-10
	2020-04	2020-04	2020-04

**RISULTATI/RESULTS**

	Level 1		Level 2	
	Risultato/ Result	Range	Risultato/ Result	Range
VES-Matic Cube 30	7	1-15	74	49-99
VES-Matic Cube 80/VES-Matic Cube 200	6	1-14	55	30-80
MINI-CUBE/Cube 30 touch	9	1-17	62	37-87

**Data/Date** 15/11/2018

**Testato da/Tested by:** Katia Tanzini

**Approvato da/Approved by:** Alessandra Cianchi (Responsabile CQ/QC Head)