

DIESE DIAGNOSTICA SENESE SpA
 Via delle Rose 10, 53035 Monteriggioni (SI), Italy
 Tel. ++39-0577-587111 fax ++39-0577-318690

Sede Legale e Amministrazione
 Via A. Solari 19 20144 Milano Italy
 Tel. ++39-02-4859121 r.a. fax ++39-02-48008530

CERTIFICATE OF ANALYSIS

CODE	PRODUCT NAME	LOT	EXPIRY DATE	PREP DATE
10439	SediCheck Cube 30 4 x 3 ml	018	2019-06	2018-06

VES-MATIC CUBE 30

Level I (NORMAL Level) lot 424	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	7	1-17

Level II (ABNORMAL Level) lot 424	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	69	47-91

MINI-CUBE

Level I (NORMAL Level) lot 424	Test conditions	VES 1st hour	Range of acceptability
MINI- CUBE	Temperature ON	4	1-17

Level II (ABNORMAL Level) lot 424	Test conditions	VES 1st hour	Range of acceptability
MINI- CUBE	Temperature ON	43	26-75

Date: 13/07/2018


 Signature (Head Quality Control)