

**DIESE DIAGNOSTICA SENESE SpA**  
 53035 Monteriggioni (SI) - Italia - Via delle Rose 10  
 Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione  
 20144 Milano (Italy) Via A. Solari 19 sc. 6  
 Tel. +39-02-4859121 r.a. Fax +39-02-48008530

**QUALITY CONTROL CERTIFICATE**

**KIT NAME: FT3 - CHORUS**

Code	86518		
Lot No.	500	2020-01	2019-01

**REAGENTS**

Label	Description	Lot.	Exp	No
D3300	FT3	500	2020-01	6X6
B3300	CALIBRATOR FT3	500	2020-01	1
C3300	CONTROL SERUM FT3	500	2020-01	1


**CONTROLS**

	Found results	Expected results
CALIBRATOR (mOD)	320	≥100
CONTROL SERUM (pmol/L)	8,7	6,1-15,5

**CALIBRATOR**

9	6	7	5	0	0	6	3	2	8	3	9	5	0	5	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader **YES**  Operator Signature *Chiara*



**CONTROL SERUM**

9	6	7	5	0	0	6	1	4	2	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES**  Operator Signature *Chiara*



The test has been carried out in accordance with the package insert

**TESTED BY:** Chiara Bruschelli

**APPROVED** RESPONSIBLE HEAD, Quality Control Reagents..... *Chiara*

**DATE** ..... 22/01/2019

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.