



VESTMMATIC EASY



Automatic
instrument
for the
determination
of the **ESR**
on 10 samples
simultaneously
or in
random access
mode.

**NIECCE
VIESSE**

VES-MATIC EASY:

*Complete with printer.
USB connector.
Ref# 10376*



technical characteristics

- Up to 10 samples using Ves-tec or Vacu-tec tubes
- Random access mode
- Printed results
- Results in 20 minutes
- External barcode reader connector
- LCD screen of 16 characters per line
- USB port for host connection
- Dimension in centimeters (15 x 15 x 22) (l x w x h)
- External power supply IN100/240 VAC OUT 9VDC
- External barcode reader (optional)

DIESSE Diagnostica Senese SpA
Via del Pozzo, 5 - Loc. S. Martino • 53035 Monteriggioni (SI) Italy
Tel.: ++39/0577/31.95.60/61/50 • Fax: ++39/0577/31.87.63