



QUALITY CONTROL CERTIFICATE

KIT NAME: DIPHTHERIA IgG - CHORUS

Code 81266/12				
Lot No.	626	2021-01	2020-04	

REAGENTS

Label	Description	Lot.	Exp	No
D9670	DIPHTHERIA IgG	626	2021-01	2X6
B9670	CALIBRATORE DIPHTHERIA IgG	626	2021-04	1
C9670	CONTROL DIPHTHERIA IgG	626	2021-04	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	360	≥ 100
Control (mIU/ml)	470,3	304,3-617,4

CALIBRATOR

9	4	0	6	2	6	5	2	3	3	5	1	5	0	6	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *SB*

940626623361606400  
Diphtheria IgG  
CALIBRATOR 0.175 mL IVD  
LOT 626  
2021-04 2°C-8°C  
DIESSE Diagnostica Senese Monteriggioni - Italy CE

CONTROL +

9	4	0	6	2	6	5	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *SB*

940626614260000001  
Diphtheria IgG  
CONTROL+ 0.425 mL IVD  
LOT 626  
2021-04 2°C-8°C  
DIESSE Diagnostica Senese Monteriggioni - Italy CE

The test has been carried out in accordance with the package insert

TESTED BY: Sara Benocci *SB*  
APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Plavet*  
DATE *05-05-20*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.