



**QUALITY CONTROL CERTIFICATE**

**KIT NAME: SARS-CoV-2 IgM - CHORUS**

<b>Code</b>	81401		
<b>Lot No.</b>	667	2021-08	2020-09

**REAGENTS**

Label	Description	Lot.	Exp	No
D4190	SARS-CoV-2 IgM	667	2021-08	6X6
B4190	CALIBRATORE SARS-CoV-2 IgM	667	2021-09	1
C4190	CONTROL SERUM SARS-CoV-2 IgM	667	2021-09	1

**CONTROLS**

	Found results	Expected results
CALIBRATOR (mOD)	526	≥100
CONTROL SERUM (Index)	2,0	1,4-2,6

**CALIBRATOR**



1	6	4	6	6	7	6	0	0	0	0	5	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES**  Operator Signature *[Signature]*

**CONTROL SERUM**

1	6	4	6	6	7	6	1	4	2	6	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES**  Operator Signature *[Signature]*



The test has been carried out in accordance with the package insert

TESTED BY: Chiara Bruschelli *[Signature]*  
APPROVED RESPONSIBLE HEAD, Quality Control Reagents... *[Signature]*  
DATE 12.11.2020

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.