

Product	REF	CPU 2003	CPU 2010			Method version
		SOFTWARE 2	SOFTWARE 3	SOFTWARE 4	SOFTWARE 5 (Only for Autosampler)	
Chorus Borrelia IgG	81046	NOT APPLICABLE	Rel. 3.2.1	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	from 1030
Chorus Borrelia IgG	81046/C03	Rel. 2.10 rev 1	NOT APPLICABLE	NOT APPLICABLE	NOT APPLICABLE	from 1030
Chorus Borrelia IgM	81047	NOT APPLICABLE	Rel. 3.2.1	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	from 1030
Chorus Borrelia IgM	81047/C03	Rel. 2.10 rev 1	NOT APPLICABLE	NOT APPLICABLE	NOT APPLICABLE	from 1030
Chorus Echinococcus IgG	81284	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Dengue IgG	81154	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Dengue IgM	81156	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Tetanus IgG	81270	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Diphtheria IgG	81266	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Tick borne encephalitis IgM	81277	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Toxocara IgG	81280	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Calprotectin K	86601	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Helicobacter pylori Ag	81063	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Clostridium difficile GDH	81168	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Clostridium difficile A/B Toxins	81170	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Tg High sensitivity	86104	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus FSH	86510	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus LH	86512	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus PRL	86514	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus TSH	86516	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus fT3	86518	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus fT4	86520	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Tick borne encephalitis IgG (up to lot 459)	81276	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038/1039
Chorus Tick borne encephalitis IgG (from lot 469)	81276	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Bordetella Pertussis Toxin IgG	81272	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Bordetella Pertussis Toxin IgA	81274	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus 25OH Vitamin D Total (up to lot 499)	86900	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	from 1034 to 1039
Chorus 25OH Vitamin D Total (from lots subsequent to 505)	86900	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Q Fever Phase II IgG	81171	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Q Fever Phase II IgM	81172	Rel. 2.10 rev 1	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus PSA	86550	NOT APPLICABLE	Rel. 3.2.9 [25]	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1039
Chorus Toxoplasma IgM Capture	81042	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	2040
Chorus Toxoplasma IgA Capture	81044	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	2040
Chorus Cytomegalovirus IgM Capture	81012	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	2040
Chorus CH-50	86102	NOT APPLICABLE	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	NOT APPLICABLE	1040-2040
Chorus Brucella IgG	81164	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1040-2040
Chorus Brucella IgM	81166	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1040-2040
Chorus Allergy	88000 88000-xx	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [22]	Rel. 5.0.0 [4]	1038
Chorus Sars-CoV-2 IgG	81400	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1041-2041
Chorus Sars-CoV-2 IgM	81401	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1041-2041
Chorus Sars-CoV-2 IgA	81402	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1041-2041
Chorus Sars-CoV-2 Screen Serum	81406	Rel. 2.10 rev 1	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1042-2042
Chorus TGE	86110	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1041-2041
Chorus TRAb	86112	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1041-2040
Chorus Phospholipid Screen-G	86106	NOT APPLICABLE	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1040-2040
Chorus Phospholipid Screen-M	86108	NOT APPLICABLE	Rel. 3.2.9 [65]	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1040-2040
Chorus Sars-CoV-2 Antigen	81410	NOT APPLICABLE	NOT APPLICABLE	Rel. 4.0.6 [65]	Rel. 5.0.0 [4]	1043-2043