

Chorus ANA-8 (REF: 86010 and REF: 86010/12) is an immunoenzymatic method for the qualitative determination of IgG class antibodies against 8 cellular and nuclear antigens (Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB) in human serum, using a disposable device applied on the chorus and Chorus TRIO instruments.

REF 86010: each package contains 6 bags with 6 devices each (sufficient for the total of 36 determinations).

REF 86010/12: each package contains 2 bags with 6 devices each (sufficient for the total of 12 determinations).

Each kit consists of:

### DEVICE

**Position 8:** Space for application of bar code label

**Position 7:** Empty

**Position 6:** MICROPLATE WELL

Coated with recombinant and/or highly purified antigens

**Position 5:** MICROPLATE WELL

Not coated

**Position 4:** TMB SUBSTRATE 0.400 ml

The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet.

Contents: Tetramethylbenzidine 0.26 mg/ml and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

**Position 3:** SAMPLE DILUENT 0.400 ml

The substance is classified as **dangerous** according to Regulation 1272/2008.

For this mixture the material safety data sheet, prepared in accordance with EC Regulation 1907/2006, is available below.

Contents: Saline proteic solution with Proclin (0.1%)

**Position 2:** CONJUGATE 0.350 ml

The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet.

Contents: Anti-human IgG monoclonal antibodies labeled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

**Position 1:** EMPTY WELL

In which the sample is transferred

## CALIBRATOR

### CALIBRATOR

0.175 ml

Contents: Diluted human serum containing IgG antibodies anti- Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB and preservative. Liquid, ready for use.

## POSITIVE CONTROL

### CONTROL +

0.425 ml

Contents: Diluted human serum containing IgG antibodies anti- Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB and preservative. Liquid, ready for use.

Both mixtures are classified as **dangerous** according to Regulation 1272/2008/EC.

For these mixtures the material safety data sheet, prepared in accordance with EC Regulation 1907/2006, is available below.



# MATERIAL SAFETY DATA SHEET

## SAMPLE DILUENT – DT20 CALIBRATOR/POSITIVE CONTROL

(According to Regulation (EC) No. 1907/2006)

### 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

#### 1.1 PRODUCT IDENTIFIERS

Product name: 1. SAMPLE DILUENT – DT20  
2. CALIBRATOR/POSITIVE CONTROL

Product code: 1. PF30212  
2. PF86810-C1 (CALIBRATOR); PF86810-C2 (POSITIVE CONTROL)

Brand: DIESSE

#### 1.2 RELEVANT IDENTIFIED USES OF THE SUBSTANCES OR MIXTURE AND USES ADVISED AGAINST

Identified use: Professional use as laboratory reagent.  
DT20: Proteic solution used for dilution of samples.  
Calibrator: diluted human serum, necessary for the calibration of the instrument.  
Positive Control: diluted human serum, used to verify the correctness of the obtained results and for the subsequent validation of the test.

#### 1.3 DETAILS OF THE SUPPLIER OF THE SAFETY DATA SHEET

Company: DIESSE Diagnostica Senese S.p.A

|  |  |
|--|--|
| <u>Registered office:</u><br>Via A. Solari 19<br>20144 Milan, Italy<br>Tel: +39 02 4859121<br>Fax: +39 02 48008530 | <u>Production plant:</u><br>Strada dei Laghi, 39<br>53035 Monteriggioni (SIENA), Italy<br>Tel: +39 0577 307100<br>e-mail: contatti@diesse.it |
|--|--|

#### 1.4 EMERGENCY TELEPHONE NUMBER

Emergency number: Centro Antiveleni, Ospedale Niguarda Ca' Granda - Milan  
Tel: +39 02 66101029

Centro Antiveleni, Azienda Ospedaliera "S.G.Battista" – Molinette di Torino - Turin  
Tel: +39 011 6637637

Centro Antiveleni – U.O. Tossicologia Medica, Azienda Ospedaliera Careggi – Florence  
Tel: +39 055 4277238

Centro Antiveleni, Policlinico A. Gemelli – Università Cattolica del Sacro Cuore - Rome  
Tel: +39 06 3054343

Centro Antiveleni, Azienda Ospedaliera A. Cardarelli – Naples  
Tel: +39 081 7472870

## 2. HAZARDS IDENTIFICATION

### 2.1 CLASSIFICATION OF THE SUBSTANCES OR MIXTURE

Classification according to Regulation (EC) No. 1272/2008: Skin Sensitizer– Category 1

Hazard statement: H317 – May cause an allergic skin reaction.

### 2.2 LABEL ELEMENTS

Pictograms: GHS07

Signal word: Warning

Hazard statement(s):  
H317 – May cause an allergic skin reaction.

Precautionary statement(s):

*Prevention:*

P261 – Avoid breathing dust/fume/gas/mist/vapors/spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

*Response:*

P302+P352 – IF ON SKIN: Wash with plenty of soap and water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P363 – Wash contaminated clothing before reuse.

*Disposal:*

P501 – Dispose of contents/container in accordance with local regulation

#### Contains:

Reaction mass of: 5-chloro-2-methyl-4-iso-thiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1);

Index No: 613-167-00-5



### 2.3 OTHER HAZARDS

None

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substances Not applicable

3.2 Mixtures

| International Chemical Identification  | Concentration | Classification<br>Regulation EC/1272/2008 |      |
|--|---------------|---|------|
| Reaction mass of: 5-chloro-2-methyl-4-iso-thiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1);<br>Cas No 55965-84-9<br>Index No 613-167-00-5 | 0.0015-0.06%  | Acute Tox. 3                              | H331 |
|  |               | Acute Tox. 3                              | H311 |
|  |               | Acute Tox. 3                              | H301 |
|  |               | Skin Corr. 1B                             | H314 |
|  |               | Skin Sens. 1                              | H317 |
|  |               | Aquatic Acute 1                           | H400 |
|  |               | Aquatic Chronic 1                         | H410 |

The entire text of Hazard Statements is reported at Section 16 of the sheet.

## 4. FIRST AID MEASURES

### 4.1 DESCRIPTION OF FIRST AID MEASURES

|                 |   |
|-----------------|---|
| General advice: | Move the person from the exposure to open air. In case of needs consult a physician immediately and show this Material Safety Datasheet. Eyewash and shower for the treatment of emergency have to be present in the workplace. |
| Skin contact:   | Wear off the contaminated clothes and wash with copious amounts of water (for at least 15 minutes). If irritation persists consult a physician.   |
| Eye contact:    | If present, remove immediately contact lenses. Wash with plenty of water (for at least 15 minutes), keeping eyelids opened. Consult an oculist if the irritation persists.  |
| Inhalation:     | Move the person from the exposure to open air. If irritation occurs consult a physician.  |
| Ingestion:      | Rinse mouth immediately and drink a copious amount of water. Call a physician immediately. Do not cause vomiting.   |

### 4.2 MOST IMPORTANT SYMPTOMS AND EFFECTS

No data available

### 4.3 INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

No data available, other than point 4.1

## 5. FIREFIGHTING MEASURES

### 5.1 EXTINGUISH MEDIA

|                                   |   |
|-----------------------------------|---|
| Suitable extinguishing media:     | Use extinguishing measures (CO <sub>2</sub> , foam, dry powder, water) that are appropriate to local circumstances and the surrounding environment. |
| Not Suitable extinguishing media: | None  |

### 5.2 SPECIAL HAZARDS ARISING FROM THE SUBSTANCE OR MIXTURE

The mixture is not flammable, however in case of fire may release irritating and/or toxic gases.

### 5.3 ADVICE FOR FIREFIGHTERS

Wear appropriate personal protective equipment and clothing. In case of fire, wear self-contained breathing apparatus and avoid that fire extinguishing water contaminates surface water and/or groundwater.

## 6. ACCIDENTAL RELEASE MEASURES

### 6.1 PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES

Avoid the contact with skin and eyes and evacuate the area, keeping people not involved in the intervention operations away. Ensure an adequate ventilation of the affected area.

Do not handle damaged containers or the leaked product before wearing appropriate protective outfit.

## 6.2 ENVIRONMENTAL PRECAUTIONS

Avoid the contamination of surface water, soil and the dispersion in the air. Do not let product enter into drains. Discharge into the environment must be avoided.

## 6.3 METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP

Contain the leakages with earth and sand. Collect the spilled material and store it in suitable containers for disposal. Use water only to remove residuals, to avoid the danger of spill of product into the sewers.

## 6.4 REFERENCE TO OTHER SECTIONS

For further information see section 8 and 13.

# 7. HANDLING AND STORAGE

## 7.1 PRECAUTION FOR SAFE HANDLING

Work in well ventilated areas and in the presence of ventilation systems or personal protective equipment. Do not inhale vapors or mists. Avoid the contact with eyes, skin and clothes. Limit repeated exposure.

## 7.2 CONDITION FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES

Keep the containers at 2-8°C in a cool, well-ventilated area, away from heat sources and humidity.

## 7.3 SPECIFIC END USE

None

# 8. EXPOSURE CONTROL/PERSONAL PROTECTION

## 8.1 CONTROL PARAMETERS

Applicable exposure limits:           None

## 8.2 EXPOSURE CONTROLS

Work and handle according to the usual precautionary measures for handling chemicals.  
Do not eat, drink or smoke while handling the product; wash hands thoroughly with soap and water before meals and after the work shift. Immediately remove all contaminated clothing.

Appropriate engineering  
controls:                           Ensure an adequate ventilation of the working area.

Personal Protective  
Equipment:                   The suggestions on the use of specific PPE are indicative. Their choice should be made according to the use of the product and the instructions given by the supplier of the device.

Hand protections:  
Chemical-resistant gloves, compliant with EN 374

Eye protections:  
Side Shields Safety Goggles compliant with EN 166

Body protections:

Work outfits

Respiratory protections:

Not required under normal work activities

### 8.3 ENVIRONMENTAL EXPOSURE CONTROLS:

Do not discard residuals in the environment.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 INFORMATION ON BASIC PHYSICAL AND CHEMICAL PROPERTIES

|   |                          |
|---|--------------------------|
| Appearance                                  | Liquid                   |
| Odor:                                       | Odorless                 |
| Odor threshold:                             | Not pertinent            |
| pH:   | See it at the first page |
| Melting point/freezing point                | Not determined           |
| Initial boiling point and boiling range     | Not determined           |
| Flash point                                 | Not applicable           |
| Evaporation rate                            | Not determined           |
| Flammability                                | Not flammable            |
| Upper/lower flammability or explosive limit | Not applicable           |
| Vapor pressure                              | Not determined           |
| Vapor density                               | Not determined           |
| Relative density                            | Not determined           |
| Water solubility                            | Soluble in water         |
| Partition coefficient: n-octano/water       | Not determined           |
| Autoignition temperature                    | Not applicable           |
| Decomposition temperature                   | Not determined           |
| Viscosity                                   | Not determined           |
| Explosive properties                        | Not explosive            |
| Oxidizing properties                        | Not determined           |

### 9.2 OTHER SAFETY INFORMATION

None

## 10. STABILITY AND REACTIVITY

### 10.1 REACTIVITY

In the normal condition, no risk of reactivity is present.

### 10.2 CHEMICAL STABILITY

Stable until the expire date under the recommended transport, handling and storage conditions

### 10.3 POSSIBILITY OF HAZARDOUS REACTIONS

The mixture could react with the basic and/or oxidizing substances

#### 10.4 CONDITIONS TO AVOID

Avoid the storage at temperature different from that are advised

#### 10.5 INCOMPATIBLE MATERIALS

Strong oxidizing and reducing agents

#### 10.6 HAZARDOUS DECOMPOSITION PRODUCTS

When heated or in case of fire, vapors potentially dangerous to health may be produced.

### 11. TOXICOLOGICAL INFORMATION

#### 11.1 INFORMATION ON TOXICOLOGICAL EFFECTS

There is no known information on toxicological effects from exposure to the mixtures. However, it's recommended to operate according to the good operational practices. Nevertheless, the mixture can cause mild health effects on sensitive people by inhalation, cutaneous absorption, contact eyes or on ingestion.

|  |                   |
|--|-------------------|
| Acute toxicity   | No data available |
| Skin irritation/corrosion                                  | No data available |
| Eye lesions/severe eye irritation                          | No data available |
| Respiratory/skin sensitization                             | No data available |
| Mutagenicity for germ cells                                | No data available |
| Carcinogenicity  | No data available |
| Reproductive toxicity                                      | No data available |
| Toxicity for target organs (single and repeated exposures) | No data available |
| Hazards in case of inhalation                              | No data available |
| Additional information                                     | No data available |

### 12. ECOLOGICAL INFORMATION

#### 12.1 TOXICITY

Use according to the good working practices, avoiding the disposal in the environment. In case the product reach waterways or sewers or contaminate soil or vegetation, inform the competent authorities.

#### 12.2 PERSISTENCE AND DEGRADABILITY

No data available

#### 12.3 BIOACCUMULATIVE POTENTIAL

No data available

#### 12.4 MOBILITY IN SOIL

No data available

#### 12.5 RESULTS OF PBT AND VPVB ASSESSMENT

No data available



## 12.6 OTHER ADVERSE EFFECTS

No data available

## 13. DISPOSAL CONSIDERATIONS

### 13.1 WASTE TREATMENT METHODS

The samples and all the used reagents have to be handled as potentially infected.  
The product and its containers should be considered special waste.  
Their transport and disposal should be performed by authorized specialized companies according to the law.

## 14. TRANSPORT INFORMATION

Not hazardous good according to the transport regulations.

## 15. REGULATORY INFORMATION

### 15.1 SAFETY, HEALTH AND ENVIRONMENTAL REGULATIONS/LEGISLATION SPECIFIC FOR THE SUBSTANCE OR MIXTURE

Hazard classification, labeling and packaging under Regulation 1272/2008 (EC) and its subsequent amendments.

#### Legislative Decree 81/2008:

The use of this chemical entails the obligation of "Risk Assessment" by the employer in accordance with the provisions of the Decree April 9, 2008 n. 81 and subsequent amendments. Workers exposed to this chemical agent must not undergo health checks whether the results of the risk assessment show that, depending on the type and quantity of chemical agent and the method and frequency of exposure to this agent, there is only one "Low risk to the safety and irrelevant to the health" of the workers and the measures provided for in the same Decree are sufficient to reduce the risk

#### Directive 96/82/EC (Seveso Directive):

Not applicable

### 15.2 CHEMICAL SAFETY ASSESSMENT

Not performed for the product

## 16. OTHER INFORMATION

This product has to be used for diagnostic use only by personnel who is qualified and trained on the hazards shown in this safety sheet.

Text of the hazard statements present at point 3

|      |  |
|------|--|
| H301 | Toxic if swallowed                                   |
| H311 | Toxic in contact with skin                           |
| H314 | Causes severe skin burns and eye damage              |
| H317 | May cause an allergic skin reaction                  |
| H331 | Toxic if inhaled                                     |
| H400 | Very toxic to aquatic life                           |
| H410 | Very toxic to aquatic life with long lasting effects |

#### REFERENCES:

1. Regulation (EC) 1907/2006 of the European Parliament (REACH) as amended
2. Regulation (EC) 1272/2008 of the European Parliament (CLP) as amended

3. ECHA European Chemicals Agency [www.echa.europa.eu](http://www.echa.europa.eu)
4. The Merck Index.
5. Handling Chemical Safety
6. NIOSH - Registry of Toxic Effects of Chemical Substances
7. INRS - Fiche Toxicologique

*DIESSE Diagnostica Senese SpA declares that the information contained in this data sheet is based on the knowledge available to us on the date of the last version. Users must verify the suitability and accuracy of the provided information according to the specific use of the product. The document has not been accepted as a guarantee of any specific product property. Because of the use of this product is not subject to our direct control, the users, under their own responsibility, must observe the laws and regulations relating to health and safety. We accept no responsibility for improper use.*

**This sheet supersedes any previous edition.**