

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI) - Italia - Via delle Rose 10  
Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione  
20144 Milano (Italy) Via A. Solari 19 sc. 6  
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

## QUALITY CONTROL CERTIFICATE

KIT NAME: FT4 - CHORUS

|         |       |         |         |
|---------|-------|---------|---------|
| Code    | 86520 |         |         |
| Lot No. | 844   | 2023-06 | 2022-06 |

## REAGENTS

| Label | Description       | Lot. | Exp     | No  |
|-------|-------------------|------|---------|-----|
| D3320 | FT4               | 844  | 2023-06 | 6X6 |
| B3320 | CALIBRATOR FT4    | 844  | 2023-06 | 1   |
| C3320 | CONTROL SERUM FT4 | 844  | 2023-06 | 1   |

## CONTROLS

| Found results | Expected results |
|---------------|------------------|
|---------------|------------------|


|                        |      |          |
|------------------------|------|----------|
| CALIBRATOR (mOD)       | 310  | ≥100     |
| CONTROL SERUM (pmol/L) | 12,7 | 2,6-29,6 |

## CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 9 | 6 | 8 | 8 | 4 | 4 | 6 | 2 | 3 | 0 | 4 | 1 | 2 | 7 | 4 | 8 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader

YES


Operator Signature *se*

## CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 9 | 6 | 8 | 8 | 4 | 4 | 6 | 2 | 6 | 4 | 8 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader

YES

Operator Signature *se*

The test has been carried out in accordance with the package insert

TESTED BY: Antonella Leo *Antonella Leo*APPROVED RESPONSIBLE HEAD, Quality Control Reagents *se*DATE *21.06.2022*

## WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.