

The kit **Enzy-Well Herpes simplex 2 IgG (REF 91026 - 91185)** is an immunoenzymatic method for the qualitative determination of IgG-class antibodies to Herpes simplex virus (type 2) in human serum.

Each	kit	is	composed	of:
------	-----	----	----------	-----

MT PLATE	MICROPLATE 12 x 8 wells coated with Herpes simplex virus (type 2) (91026) 10 x 96 wells coated with Herpes simplex virus (type 2) (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.
CONTROL +	POSITIVE CONTROL 1 x 1.6 mL (91026) 4 x 1.6 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.
	<u>Contents</u> : diluted human serum, containing anti-HSV 2 IgG antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.
CONTROL CUT-OFF	
	$5 \times 2 \text{ mL}$ (91185) The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive
	1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material
	safety data sheet. <u>Contents</u> : diluted human serum, containing anti-HSV 2 IgG antibodies, in phosphate buffer 0.01 mol/L with BSA
	1% and sodium azide 0.09%.
CONJ	CONJUGATE 1 x 16 mL (91026) 8 x 16 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations
	above the established limits. For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material
	safety data sheet.
	<u>Contents</u> : anti-human IgG monoclonal antibodies labelled with peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.
SORBENT 1	HSV Type 1 Freeze-dried powder x 6 vials (91026) Freeze-dried powder x 40 vials (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive
	<u>1999/45/CE and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.
	<u>Contents</u> : purified and inactivated HSV 1, in phosphate buffer with lactose.
CONTROL IgG -	IgG NEGATIVE CONTROL 1 x 1.6 mL (91026) 4 x 1.6 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and / or with limit values for exposure in
	concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.
	Contents: human serum diluted in phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%.
WASH BUF 10x	WASHING BUFFER 10X 1 x 100 mL (91026) 5 x 100 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive
	<u>1999/45/CE and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.
	<u>Contents</u> : Phosphate buffered saline (PBS), concentrated 10 times; contains Brij 0.5%.



SAFETY INFORMATION ENZY-WELL HERPES SIMPLEX 2 IgG

SAMP DIL 2	DILUENT 2 1 x 100 mL (91026)
	6 x 100 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive
	1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations
	above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material
	safety data sheet.
	Contents: Proteic solution in phosphate buffer with sodium azide 0.09% containing methyl orange.
SUBS TMB	SUBSTRATE 1 x 12 mL (91026)
	1 x 120 mL (91185)
	The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive
	1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations
	above the established limits.
	For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material
	safety data sheet.
	Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05
	mol/L (pH 3.8).
H ₂ SO ₄ 0.3 M	STOP SOLUTION 1 x 16 mL (91026)
	1 x 120 mL (91185)
	The substance is classified as dangerous according to Regulation 1272/2008 and European Directive
	<u>1999/45/CE.</u>
	The material safety data sheet, prepared in accordance with EC regulations 1907/2006 and 453/2010, is
	available below.
	<u>Contents</u> : solution of sulphuric acid, 0.3 mol/L.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

1. IDENTIFICATION OF THE SUBST	ANCE/MIXTURE AND THE COMPANY/	UNDERTAKING
1.1 IDENTIFICATION OF THE PRODU	ст	
Product name:	ENZY-WELL STOP SOLUTION	
Product number:	PF 93602	
Brand:	DIESSE	
1.2 IDENTIFIED USES OF THE SUBST	ANCE OR MIXTURE AND USES ADVISED AGAIN	IST
Identified use:	Professional use as laboratory rea	gent.
1.3 DETAILS OF THE SUPPLIER OF T	HE SAFETY DATA SHEET	
Company:	DIESSE Diagnostica Senese SpA	
	<u>Registered office:</u> Via A. Solari 19 20144 Milan, Italy Tel: +39 02 4859121 Fax: +39 02 48008530	<u>Production plant</u> : Via delle Rose, 10 – Loc. Tognazza 53035 Monteriggioni (SI), Italy Tel: +39 0577 587111 Fax: +39 0577 318690
1.4 EMERGENCY TELEPHONE NUMBE	R	
Emergency number:	Centro Antiveleni, Ospedale Nigua Tel: +39 02 66101029	rda Ca' Granda - Milano
	Centro Antiveleni, Azienda Ospeda Tel: +39 011 6637637	aliera "S.G.Battista" – Molinette di Torino - Torino
	Centro Antiveleni – U.O. Tossicolo Tel: +39 055 4277238	gia Medica, Azienda Ospedaliera Careggi – Firenze

Centro Antiveleni, Policlinico A. Gemelli – Università Cattolica del Sacro Cuore - Roma Tel: +39 06 3054343

Centro Antiveleni, Azienda Ospedaliera A. Cardarelli – Napoli Tel: +39 081 7472870

web: http://www.salute.gov.it/servizio/documenti/centri antiveleni.pdf

2. HAZARDS IDENTIFICATION

2.1 CLASSIFICATION OF THE SUBSTANCES OR MIXTURE

Classification according to Regulation (EC) No. 1272/2008: Skin Corrosion – Category 1A

Dangerous components:

H314 - Causes severe skin burns and eye damage

Classificazione ai sensi delle Direttive 67/548/CEE e 1999/45/CE: C - Corrosive

R-phrases:

R35 - Causes severe skin burns



(According to Regulation (EC) No. 1907/2006 e 453/2010)

2.2 LABEL ELEMENTS

Pictograms: GHS05

Signal word: Danger

Hazard statement:

H314 - Causes severe skin burns and eye damage

Precautionary statement:

Prevention P260 - Do not breathe dust/fume/gas/mist/vapours/spray P264 - Wash ... thoroughly after handling P280 - Wear protective gloves/protective clothing/eye protection/face protection *Reaction* P301 + P330 + P331 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting P303 + P361 + P353 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower P363 - Wash contaminated clothing before reuse P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing P310 - Immediately call a POISON CENTER or doctor/physician P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing *Storage* P405 - Store locked up.

Contains:

Sulphuric Acid INDEX: 016-020-00-8

2.3 OTHER HAZARDS

None

3. COMPOSITION/INFORMATION ON INGREDIENTS

3.1 SUBSTANCE

Non applicable

3.2 MIXTURE

Substance name	Conc.	Directive	Classi 67/548/CEE	fication EC Regulation 1272/2008
Sulphuric Acid Cas No 7664-93-9 EC No 231-639-5 Index No 016-020-00-8	<3%	С	R35	Skin Corr. 1A H314

4. FIRST AID MEASURES

4.1 DESCRIPTION OF FIRST AID MEASURES

General advice:	Move the person from the exposure to open air. In case of severe condition consult a physician immediately and present this safety data sheet. Eyewash and shower for the treatment of emergency have to be present in the workplace.
Skin contact:	Take off the contaminated clothes and wash with copious amounts of water. If irritation persist consult a physician.
Eye contact:	Take off lenses (if present). Flush eyes with plenty of water for at least 15 minutes keeping eyelids opened. Consult immediately a physician.
Inhalation:	In case of severe condition consult a physician.
Ingestion:	Clean mouth with water and drink afterwards plenty of water. Do not induce vomiting. Do not administer any medicine if the person is unconscious. Consult immediately a physician.

4.2 MOST IMPORTANT SYMPTOMS AND EFFECTS

No data available

4.3 INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

No data available



(According to Regulation (EC) No. 1907/2006 e 453/2010)

5. FIREFIGHTING MEASURES

5.1 EXTINGUISHING MEDIA

 Suitable extinguishing media:
 Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

 Unsuitable extinguishing media:
 None

5.2 SPECIAL HAZARDS ARISING FROM THE SUBSTANCE OR MIXTURE

The mixture is not inflammable. However, in case of fire, irritating and/or toxic fumes can be released.

5.3 Advice for firefighters

If necessary, wear self-contained breathing apparatus. Avoid that fire extinguishing water contaminates surface water and/or groundwater.

6. ACCIDENTAL RELEASE MEASURES

6.1 PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES

Keep people not involved in the intervention operations away. Ensure an adequate ventilation of the affected area. Wash with plenty of water all the body parts entered in contact with the mixture immediately. Do not handle damaged containers or the leaked product before wearing appropriate protective outfit.

6.2 ENVIRONMENTAL PRECAUTIONS

Avoid the contamination of surface water and soil.

6.3 METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP

Absorb on inert material (sand/earth). Collect the spilled material with mechanical means and store it in suitable containers for disposal. Use water only to remove residuals, to avoid the danger of spill of product into the sewers.

6.4 REFERENCE TO OTHER SECTIONS

For more information on disposal see section 8 and 13.

7. HANDLING AND STORAGE

7.1 PRECAUTION FOR SAFE HANDLING

Work in well-ventilated area. Avoid the contact with eyes and skin and the inhalation of vapors.

7.2 CONDITION FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES

Store at 2-30°C, away from basic substances.

7.3 SPECIFIC END USE

None

8. EXPOSURE CONTROL/PERSONAL PROTECTION

8.1 CONTROL PARAMETERS

Applicable exposure limits	D	-	TWA/8h		STEL/15min	
	Description	Туре	mg/m3	ppm	mg/m3	ppm
	Sulphuric Acid	TLV-ACGIH 2008	0.2			
8.2 EXPOSURE CONTROLS	Do not eat, drink	according to the usual prec or smoke while handling th meals and after the work s	he product;			
Appropriate engineering controls:	Ensure an adequa	ate ventilation of the workin	ig area.			



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3 Date: 08-Jul-2014 Page: 6/9

Personal Protective Equipment:	Suggestions on the use of specific PPE are indicated below. Their choice should be made according to the use of the product and the instructions given by the supplier of the device.
	Hand protections: Chemical-resistant gloves, comply with EN 374.
	Eye protections: Safety goggles side shields
	Body protections: Work outfits for complete skin protection
	Respiratory protection: Not required under normal work activities
3.3 ENVIRONMENTAL EXPOSURE	
CONTROLS	Do not discard residuals in the environment. Discard the empty containers as hazardous waste.

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 INFORMATION ON BASIC PHYSICAL AND CHEMICAL PROPERTIES

Form Liquid Odour No data available Odour Threshold No data available < 0.5 pН Freezing/Melting point No data available Intital boiling and melting point range No data available No data available Flash point Evaporation rate No data available Flammability (solid, gas) Not applicable Upper/lower flammability or explosive limit Not applicable Vapour pressure No data available Vapour density No data available Relative density No data available Water solubility Soluble in water Partition coefficient: n-octano/water No data available Autoignition temperature Not applicable Decomposition temperature No data available No data available Viscosity Explosive properties No data available Oxidizing properties No data available

9.2 OTHER SAFETY INFORMATION

None

10. STABILITY E REACTIVITY

10.1 REACTIVITY

No data available

10.2 CHEMICAL STABILITY

Stable under the recommended transport, handling and storage conditions.

10.3 POSSIBILITY OF HAZARDOUS REACTIONS

The mixture reacts with the basic and/or oxidizing substances with formation of heat.

10.4 CONDITIONS TO AVOID

Avoid contact with basic and/or oxidizing substances.

10.5 INCOMPATIBLE MATERIALS

Basic substances.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3 Date: 08-Jul-2014 Page: 7/9

10.6 HAZARDOUS DECOMPOSITION PRODUCTS

When heated or in case of fire, vapors that are potentially dangerous to health (SOx) may be produced.

11. TOXICOLOGICAL INFORMATION

11.1 INFORMATION ON TOXICOLOGICAL EFFECTS

The mixture is corrosive for skin and eyes.

Acute toxicity:	No data available
Skin irritation:	Corrosive for skin and mucous membranes
Eye irritation:	Severe eye damages
Subcutaneous sensitization:	No data available
Germ cell mutagenesis:	No data available
Carcinogenicity:	No data available
Reproduction toxicity:	No data available
Organ specific toxicity (on single or repeated exposure):	No data available
Acute inhalation toxicity:	No data available
, Other information	No data available

12. ECOLOGICAL INFORMATION

12.1 TOXICITY

Use according to good working practices, avoiding release to the environment. In case the product reach waterways or sewers or contaminate soil or vegetation, inform the competent authorities.

12.2 PERSISTENCE AND DEGRADABILITY

No data available

12.3 BIOACCUMULATIVE POTENTIAL

No data available

12.4 MOBILITY IN SOIL

No data available

12.5 RESULTS OF PBT AND VPVB ASSESSMENT

No data available

12.6 OTHER ADVERSE EFFECTS

The disposal into drains or water bodies can result in a decreasing of the pH value. A low pH value harms aquatic organisms.

13. DISPOSAL CONSIDERATIONS

13.1 WASTE TREATMENT METHODS

Serum samples and reagents used should be treated as infectious waste. The product and its containers should be considered special waste. Their transport and disposal should be performed by authorized specialized companies according to the law.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3 Date: 08-Jul-2014 Page: 8/9

14. TRANSPORT INFORMATION

Even if the mixture is classified as hazardous by the transport regulations (see table below), the exemptions to the classification of ADR Chapter 3.5 - *Dangerous goods packed in excepted quantities* – can be applied.

These goods must be transported by vehicles authorized to the transport of hazardous goods according to the requirements of the current edition of ADR and to the national regulations.

People, responsible for the loading and unloading of this dangerous good must be trained on all the risks deriving from the substance and on procedures to be taken in case of emergency.

ADR/RID:			
UN N°	2796	Class:	8
Classification:	C1	Packaging group:	II
Tunnel restriction code:	(E)		
If subject to the following character	istics is ADR exempt:		
Exemption for limited quantities: 1L			
Exemption for expected quantities:			
 Maximum quantity for the 	internal package: 30 m		
 Maximum quantity for the 	external package: 500	nl	
Sea transport (IMDG Code):			
UN N°	2796	Class:	8
Packaging group:	II	EmS:	F-A, S-B
Marine polluting agent	no		
Air transport (ICAO-IATA):			
Class:	8	Packaging group: II	
Name:	Sulphuric	Acid (not containing more than	51% of acid)
Label:	n.8		

15. REGULATORY INFORMATION

15.1 SAFETY, HEALTH AND ENVIRONMENTAL REGULATIONS/LEGISLATION SPECIFIC FOR THE SUBSTANCE OR MIXTURE

Classification and labeling is done in accordance with the directives 1272/2008 (EC) and its subsequent amendments.

Legislative Decree 81/2008:

The use of this chemical entails the obligation of "Risk Assessment" by the employer in accordance with the provisions of the Decree April 9, 2008 n. 81 and subsequent amendments. Workers exposed to this chemical agent must not undergo health checks whether the results of the risk assessment show that, depending on the type and quantity of chemical agent and the method and frequency of exposure to this agent, there is only one "Low risk to the safety and irrelevant to the health" of the workers and the measures provided for in the same Decree are sufficient to reduce the risk.

European Standards:1907/2006/EC:

Not applicable

<u>Directive 96/82/EC (Seveso Directive)</u>: Not applicable

15.2 CHEMICAL SAFETY ASSESSMENT

Not performed for the product.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

16. OTHER INFORMATION

This product has to be used for diagnostic use only by personnel who is qualified and trained on the hazards shown in this safety sheet.

Text of hazard statements set out in paragraph 3

H314 Causes severe skin burns and eye damage

Text of risk phrases set out in paragraph 3

R35 Causes severe burns.

REFERENCES:

- 1. Directive 1999/45/EC and following amendments
- 2. Directive 67/548/EEC and following amendments and technical adjustments
- 3. Regulation (EC) 1907/2006 of the European Parliament (REACH) as amended
- 4. Regulation (EC) 1272/2008 of the European Parliament (CLP) as amended
- 5. Commission Regulation (EU) 453/2010
- 6. ECHA European Chemicals Agency www.echa.europa.eu

7. The Merck Index.

8. Handling Chemical Safety

- 9. NIOSH Registry of Toxic Effects of Chemical Substances
- 10. INRS Fiche Toxicologique

DIESSE Diagnostica Senese SpA declares that the information contained in this data sheet is based on the knowledge available to us on the date of the last version. Users must verify the suitability and accuracy of the provided information according to the specific use of the product. The document has not be accepted as a guarantee of any specific product property. Because of the use of this product is not subject to our direct control, the users, under their own responsibility, must observe the laws and regulations relating to health and safety. We accept no responsibility for improper use.

This data sheet supersedes each previous edition.