



CATALOGUE **2020**

CATALOGUE 2020

**DIESSE
DIAGNOSTICA SENESE SPA**

Sales Office:
Strada dei Laghi 39
53035 Monteriggioni (SI), Italy

T. +39 0577 319574/75
F. +39 0577 319020

venditeitalia@diesse.it
www.diesse.it



Company with quality management system certified
according to the standards UNI EN ISO 9001:2015 e EN ISO
13485:2016





INTRODUCTION

DIESSE Diagnostica Senese SpA is an Italian manufacturer of innovative in vitro diagnostics systems. DIESSE was founded in 1980, DIESSE has focused its efforts on developing a solid position in the IVD market segments of hematology, seroimmunology and bacteriology.

A mainstay of DIESSE strategy is the in-house production of most of the basic biologic components of the diagnostic kits (polyclonal and monoclonal antibodies, native and recombinant antigens), along with the capacity of designing instruments and dedicated plastic devices in order to develop close automated systems.

CATALOGUE

CUBE NEXT	8
LINE	
VES-MATIC CUBE	10
LINE	
VES-MATIC ORIGINAL	12
LINE	
CHORUS	14
LINE	
ENZYWELL	34
LINE	
AUTO-DAT	38
LINE	



CUBE NEXT LINE instruments

CUBE 30 touch

Automatic continuous loading instrument for the determination of the ESR on 30 samples simultaneously directly on the hematocrit tube. Automatic barcode reading, internal mixing and internal printer.

1 pc

10395

MINI-CUBE

Automatic continuous loading instrument for the determination of the ESR on 4 samples simultaneously directly on the hematocrit tube. USB and Bluetooth connector for personal computer, smartphone and tablet. External bar-code reader and printer (optional).

1 pc

10392

accessories

Bluetooth MINI-CUBE Printer	1 pc	R21430320
Customized Bluetooth thermal printer for MINI-CUBE for automatic printing of results.		
ESR CONTROL CUBE	4x9 ml	10435
	2x9 ml	10436
Control blood for the ESR tests to be used with the VES-Matic Cube systems. Contains normal and pathological level vials and dedicated bar-code labels.		
EXTERNAL BARCODE READER	1 pc	P20550510
External bar-code reader for all instruments of Ves-matic Original, Ves-matic Cube line, MINI-CUBE and CUBE 30 touch.		
MINI-CUBE Thermal paper	1 pc	10404
Thermal paper for Bluetooth Printer MINI-CUBE		
Thermal paper	1 pc	10403
Paper rolls for thermal printer in Ves-Matic Cube, Ves-Matic Original and CUBE 30 touch instruments.		
TEST DEVICE NEXT 500	500 test	10293
TEST DEVICE NEXT 1K	1.000 test	10294
TEST DEVICE NEXT 5K	5.000 test	10296
TEST DEVICE NEXT 10K	10.000 test	10297
Dedicated tube with transponder RF. To be used only once to charge the MINI-CUBE and CUBE 30 touch systems with the exact number of tests (not compatible with the VES-MATIC CUBE instrument)		



VES-MATIC CUBE LINE

instruments

VES-MATIC CUBE 200

Automatic instrument for the determination of the ESR directly on the hematocrit tube (190 samples/hour), using racks from the main cell counters. The first result is performed after 24 minutes. Host query interface.

1 pc

10370

Advia	1 pc	10370/A
ABX	1 pc	10370/AB
Beckman Coulter	1 pc	10370/BC
Beckman Coulter (5 positions)	1 pc	10370/BC5
Sysmex	1 pc	10370/S

VES-MATIC CUBE 80

Automatic instrument for the determination of the ESR directly on the hematocrit tube (95 samples/hour), using dedicated racks. The first result is performed after 26 minutes. Host query interface.

1 pc

10380

VES-MATIC CUBE TRACK

Automatic instrument for the determination of the ESR directly on the hematocrit tube (95 samples/hour), to be connected to clinical laboratory total automation tracks.

1 pc

10385

accessories

ESR CONTROL CUBE

Control blood for the ESR tests to be used with the VES-Matic Cube and CUBE NEXT line systems. Contains normal and pathological level vials and dedicated bar-code labels.

4x9 ml
2x9 ml

10435
10436

EXTERNAL BARCODE READER

External bar-code reader for all instruments of Ves-matic Original, Ves-matic Cube line, MINI-CUBE and CUBE 30 touch.

1 pc **P20550510**

TEST DEVICE 1K

1.000 test

10292

TEST DEVICE 5K

5.000 test

10291

TEST DEVICE 10K

10.000 test

10290

Dedicated tube with transponder RF. To be used only once to charge the VES-Matic Cube systems with the exact number of tests (not compatible with the MINI-CUBE and CUBE 30 touch instrument).

THERMAL PAPER

Paper rolls for thermal printer in Ves-Matic Cube, Ves-Matic Original and CUBE 30 touch instruments.

1 pc **10403**



VES-MATIC ORIGINAL LINE instruments

VES-MATIC EASY

Automatic instrument for the determination of the ESR on 10 samples simultaneously and random access. Complete with printer. USB connector and external bar code reader (optional) connector.

1 pc

10376

accessories

ESR CONTROL

Control blood for the ESR tests to be used with the Ves-matic system. Contains normal and pathological level vials.

**4x9 ml
2x9 ml**

**10430
10434**

EXTERNAL BARCODE READER

External bar-code reader for all instruments of Ves-matic Original, Ves-matic Cube line, MINI-CUBE and CUBE 30 touch.

1 pc

P20550510

VES-MATIC ORIGINAL LINE

accessories

THERMAL PAPER

1 pz.

10403

Paper rolls for thermal printer in Ves-Matic Cube, Ves-Matic Original and CUBE 30 touch instruments.

blood drawing and analysis system

VACU-CODE

100 test

10230

Vacuum tubes for blood drawing, without outertubes, containing sodium citrate solution for direct determination of the ESR. Suitable for bar code label on primary tube. To be used with VES-MATIC "PLUS". Sample volume: 1 ml.

VACU-TEC

100 test

10200

Vacuum tubes for blood drawing, with outertubes, containing sodium citrate solution for direct determination of the ESR. Sample volume: 1 ml.

VACU-TEC CUVETTE-L

200 test

10208

Vacuum tubes for blood drawing, without outertubes, containing sodium citrate solution for direct determination of the ESR. Complete with identification label. Sample volume: 1 ml.

VACU-TEC DOUBLE LABEL

100 test

10200/W

Vacuum tubes for blood drawing, with outertubes, containing sodium citrate solution for direct determination of the ESR and for a positive identification of the sample. Sample volume: 1 ml.

VACU-TEC FLAG

200 tests

10235

Vacuum tubes for blood drawing, without outertubes, containing sodium citrate solution for direct determination of the ESR. Complete with identification label. Labels for bar code support included in the package. Sample volume: 1 ml.

VACU-TEC S.C.

100 tests

10600

100 sterile tests

10602

Vacuum tubes for blood drawing, with outertubes, joined to a safety collar where a double label is applied, containing sodium citrate solution for direct determination of the ESR. Sample volume: 1 ml.

VES-TEC

100 tests

10201/A

Tubes for blood drawing, with outertubes, containing sodium citrate solution for direct determination of the ESR. Sample volume: 1 ml.

VES-TEC CUVETTE-L

200 tests

10214

Tubes for blood drawing, without outertubes, containing sodium citrate solution for direct determination of the ESR. Complete with identification label. Sample volume: 1 ml.



CHORUS LINE instrument

Chorus TRIO

Easy-to-use single test multiparametric system for immunometric assay.

1 pc

81200

Chorus TRIO RAS

Easy-to-use single test multiparametric system for immunometric assay. Prepared for the Autosampler accessory (81210).

1 pc

81220

accessories

Chorus AUTOSAMPLER

Integrated automated pipettor instrument for dispensing of the samples in diagnostic devices, to be associated with the Chorus Trio RAS (81220).

1 pc

81210



Chorus AUTOSAMPLER UPGRADE KIT

Integrated automated pipettor instrument for dispensing of the samples in diagnostic devices, to be associated with the Chorus Trio (81200).

1 pc

81215

THERMAL PAPER

Paper rolls for thermal printer.

4 pcs

83610

TIPS FOR CHORUS AUTOSAMPLER

Tips for automation in Autosampler.

10X96 pcs

83611

common reagents

CHORUS NEGATIVE CONTROL/SAMPLE DILUENT	1x5 ml	83607
Solution to be used as negative control sample and as diluent for highly reactive samples.		
CLEANING SOLUTION 2000	4x20 ml	83609
Cleaning solution for washing tubing system and sampling needles. 2000 X concentrated.		
SANITIZING SOLUTION	1x20 ml	83604
Sanitizing solution for sanitizing tubing system and sampling needles. 1000 X concentrated.		
WASHING BUFFER	4x100 ml	83606
Cleaning solution for washing of reaction wells. 10 X concentrated.		
WASHING BUFFER AUTOIMMUNITY	12x20 ml	86004
Cleaning solution for washing of reaction wells. 50 X concentrated.		

infectious diseases

ToRCH

Chorus CYTOMEGALOVIRUS IgG **36 tests** **81010**

Kit for the quantitative determination of anti-Cytomegalovirus IgG antibodies in human serum. Contains calibrator and positive control.

Chorus CYTOMEGALOVIRUS IgG AVIDITY **12 tests** **81092**

Kit for the qualitative determination of anti-Cytomegalovirus IgG antibodies avidity in human serum. Contains avidity controls.

Chorus CYTOMEGALOVIRUS IgM **36 tests** **81011**

Kit for the qualitative determination of anti-Cytomegalovirus IgM antibodies in human serum. Contains calibrator and positive control.

NEW! **Chorus CYTOMEGALOVIRUS IgM CAPTURE** **36 tests** **81012**

Kit for the qualitative determination by capture method of anti Cytomegalovirus IgM antibodies in human serum. Contains calibrator and positive control.

Chorus HERPES SIMPLEX 1+2 IgG **36 tests** **81020**

Kit for the qualitative determination of anti-Herpes simplex (type 1 and 2) IgG antibodies in human serum. Contains calibrator and positive control.

Chorus HERPES SIMPLEX 1+2 IgM **36 tests** **81021**

Kit for the qualitative determination of anti-Herpes simplex (type 1 and 2) IgM antibodies in human serum. Contains calibrator and positive control.

Chorus HERPES SIMPLEX 1 IgG Recombinant **36 tests** **81023**

Kit for the qualitative determination of anti-Herpes simplex type 1 IgG antibodies in human serum, using recombinant antigens. Contains calibrator and positive control.

Chorus HERPES SIMPLEX 2 IgG Recombinant **36 tests** **81029**

Kit for the qualitative determination of anti-Herpes simplex type 2 IgG antibodies in human serum, using recombinant antigens. Contains calibrator and positive control.

Chorus HSV 1 SCREEN **36 tests** **81018**

Kit for the qualitative determination of anti-Herpes simplex type 1 immunoglobulins in human serum. Immunoenzymatic competitive method. Contains calibrator and positive control.

* available soon

infectious diseases

	<p>Chorus HSV 2 SCREEN</p> <p>Kit for the qualitative determination of anti-Herpes simplex type 2 immunoglobulins in human serum. Immunoenzymatic competitive method. Contains calibrator and positive control.</p>	36 tests	81019
	<p>Chorus RUBELLA IgG</p> <p>Kit for the qualitative determination of anti-Rubella IgG antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81030
	<p>Chorus RUBELLA IgG AVIDITY</p> <p>Kit for the qualitative determination of anti-Rubella IgG antibodies avidity in human serum. Contains avidity controls.</p>	12 tests	81095
	<p>Chorus RUBELLA IgM</p> <p>Kit for the qualitative determination of anti-Rubella IgM antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81031
	<p>Chorus TOXOPLASMA IgA</p> <p>Kit for the qualitative determination of anti-Toxoplasma gondii IgA antibodies in human serum. Contains calibrator and positive control.</p>	36 tests 12 tests	81043 81043/12 Upon request
NEW!	<p>Chorus TOXOPLASMA IgA CAPTURE</p> <p>Kit for the qualitative determination by capture method of anti-Toxoplasma gondii IgA antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81044
	<p>Chorus TOXOPLASMA IgG</p> <p>Kit for the quantitative determination of anti-Toxoplasma gondii IgG antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81040
	<p>Chorus TOXOPLASMA IgG AVIDITY</p> <p>Kit for the qualitative determination of anti-Toxoplasma gondii IgG antibodies avidity in human serum. Contains avidity controls.</p>	12 tests	81098
	<p>Chorus TOXOPLASMA IgM</p> <p>Kit for the qualitative determination of anti-Toxoplasma gondii IgM antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81041
NEW!	<p>Chorus TOXOPLASMA IgM CAPTURE</p> <p>Kit for the qualitative determination by capture method of anti-Toxoplasma gondii IgM antibodies in human serum. Contains calibrator and positive control.</p>	36 tests	81042

CHILDHOOD DISEASES

Chorus MEASLES IgG Kit for the qualitative determination of anti-Measles IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81072
Chorus MEASLES IgM Kit for the qualitative determination of anti-Measles IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81073
Chorus MUMPS IgG Kit for the qualitative determination of anti-Mumps IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81075
Chorus MUMPS IgM Kit for the qualitative determination of anti-Mumps IgM antibodies in human serum. Contains calibrator and positive control.	36 tests 12 tests	81076 81076/12
Chorus PARVOVIRUS B19 IgG Kit for the qualitative determination of anti-Parvovirus B19 IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81150
Chorus PARVOVIRUS B19 IgM Kit for the qualitative determination of anti- Parvovirus B19 IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81152
Chorus VARICELLA IgG Kit for the qualitative determination of anti-Varicella Zoster IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81078
Chorus VARICELLA IgM Kit for the qualitative determination of anti-Varicella Zoster IgM antibodies in human serum. Contains calibrator and positive control.	36 tests 12 tests	81079 81079/12

VACCINABLES

Chorus DIPHTHERIA IgG Kit for the quantitative determination of IgG class antibodies against Diphtheria toxin in human serum and plasma. Contains calibrator and positive control.	36 tests	81266
Chorus POLIOVIRUS IgG Kit for the quantitative determination of IgG class antibodies against Poliovirus in human serum and plasma. Contains calibrator and positive control.	36 tests	81264

infectious diseases

Chorus TETANUS IgG

Kit for the quantitative determination of IgG antibodies against Tetanus toxin in human serum and plasma. Contains calibrator and positive control.

36 tests
12 tests

81270
81270/12
Upon request

BACTERIAL RESPIRATORY INFECTIONS

NEW!

Chorus BORDETELLA PERTUSSIS TOXIN IgA

Kit for the quantitative determination of anti-Bordetella pertussis toxin IgA antibodies.

36 tests

81274

NEW!

Chorus BORDETELLA PERTUSSIS TOXIN IgG

Kit for the quantitative determination of anti-Bordetella pertussis toxin IgG antibodies.

36 tests

81272

Chorus CHLAMYDOPHILA PNEUMONIAE IgA

Kit for the semiquantitative determination of anti-Chlamydomphila Pneumoniae IgA antibodies in human serum. Contains calibrator and positive control.

36 tests

81252

Chorus CHLAMYDOPHILA PNEUMONIAE IgG

Kit for the semiquantitative determination of anti-Chlamydomphila pneumoniae IgG antibodies in human serum. Contains calibrator and positive control.

36 tests

81250

Chorus CHLAMYDOPHILA PNEUMONIAE IgM

Kit for the semiquantitative determination of anti-Chlamydomphila pneumoniae IgM antibodies in human serum. Contains calibrator and positive control.

36 tests

81251

Chorus LEGIONELLA PNEUMOPHILA 1-6 IgG

Kit for the qualitative determination of anti-Legionella pneumophila serogroups 1-6 IgG antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81090
81090/12
Upon request

Chorus LEGIONELLA PNEUMOPHILA 1 IgG

Kit for the qualitative determination of anti-Legionella pneumophila serogroup 1 IgG antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81091
81091/12
Upon request

Chorus LEGIONELLA PNEUMOPHILA IgM

Kit for the qualitative determination of anti-Legionella pneumophila IgM antibodies in human urine. Contains calibrator and positive control.

36 tests
12 tests

81093
81093/12
Upon request

Chorus MYCOPLASMA PNEUMONIAE IgA **36 tests** **81033**

Kit for the semiquantitative determination of anti-Mycoplasma pneumoniae IgA antibodies in human serum. Contains calibrator and positive control.

Chorus MYCOPLASMA PNEUMONIAE IgG **36 tests** **81034**

Kit for the semiquantitative determination of anti-Mycoplasma pneumoniae IgG antibodies in human serum. Contains calibrator and positive control.

Chorus MYCOPLASMA PNEUMONIAE IgM **36 tests** **81035**

Kit for the qualitative determination of anti-Mycoplasma pneumoniae IgM antibodies in human serum. Contains calibrator and positive control.

NEW!

Chorus Q FEVER PHASE II IgG **36 tests** **81171**

Kit for the qualitative determination of anti-Coxiella burnetii IgG antibodies in human serum. Contains calibrator and positive control.

NEW!

Chorus Q FEVER PHASE II IgM **36 tests** **81172**

Kit for the qualitative determination of anti-Coxiella burnetii IgM antibodies in human serum. Contains calibrator and positive control.

VIRAL RESPIRATORY INFECTIONS

Chorus ADENOVIRUS IgA **36 tests** **81198**
12 tests **81198/12**
 Kit for the qualitative determination of anti-Adenovirus IgA antibodies in human serum. Contains calibrator and positive control. Upon request

Chorus ADENOVIRUS IgG **36 tests** **81196**
12 tests **81196/12**
 Kit for the qualitative determination of anti-Adenovirus IgG antibodies in human serum. Contains calibrator and positive control. Upon request

Chorus INFLUENZA A IgA **36 tests** **81192**
12 tests **81192/12**
 Kit for the qualitative determination of anti-Influenza A IgA antibodies in human serum. Contains calibrator and positive control. Upon request

Chorus INFLUENZA A IgG **36 tests** **81190**
12 tests **81190/12**
 Kit for the qualitative determination of anti-Influenza A IgG antibodies in human serum. Contains calibrator and positive control. Upon request

Chorus INFLUENZA B IgA **36 tests** **81195**
12 tests **81195/12**
 Kit for the qualitative determination of anti-Influenza B IgA antibodies in human serum. Contains calibrator and positive control. Upon request

infectious diseases**Chorus INFLUENZA B IgG**

Kit for the qualitative determination of anti-Influenza B IgG antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81193
81193/12
Upon request

Chorus RESPIRATORY SYNCYTIAL VIRUS IgA

Kit for the qualitative determination of anti-Respiratory Syncytial Virus IgA antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81037
81037/12
Upon request

Chorus RESPIRATORY SYNCYTIAL VIRUS IgG

Kit for the qualitative determination of anti-Respiratory Syncytial Virus IgG antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81036
81036/12
Upon request

SEXUALLY TRANSMITTED DISEASES**Chorus CHLAMYDIA TRACHOMATIS IgG**

Kit for the qualitative determination of anti- Chlamydia Trachomatis IgG antibodies in human serum. Contains calibrator and positive control.

36 tests

81261

Chorus CHLAMYDIA TRACHOMATIS IgA

Kit for the qualitative determination of anti- Chlamydia Trachomatis IgA antibodies in human serum. Contains calibrator and positive control.

36 tests

81262

Chorus HERPES SIMPLEX 2 IgG Recombinant

Kit for the qualitative determination of anti-Herpes simplex type 2 IgG antibodies in human serum, using recombinant antigens. Contains calibrator and positive control.

36 tests

81029

Chorus HSV 2 SCREEN

Kit for the qualitative determination of anti-Herpes simplex type 2 immunoglobulins in human serum. Immunoenzymatic competitive method. Contains calibrator and positive control.

36 tests

81019

Chorus SYPHILIS SCREEN RECOMBINANT

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) immunoglobulins in human serum, using recombinant antigens. Immunoenzymatic method. Contains calibrator and positive control.

36 tests

81100

Chorus TREPONEMA IgG

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) IgG antibodies in human serum. Contains calibrator and positive control.

36 tests
12 tests

81050
81050/12
Upon request

Chorus TREPONEMA IgM

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) IgM antibodies in human serum. Capture method. Contains calibrator and positive control.

36 tests
12 tests

81051
81051/12
Upon request

GASTROINTESTINAL DISEASES**Chorus HELICOBACTER PYLORI IgA**

Kit for the qualitative determination of anti-Helicobacter pylori IgA antibodies in human serum. Contains calibrator and positive control.

36 tests

81062

Chorus HELICOBACTER PYLORI IgG

Kit for the semiquantitative determination of anti-Helicobacter pylori IgG antibodies in human serum. Contains calibrator and positive control.

36 tests

81060

INFECTIOUS MONONUCLEOSIS**Chorus EPSTEIN-BARR EARLY ANTIGEN IgG**

Kit for the qualitative determination of IgG antibodies against Early Antigen of Epstein-Barr Virus in human serum. Contains calibrator and positive control.

36 tests
12 tests

81058
81058/12
Upon request

Chorus EPSTEIN-BARR EARLY ANTIGEN IgM

Kit for the qualitative determination of IgM antibodies against Early Antigen of Epstein-Barr Virus in human serum. Contains calibrator and positive control.

36 tests
12 tests

81059
81059/12
Upon request

Chorus EPSTEIN-BARR EBNA IgG

Kit for the qualitative determination of IgG antibodies against Nuclear Antigen of Epstein-Barr Virus in human serum. Contains calibrator and positive control.

36 tests

81057

Chorus EPSTEIN-BARR VCA IgG

Kit for the qualitative determination of IgG antibodies against Viral Capsid Antigen of Epstein-Barr Virus in human serum. Contains calibrator and positive control.

36 tests

81055

Chorus EPSTEIN-BARR VCA IgM II

Kit for the qualitative determination of IgM antibodies against Viral Capsid Antigen of Epstein-Barr Virus in human serum. Contains calibrator and positive control.

36 tests

81054

infectious diseases

VECTOR TRANSMITTED DISEASES

Chorus CHIKUNGUNYA IgG	36 tests	81133
Kit for the qualitative determination of anti-Chikungunya IgG antibodies in human serum. Contains calibrator and positive control.		
Chorus CHIKUNGUNYA IgM	36 tests	81134
Kit for the qualitative determination of anti-Chikungunya IgM antibodies in human serum. Contains calibrator and positive control.		
Chorus DENGUE IgG	36 tests	81154
Kit for the qualitative determination of anti-Dengue IgG antibodies in human serum. Contains calibrator and positive control.		
Chorus DENGUE IgM	36 tests	81156
Kit for the qualitative determination of anti-Dengue IgM antibodies in human serum. Contains calibrator and positive control.		
Chorus TOSCANA VIRUS IgG	36 tests	81067
Kit for the qualitative determination of anti-Toscana Virus IgG antibodies in human serum. Contains calibrator and positive control.	12 tests	81067/12
		Upon request
Chorus TOSCANA VIRUS IgM	36 tests	81068
Kit for the qualitative determination of anti-Toscana Virus IgM antibodies in human serum. Capture method. Contains calibrator and positive control.	12 tests	81068/12
		Upon request

ZOONOSIS

Chorus BORRELIA IgG	36 tests	81046
Kit for the qualitative determination of anti-Borrelia bourgdorferi IgG antibodies in human serum. Contains calibrator and positive control.		
Chorus BORRELIA IgM	36 tests	81047
Kit for the qualitative determination of anti-Borrelia bourgdorferi IgM antibodies in human serum. Contains calibrator and positive control.		
Chorus BORRELIA IgG CSF	12 tests	81048
Kit for the qualitative determination of anti-Borrelia bourgdorferi IgG in liquor. Contains calibrator and positive control.		
NEW! Chorus BRUCELLA IgG	36 tests	81164
Kit for the qualitative determination of anti-Brucella IgG antibodies in human serum. Contains calibrator and positive control.		

* available soon

NEW!	Chorus BRUCELLA IgM Kit for the qualitative determination of anti-Brucella IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81166
NEW!	Chorus ECHINOCOCCUS IgG Kit for the qualitative determination of anti-Echinococcus IgG antibodies in human serum. Contains calibrator and positive control.	12 tests	81284
	Chorus LEPTOSPIRA IgG* Kit for the qualitative determination of anti-Leptospira IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81288
	Chorus LEPTOSPIRA IgM* Kit for the qualitative determination of anti-Leptospira IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81289
NEW!	Chorus Q FEVER PHASE II IgG Kit for the qualitative determination of anti-Coxiella burnetii IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81171
NEW!	Chorus Q FEVER PHASE II IgM Kit for the qualitative determination of anti-Coxiella burnetii IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81172
	Chorus TICK-BORNE ENCEPHALITIS VIRUS IgG Kit for the semiquantitative determination of anti-Tick-Borne Encephalitis Virus IgG antibodies in human serum. Contains calibrator and positive control.	36 tests	81276
	Chorus TICK-BORNE ENCEPHALITIS VIRUS IgM Kit for the semiquantitative determination of anti-Tick-Borne Encephalitis Virus IgM antibodies in human serum. Contains calibrator and positive control.	36 tests	81277
	Chorus TOxocara IgG Kit for the qualitative determination of anti-Toxocara IgG antibodies in human serum. Contains calibrator and positive control.	12 tests	81280

URINARY ANTIGENS

	Chorus LEGIONELLA URINARY ANTIGEN Kit for the qualitative determination of Legionella pneumophila urinary antigen in human urine. Contains calibrator and positive control.	12 tests	81301
--	---	-----------------	--------------

* available soon

infectious diseases

STOOL ANTIGENS

Chorus CALPROTECTIN Kit for the quantitative determination of Calprotectin in human stool. Contains calibrator and positive control.	12 tests	86600 86600/C03
Chorus CALPROTECTIN K Kit for the quantitative determination of Calprotectin in human stool. Contains calibrator and positive control. Range of calibration up to 1.000 mg/Kg.	12 tests	86601
Chorus CLOSTRIDIUM DIFFICILE A/B TOXINS Kit for the qualitative determination of A/B Toxins of Clostridium difficile in human stool. Contains calibrator and positive control.	12 tests	81170
Chorus CLOSTRIDIUM DIFFICILE GDH Kit for the qualitative determination of the Clostridium difficile GDH antigen in human stool. Contains calibrator and positive control.	12 tests	81168
Chorus COPROCOLLECT Device for the collection, transport and extraction of stool human samples. (For Calprotectin kit 86600, 86600/C03 and 86601).	20 pcs.	86602
Chorus HELICOBACTER PYLORI Ag Kit for the qualitative determination of Helicobacter pylori antigen in human stool. Contains calibrator and positive control.	12 tests	81063

BONE METABOLISM

Chorus 25 OH VITAMIN D TOTAL Kit for the quantitative determination of 25 (OH) vitamin D in human serum. Contains calibrator and positive control.	36 test	86900
--	----------------	--------------

autoimmunity

REUMATOLOGY

Chorus ANA-8 Kit for the quantitative determination of IgG autoantibodies to 8 different nuclear antigens: U1-snRNP 70 kDa, SS-B, SS-A 60/52 kDa, Scl-70, Cenp-B, Jo-1, snRNP complex (snRNP/Sm), Sm (screening) in human serum. Contains liquid calibrator and positive control.	36 tests 12 tests	86010 86010/12 Upon request
---	------------------------------------	---

Chorus ANA SCREEN

Kit for the qualitative determination of IgG autoantibodies against HEp2 cells extract enriched with purified nuclear antigens (dsDNA, SS-B, SS-A 60/52 kDa, Scl-70, Jo-1, snRNP complex (snRNP/Sm), Sm) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86014
86014/12
Upon request

Chorus anti-CCP

Kit for the semiquantitative determination of IgG autoantibodies against cyclic citrullinated peptides in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86094
86094/12
Upon request

Chorus Cenp-B

Kit for the semiquantitative determination of IgG autoantibodies to the 80 kDa centromere protein B (Cenp-B) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86028
86028/12
Upon request

Chorus CH50*

Kit for the quantitative determination of the serum complement haemolytic activity and of the anticomplementary power. Contains freeze-dried calibrator and positive control.

12 tests

86102

Chorus dsDNA-G

Kit for the quantitative determination of IgG autoantibodies to double stranded DNA (dsDNA) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86032
86032/12
Upon request

Chorus dsDNA-M

Kit for the semiquantitative determination of IgM autoantibodies to double stranded DNA (dsDNA) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86034
86034/12
Upon request

Chorus ENA-6 S

Kit for the qualitative determination of IgG autoantibodies to 6 different nuclear antigens: SS-B, SS-A 60/52 kDa, Scl-70, Jo-1, snRNP complex (snRNP/Sm), Sm (screening) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86012
86012/12
Upon request

Chorus Jo-1

Kit for the semiquantitative determination of IgG autoantibodies to human histidyl-tRNA synthetase (Jo-1) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86030
86030/12
Upon request

Chorus RF-G

Kit for the semiquantitative determination of IgG rheumatoid factors in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86038
86038/12
Upon request

* available soon

CHORUS LINE

autoimmunity

Chorus RF-M

Kit for the quantitative determination of IgM rheumatoid factors in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86040
86040/12

Upon request

Chorus Scl-70

Kit for the semiquantitative determination of IgG autoantibodies to the 70 kDa fragment of human DNA-topoisomerase I (Scl-70) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86026
86026/12

Upon request

Chorus SM

Kit for the semiquantitative determination of IgG autoantibodies to Sm proteins in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86016
86016/12

Upon request

Chorus snRNP-C

Kit for the semiquantitative determination of IgG autoantibodies to snRNP Complex in human serum (U1-snRNP Complex). Contains liquid calibrator and positive control.

36 tests
12 tests

86092
86092/12

Upon request

Chorus SS-A

Kit for the semiquantitative determination of IgG autoantibodies to the SS-A 60/52 kDa proteins (formerly called Ro) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86018
86018/12

Upon request

Chorus SS-B

Kit for the semiquantitative determination of IgG autoantibodies to the 48 kDa SS-B protein (formerly called La) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86024
86024/12

Upon request

Chorus U1-70 RNP

Kit for the semiquantitative determination of IgG autoantibodies to 70 kDa fragment of ribonucleoprotein U1 in human serum (U1-70 RNP). Contains liquid calibrator and positive control.

36 tests
12 tests

86088
86088/12

Upon request

THYROID

Chorus a-TG

Kit for the quantitative determination of IgG autoantibodies to thyroglobulin (Tg) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86072
86072/12

Upon request

Chorus a-TPO

Kit for the quantitative determination of IgG autoantibodies to thyroid peroxidase (TPO) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86074
86074/12

Upon request

THROMBOSIS

Chorus BETA 2-GLYCOPROTEIN-G

Kit for the semiquantitative determination of IgG autoantibodies to β 2-glycoprotein I in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86050
86050/12
Upon request

Chorus BETA 2-GLYCOPROTEIN-M

Kit for the semiquantitative determination of IgM autoantibodies to β 2- glycoprotein I in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86052
86052/12
Upon request

Chorus CARDIOLIPIN-G

Kit for the semiquantitative determination of IgG autoantibodies to the cardiolipin- β 2-glycoprotein I complex in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86046
86046/12
Upon request

Chorus CARDIOLIPIN-M

Kit for the semiquantitative determination of IgM autoantibodies to the cardiolipin- β 2-glycoprotein I complex in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86048
86048/12
Upon request

Chorus PHOSPHOLIPID SCREEN-G*

Kit for the semiquantitative determination of anti-phospholipid IgG autoantibodies in human serum. Contains liquid calibrator and positive control.

36 tests

86106

Chorus PHOSPHOLIPID SCREEN-M*

Kit for the semiquantitative determination of anti-phospholipid IgM autoantibodies in human serum. Contains liquid calibrator and positive control.

36 tests

86108

HEPATOLOGY

Chorus AMA-M2

Kit for the semiquantitative determination of anti-mitochondrial autoantibodies (IgG) to M2 in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86084
86084/12
Upon request

Chorus LKM-1

Kit for the semiquantitative determination of LKM-1 autoantibodies (IgG) to liver-kidney microsomes in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86080
86080/12
Upon request

* available soon

**CHORUS
LINE****autoimmunity****VASCULITIS****Chorus GBM**

Kit for the semiquantitative determination of IgG autoantibodies to human glomerular basement membrane protein (GBM) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86070
86070/12
Upon request

Chorus MPO

Kit for the semiquantitative determination of IgG autoantibodies to myeloperoxidase (MPO) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86068
86068/12
Upon request

Chorus PR3

Kit for the semiquantitative determination of IgG autoantibodies to proteinase 3 (PR3) in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86066
86066/12
Upon request

DIABETES**Chorus INSULIN**

Kit for the semiquantitative determination of IgG autoantibodies to insulin in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86086
86086/12
Upon request

GASTROENTEROLOGY**Chorus ASCA-A**

Kit for the semiquantitative determination of IgA autoantibodies to cell wall proteoglycans of *Saccharomyces cerevisiae* in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86062
86062/12
Upon request

Chorus ASCA-G

Kit for the semiquantitative determination of IgG autoantibodies to cell wall proteoglycans of *Saccharomyces cerevisiae* in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86064
86064/12
Upon request

Chorus DEAMIDATED GLIADIN-A

Kit for the semiquantitative determination of IgA autoantibodies against deamidated gliadin peptides in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86098
86098/12
Upon request

Chorus DEAMIDATED GLIADIN-G

Kit for the semiquantitative determination of IgG autoantibodies against deamidated gliadin peptides in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86100
86100/12
Upon request

Chorus GLIADIN-A

Kit for the semiquantitative determination of IgA autoantibodies to alpha gliadin in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86054
86054/12
Upon request

Chorus GLIADIN-G

Kit for the semiquantitative determination of IgG autoantibodies to alpha gliadin in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86056
86056/12
Upon request

Chorus INTRINSIC FACTOR

Kit for the semiquantitative determination of IgG autoantibodies to human intrinsic factor in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86090
86090/12
Upon request

Chorus tTg-A

Kit for the semiquantitative determination of IgA autoantibodies to tissue-transglutaminase in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86058
86058/12
Upon request

Chorus tTg-G

Kit for the semiquantitative determination of IgG autoantibodies to tissue-transglutaminase in human serum. Contains liquid calibrator and positive control.

36 tests
12 tests

86060
86060/12
Upon request

endocrinology

THYROID

NEW!**Chorus FT3**

Kit for the quantitative determination of free triiodothyronine (fT3) in human serum. Contains calibrator and positive control.

36 tests**86518****NEW!****Chorus FT4**

Kit for the quantitative determination of free thyroxine (fT4) in human serum. Contains calibrator and positive control.

36 tests**86520****NEW!****Chorus TSH**

Kit for the quantitative determination of human thyroid-stimulating hormone (hTSH) in human serum. Contains calibrator and positive control.

36 tests**86516**

CHORUS LINE

endocrinology

FERTILITY

NEW!

Chorus FSH

36 tests

86510

Kit for the quantitative determination of Follicle-Stimulating hormone (FSH) in human serum. Contains calibrator and positive control.

NEW!

Chorus LH

36 tests

86512

Kit for the quantitative determination of Luteneizing hormone (LH) in human serum. Contains calibrator and positive control.

NEW!

Chorus PRL

36 tests

86514

Kit for the quantitative determination of Prolactine hormone (PRL) in human serum. Contains calibrator and positive control.

tumor markers

NEW!

Chorus PSA

36 tests

86550

Kit for the quantitative determination of the total Prostate Specific Antigen (tPSA) in human serum. Contains calibrator and positive control.

Chorus TG

36 tests

86076

Kit for the quantitative determination of thyroglobulin (Tg) in human serum. Contains liquid calibrator and positive control.

12 tests

86076/12

Upon request

NEW!

Chorus TG HIGH SENSITIVITY

36 tests

86104

Kit for the high sensitivity quantitative determination of thyroglobulin (Tg) in human serum. Contains calibrator and positive control.

12 tests

86104/12

Upon request

allergy

NEW!

Chorus TOTAL IgE

36 tests

88000

Kit for the quantitative determination of total IgE in human serum, complete with calibrator and positive control.

complement fixation test

BACTERIAL SEROLOGY ANTIGENS

Chorus BRUCELLA	12 tests	85358
Chorus LEPTOSPIRA MIX	12 tests	85367
Chorus LISTERIA MONOCYTOGENES	12 tests	85344

NEUROTROPIC ANTIGENS

Chorus ECHO VIRUS N MIX	12 tests	85275
--------------------------------	-----------------	--------------

PNEUMOTROPIC ANTIGENS

Chorus ECHO VIRUS P MIX	12 tests	85221
Chorus PARAINFLUENZA MIX	12 tests	85225
Chorus Q-FEVER	12 tests	85330

COXSACKIE ANTIGEN

Chorus COXSACKIE VIRUS A MIX	12 tests	85313
Chorus COXSACKIE VIRUS B MIX	12 tests	85251

control sera

Chorus COXSACKIE B CONTROL SERUM	0,5 ml	85903
Chorus LEPTOSPIRA CONTROL SERUM	0,5 ml	85896
Chorus LISTERIA MONOCYTOGENES CONTROL SERUM	0,5 ml	85875
Chorus Q-FEVER CONTROL SERUM	0,5 ml	85886
Chorus PARAINFLUENZA CONTROL SERUM	0,5 ml	85893

infectious diseases

ToRCH

CMV SCREEN Kit for the qualitative determination of anti-Cytomegalovirus immunoglobulins in human serum. Immunoenzymatic competitive method.	96 tests	91013
CYTOMEGALOVIRUS IgG Kit for the qualitative and quantitative determination of anti-Cytomegalovirus IgG antibodies in human serum. Contains ready for use calibrators for the preparation of the calibration curve in International Units (IU/ml).	96 tests	91010
CYTOMEGALOVIRUS IgG AVIDITY Kit for the determination of anti-Cytomegalovirus IgG antibodies avidity in human serum.	24 tests	91092
CYTOMEGALOVIRUS IgM Kit for the determination of anti-Cytomegalovirus IgM antibodies in human serum. Capture method.	96 tests	91011
HERPES SIMPLEX 1+2 IgG Kit for the qualitative determination of anti-Herpes simplex (type 1 and 2) IgG antibodies in human serum.	96 tests	91020
HERPES SIMPLEX 1+2 IgM Kit for the qualitative determination of anti-Herpes simplex (type 1 and 2) IgM antibodies in human serum. Capture method.	96 tests	91021
HERPES SIMPLEX 1 IgG Kit for the qualitative determination of anti-Herpes simplex type 1 IgG antibodies in human serum.	96 tests	91027
HERPES SIMPLEX 1 IgG Recombinant Kit for the qualitative determination of IgG-class antibodies to Herpes Simplex Virus (type 1)	96 tests	91023
HERPES SIMPLEX 1 IgM Kit for the qualitative determination of anti-Herpes simplex type 1 IgM antibodies in human serum. Capture method.	96 tests	91028
HERPES SIMPLEX 2 IgG Kit for the qualitative determination of anti-Herpes simplex type 2 IgG antibodies in human serum.	96 tests	91026
HERPES SIMPLEX 2 IgG Recombinant Kit for the qualitative determination of IgM-class antibodies to Herpes Simplex Virus (type 2)	96 tests	91029

HERPES SIMPLEX 2 IgM	96 tests	91025
Kit for the qualitative determination of anti-Herpes simplex type 2 IgM antibodies in human serum. Capture method.		
TOXOPLASMA IgG	96 tests	91040
Kit for the qualitative and quantitative determination of anti-Toxoplasma gondii IgG antibodies in human serum. Contains liquid ready for use calibrators for the preparation of the calibration curve in International Units (IU/ml).		
TOXOPLASMA IgG AVIDITY	24 tests	91098
Kit for the determination of anti-Toxoplasma gondii IgG antibodies avidity in human serum.		
TOXOPLASMA IgA	96 tests	91043
Kit for the qualitative determination of anti-Toxoplasma gondii IgA antibodies in human serum.		
TOXOPLASMA IgM	96 tests	91041
Kit for the qualitative determination of anti-Toxoplasma gondii IgM antibodies in human serum. Capture method.		
RUBELLA IgG	96 tests	91030
Kit for the qualitative and quantitative determination of anti-Rubella IgG antibodies in human serum. Contains liquid ready for use calibrators for the preparation of the calibration curve in International Units (IU/ml).		
RUBELLA IgG AVIDITY	24 tests	91095
Kit for the determination of anti-Rubella IgG antibodies avidity in human serum.		
RUBELLA IgM	96 tests	91031
Kit for the qualitative determination of anti-Rubella IgM antibodies in human serum. Capture method.		
CHILDHOOD INFECTIONS		
MEASLES IgG	48 tests	91072
Kit for the qualitative determination of anti-Measles IgG antibodies in human serum.		
MEASLES IgM	48 tests	91073
Kit for the qualitative determination of anti-Measles IgM antibodies in human serum. Capture method.		
MUMPS IgG	48 tests	91075
Kit for the qualitative determination of anti-Mumps IgG antibodies in human serum.		

infectious diseases

MUMPS IgM

Kit for the qualitative determination of anti-Mumps IgM antibodies in human serum. Capture method.

48 tests
91076

Upon request

VARICELLA IgG

Kit for the qualitative determination of anti-Varicella Zoster IgG antibodies in human serum.

48 tests
91078

Upon request

VARICELLA IgM

Kit for the qualitative determination of anti-Varicella Zoster IgM antibodies in human serum. Capture method.

48 tests
91079

Upon request

SEXUALLY TRANSMITTED DISEASES

HERPES SIMPLEX 2 IgG

Kit for the qualitative determination of anti-Herpes simplex type 2 IgG antibodies in human serum.

96 tests
91026

HERPE SIMPLEX 2 IgG RECOMBINANT

Kit for the qualitative determination of anti-Herpes simplex type 2 IgG antibodies in human serum, using recombinant antigens.

96 tests
91029

HERPES SIMPLEX 2 IgM

Kit for the qualitative determination of anti-Herpes simplex type 2 IgM antibodies in human serum. Capture method.

96 tests
91025

SYPHILIS SCREEN RECOMBINANT

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) immunoglobulins in human serum, using recombinant antigens. Immunoenzymatic method.

2x96 tests
91100

TREPONEMA IgG

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) IgG antibodies in human serum. Recombinant antigen

96 tests
91050

TREPONEMA IgM

Kit for the qualitative determination of anti-Treponema pallidum (Nichols strain) IgM antibodies in human serum. Recombinant antigen. Capture method.

96 tests
91051

GASTROINTESTINAL DISEASES

HELICOBACTER PYLORI IgA **96 tests** **91062**
 Kit for the qualitative determination of anti-Helicobacter pylori IgA antibodies in human serum. Upon request

HELICOBACTER PYLORI IgG **96 tests** **91060**
 Kit for the qualitative and semiquantitative determination of anti-Helicobacter pylori IgG antibodies in human serum. Contains ready for use calibrators for the preparation of the calibration curve in Arbitrary Units (AU/ml).

INFECTIOUS MONONUCLEOSIS

EPSTEIN-BARR EBNA IgG **96 tests** **91057**
 Kit for the semiquantitative determination of IgG antibodies against Nuclear Antigen of Epstein-Barr Virus in human serum.

EPSTEIN-BARR VCA IgG **96 tests** **91055**
 Kit for the qualitative determination of IgG antibodies against Viral Capsid Antigen of Epstein-Barr Virus in human serum. Upon request

EPSTEIN-BARR VCA IgM **96 tests** **91056**
 Kit for the semiquantitative determination of IgM antibodies against Viral Capsid Antigen of Epstein-Barr Virus in human serum.

VECTOR TRANSMITTED DISEASES

CHIKUNGUNYA IgG **96 tests** **91133**
 Kit for the qualitative determination of anti-Chikungunya IgG antibodies in human serum. Upon request

CHIKUNGUNYA IgM **96 tests** **91134**
 Kit for the qualitative determination of anti-Chikungunya IgM antibodies in human serum.

TOSCANA VIRUS IgG-IgM (Sandfly Fever Virus) **48+48 tests** **91067**
 Kit for the qualitative determination of anti-Toscana Virus IgG and IgM antibodies in human serum.



AUTO-DAT LINE instrument

AUTO-DAT*

Instrument for the automation of the incubation phase of sero agglutination tests. Semiquantitative interpretation of results and through image analysis software.

1 pc

26000

reagents

AUTO-DAT MICROPLATE*	20 pcs.	26005
To be used for sero agglutination reaction with AUTO-DAT instrument.		
RPR-DAT*	500 tests	26030
Non Treponemal agglutination test for syphilis screening. Automated method. Includes controls and dedicated microplates.		
TPHA-DAT*	50 tests	26035
Kit for the semi-quantitative determination of anti-Treponema pallidum antibodies in human serum or plasma, by indirect hemagglutination. Automated method, complete with dedicated controls and microplates.		
W.W.-DAT*	700 tests	26010
For serological diagnosis of bacterial infection (Widal-Wright) Panel of 7 coloured suspensions (TO, TH, AO, AH, BO, BH, Brucella)		
W.F.-DAT*	180 tests	26020
For serological diagnosis of Rickettsial infection.		

* available soon

QUALITY CERTIFICATIONS

DIESSE Diagnostica Senese SpA has been operating for many years in accordance with leading standards for quality management systems, ensuring traceability, transparency, highest quality of its processes, products and services and working on continuous improvement.

The Company is certified in accordance with UNI EN ISO 9001: 2015 and EN ISO 13485: 2016, also in compliance with the MDSAP (Medical Device Single Audit Program) program.

As part of this program, DIESSE Diagnostica Senese SpA has obtained, starting from 2019, the MDSAP certification in accordance with the Regulations of all 5 countries participating to the program: the United States, Canada, Australia, Japan and Brazil. All products are CE marked in accordance with the European Directive 98/79/EC concerning In Vitro Diagnostic Medical Devices and the relative applicable harmonized standards.

Starting from June 2020 DIESSE Diagnostica Senese S.p.A. has adopted an Organization, Management and Control Model in accordance with Legislative Decree 231/01 (Italian law for crime-risk prevention within business activities) and undertakes to conduct, execute and monitor any activity in accordance with its new Ethical Code

Reference standard	Certificate Number	Validity – Expire date
ISO 13485:2016 MDSAP	QS6 049539 0016 Rev. 00	2022-05-17
EN ISO 13485:2016	Q5 049539 0013 Rev. 01	2022-05-07
UNI EN ISO 9001:2015	50 100 5099 Rev. 008	2022-04-16
EC Certificate (according to Directive 98/79/EC, Annex IV excluding (4, 6)	V1 049539 0014 Rev. 01	2023-08-12

Corporate Social Responsibility

DIESSE Diagnostica Senese S.p.A. (DIESSE) is committed to conducting, executing, and monitoring its activities in an ethical, socially responsible and transparent way to safeguard the respect and preservation of the environment, as well as cultural and social heritage, in order to protect this generation and those to come.

The principles that DIESSE intends to follow are described in the Corporate Social Responsibility policy and are in agreement with its Code of Ethics and its Organization, Management and Control Model (in compliance with the Italian Legislative Decree n. 231/2001, Italian law on crime-risk prevention within business activities).

DIESSE wants to express the precise intention to inspire its behaviors, not only to respect laws and regulations, but also to respect moral principles, which regulate all the aspects of civil coexistence, all relationships between individuals, private and public companies, and institutions.

The goals that DIESSE, with this policy, aims to achieve are taken from the 17 United Nation Sustainable Development Goals.

In particular, DIESSE focuses its interest on five goals:



This new culture and civic sense responsibility will be integrated into the Quality Management System, with the aim of monitoring all the company processes by promoting continuous improvement.

It is on this foundation that DIESSE builds and develops its strategies and competitiveness on the market.

GENERAL PRODUCT INDEX

CODE	DESCRIPTION	PAGE
10200	VACU-TEC	13
10200/W	VACU-TEC DOUBLE LABEL	13
10201/A	VES-TEC	13
10208	VACU-TEC CUVETTE-L	13
10214	VES-TEC CUVETTE-L	13
10230	VACU-CODE	13
10235	VACU-TEC FLAG	13
10290	TEST DEVICE 10K	11
10291	TEST DEVICE 5K	11
10292	TEST DEVICE 1K	11
10293	TEST DEVICE NEXT 500	9
10294	TEST DEVICE NEXT 1K	9
10296	TEST DEVICE NEXT 5K	9
10297	TEST DEVICE NEXT 10K	9
10370/A	VES-MATIC CUBE 200 for Advia Bayer	10
10370/AB	VES-MATIC CUBE 200 for ABX	10
10370/BC	VES-MATIC CUBE 200 for Beckman Coulter	10
10370/BC5	VES-MATIC CUBE 200 PRED.BECK.COULTER 5 P	10
10370/S	VES-MATIC CUBE 200 for Sysmex	10
10376	VES-MATIC EASY	12
10380	VES-MATIC CUBE 80	10
10385	VES-MATIC CUBE TRACK	10
10392	MINI-CUBE	8
10395	CUBE 30 TOUCH	8
10403	THERMAL PAPER	9,11,13
10404	MINI-CUBE THERMAL PAPER	9
10430	ESR CONTROL	12
10434	ESR CONTROL	12
10435	ESR CONTROL CUBE	9,11
10436	ESR CONTROL CUBE	9,11
10600	VACU-TEC DS.C.	13
10602	VACU-TEC S.C. STERILE	13

CODE	DESCRIPTION	PAGE
26000	AUTO-DAT*	38
26005	AUTO-DAT MICROPLATE*	39
26010	WW-DAT*	39
26020	WF-DAT*	39
26030	RPR-DAT*	39
26035	TPHA DAT*	39
81010	CHORUS CYTOMEGALOVIRUS IgG	17
81011	CHORUS CYTOMEGALOVIRUS IgM	17
81012	CHORUS CYTOMEGALOVIRUS IgM CAPTURE (NEW)	17
81018	CHORUS HSV 1 SCREEN	17
81019	CHORUS HSV 2 SCREEN	18,22
81020	CHORUS HERPES SIMPLEX 1+2 IgG	17
81021	CHORUS HERPES SIMPLEX 1+2 IgM	17
81023	CHORUS HERPES SIMPLEX 1 IgG RECOMBINANT	17
81029	CHORUS HERPES SIMPLEX 2 IgG RECOMBINANT	17,22
81030	CHORUS RUBELLA IgG	18
81031	CHORUS RUBELLA IgM	18
81033	CHORUS MYCOPLASMA PNEUMONIAE IgA	21
81034	CHORUS MYCOPLASMA PNEUMONIAE IgG	21
81035	CHORUS MYCOPLASMA PNEUMONIAE IgM	21
81036	CHORUS RESPIRATORY SYNCYTIAL VIRUS IgG	22
81036/12	CHORUS RESPIRATORY SYNCYTIAL VIRUS IgG	22
81037	CHORUS RESPIRATORY SYNCYTIAL VIRUS IgA	22
81037/12	CHORUS RESPIRATORY SYNCYTIAL VIRUS IgA	22
81040	CHORUS TOXOPLASMA IgG	18
81041	CHORUS TOXOPLASMA IgM	18
81042	CHORUS TOXOPLASMA IgM CAPTURE (NEW)	18
81043	CHORUS TOXOPLASMA IgA	18
81043/12	CHORUS TOXOPLASMA IgA	18
81044	CHORUS TOXOPLASMA IgA CAPTURE (NEW)	18
81046	CHORUS BORRELIA IgG	24
81047	CHORUS BORRELIA IgM	24

GENERAL PRODUCT INDEX

CODE	DESCRIPTION	PAGE
81048	CHORUS BORRELIA IgG CSF	24
81050	CHORUS TREPONEMA IgG	22
81050/12	CHORUS TREPONEMA IgG	22
81051	CHORUS TREPONEMA IgM	23
81051/12	CHORUS TREPONEMA IgM	23
81054	CHORUS EPSTEIN BARR VCA IgM II	23
81055	CHORUS EPSTEIN BARR VCA IgG	23
81057	CHORUS EPSTEIN BARR EBNA IgG	23
81058	CHORUS EPSTEIN BARR EARLY ANTIGEN IgG	23
81058/12	CHORUS EPSTEIN BARR EARLY ANTIGEN IgG	23
81059	CHORUS EPSTEIN BARR EARLY ANTIGEN IgM	23
81059/12	CHORUS EPSTEIN BARR EARLY ANTIGEN IgM	23
81060	CHORUS HELICOBACTER PYLORI IgG	23
81062	CHORUS HELICOBACTER PYLORI IgA	23
81063	CHORUS HELICOBACTER PYLORI Ag	26
81067	CHORUS TOSCANA VIRUS IgG	24
81067/12	CHORUS TOSCANA VIRUS IgG	24
81068	CHORUS TOSCANA VIRUS IgM	24
81068/12	CHORUS TOSCANA VIRUS IgM	24
81072	CHORUS MEASLES IgG	19
81073	CHORUS MEASLES IgM	19
81073/12	CHORUS MEASLES IgM	19
81075	CHORUS MUMPS IgG	19
81076	CHORUS MUMPS IgM	19
81076/12	CHORUS MUMPS IgM	19
81078	CHORUS VARICELLA IgG	19
81079	CHORUS VARICELLA IgM	19
81079/12	CHORUS VARICELLA IgM	19
81090	CHORUS LEGIONELLA PNEUMOPHILA 1-6 IgG	20
81090/12	CHORUS LEGIONELLA PNEUMOPHILA 1-6 IgG	20
81091	CHORUS LEGIONELLA PNEUMOPHILA 1 IgG	20
81091/12	CHORUS LEGIONELLA PNEUMOPHILA 1 IgG	20

CODE	DESCRIPTION	PAGE
81092	CHORUS CYTOMEGALOVIRUS IgG AVIDITY	17
81093	CHORUS LEGIONELLA PNEUMOPHILA IgM	20
81093/12	CHORUS LEGIONELLA PNEUMOPHILA IgM	20
81095	CHORUS RUBELLA IgG AVIDITY	18
81098	CHORUS TOXOPLASMA IgG AVIDITY	18
81100	CHORUS SYPHILIS SCREEN RECOMBINANT	22
81133	CHORUS CHIKUNGUNYA IgG	24
81134	CHORUS CHIKUNGUNYA IgM	24
81150	CHORUS PARVOVIRUS B19 IgG	19
81152	CHORUS PARVOVIRUS B19 IgM	19
81154	CHORUS DENGUE IgG	23
81156	CHORUS DENGUE IgM	24
81164	CHORUS BRUCELLA IgG*	24
81166	CHORUS BRUCELLA IgM*	25
81168	CHORUS CLOSTRIDIUM DIFFICILE GDH	26
81170	CHORUS CLOSTRIDIUM DIFFICILE A/B TOXINS	26
81171	CHORUS Q FEVER FASE II IgG (NEW)	21,25
81172	CHORUS Q FEVER FASE II IgM (NEW)	21,25
81190	CHORUS INFLUENZA A IgG	21
81190/12	CHORUS INFLUENZA A IgG	21
81192	CHORUS INFLUENZA A IgA	21
81192/12	CHORUS INFLUENZA A IgA	21
81193	CHORUS INFLUENZA B IgG	22
81193/12	CHORUS INFLUENZA B IgG	22
81195	CHORUS INFLUENZA B IgA	21
81195/12	CHORUS INFLUENZA B IgA	21
81196	CHORUS ADENOVIRUS IgG	21
81196/12	CHORUS ADENOVIRUS IgG	21
81198	CHORUS ADENOVIRUS IgA	21
81198/12	CHORUS ADENOVIRUS IgA	21
81200	CHORUS TRIO	14
81210	CHORUS AUTOSAMPLER	14

GENERAL PRODUCT INDEX

CODE	DESCRIPTION	PAGE
81215	AUTOSAMPLER UPGRADE KIT	14
81220	CHORUS TRIO RAS	14
81250	CHORUS CHLAMYDOPHILA PNEUMONIAE IgG	20
81251	CHORUS CHLAMYDOPHILA PNEUMONIAE IgM	20
81252	CHORUS CHLAMYDOPHILA PNEUMONIAE IgA	20
81261	CHORUS CHLAMYDIA TRACHOMATIS IgG	22
81262	CHORUS CHLAMYDIA TRACHOMATIS IgA	22
81264	CHORUS POLIOVIRUS IgG	19
81264/12	CHORUS POLIOVIRUS IgG	19
81266	CHORUS DIPHTHERIA IgG	19
81266/12	CHORUS DIPHTHERIA IgG	19
81270	CHORUS TETANUS IgG	20
81270/12	CHORUS TETANUS IgG	20
81272	CHORUS BORDETELLA PERTUSSIS IgG (NEW)	20
81274	CHORUS BORDETELLA PERTUSSIS IgA (NEW)	20
81276	CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgG	25
81277	CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgM	25
81280	CHORUS TOXOCARA IgG	25
81284	CHORUS ECHINOCOCCUS IgG (NEW)	25
81288	CHORUS LEPTOSPIRA IgG*	25
81289	CHORUS LEPTOSPIRA IgM*	25
81301	CHORUS LEGIONELLA URINARY ANTIGEN	25
83604	SANITIZING SOLUTION	16
83606	WASHING BUFFER	16
83607	CHORUS NEGATIVE CONTROL/SAMPLE DILUENT	16
83609	CLEANING SOLUTION 2000	16
83610	THERMAL PAPER	15
83611	TIPS FOR CHORUS AUTOSAMPLER	15
85221	CHORUS ECHOVIRUS P MIX	33
85225	CHORUS PARAINFLUENZA MIX	33
85251	CHORUS COXSACKIE B MIX	33
85275	CHORUS ECHO VIRUS N MIX	33

CODE	DESCRIPTION	PAGE
85313	CHORUS COXSACKIE A MIX	33
85330	CHORUS Q-FEVER	33
85344	CHORUS LISTERIA MONOCYTOGENES	33
85358	CHORUS BRUCELLA	33
85367	CHORUS LEPTOSPIRA MIX	33
85875	CHORUS LISTERIA MONOCYTOGENES CONTROL SERUM	33
85886	CHORUS Q-FEVER CONTROL SERUM	33
85893	CHORUS PARAINFLUENZA CONTROL SERUM	33
85896	CHORUS LEPTOSPIRA CONTROL SERUM	33
85903	CHORUS COXSACKIE B CONTROL SERUM	33
86004	WASHING BUFFER AUTOIMMUNITY	16
86010	CHORUS ANA-8	26
86010/12	CHORUS ANA-8	26
86012	CHORUS ENA-6 S	27
86012/12	CHORUS ENA-6 S	27
86014	CHORUS ANA-Screen	27
86014/12	CHORUS ANA-Screen	27
86016	CHORUS Sm	28
86016/12	CHORUS Sm	28
86018	CHORUS SS-A	28
86018/12	CHORUS SS-A	28
86024	CHORUS SS-B	28
86024/12	CHORUS SS-B	28
86026	CHORUS Scl-70	28
86026/12	CHORUS Scl-70	28
86028	CHORUS Cenp-B	27
86028/12	CHORUS Cenp-B	27
86030	CHORUS Jo-1	27
86030/12	CHORUS Jo-1	27
86032	CHORUS dsDNA-G	27
86032/12	CHORUS dsDNA-G	27
86034	CHORUS dsDNA-M	27

GENERAL PRODUCT INDEX

CODE	DESCRIPTION	PAGE
86034/12	CHORUS dsDNA-M	27
86038	CHORUS RF-G	27
86038/12	CHORUS RF-G	27
86040	CHORUS RF-M	28
86040/12	CHORUS RF-M	28
86046	CHORUS CARDIOLIPIN-G	29
86046/12	CHORUS CARDIOLIPIN-G	29
86048	CHORUS CARDIOLIPIN-M	29
86048/12	CHORUS CARDIOLIPIN-M	29
86050	CHORUS BETA 2-GLYCOPROTEIN-G	29
86050/12	CHORUS BETA 2-GLYCOPROTEIN-G	29
86052	CHORUS BETA 2-GLYCOPROTEIN-M	29
86052/12	CHORUS BETA 2-GLYCOPROTEIN-M	29
86054	CHORUS GLIADIN-A	31
86054/12	CHORUS GLIADIN-A	31
86056	CHORUS GLIADIN-G	31
86056/12	CHORUS GLIADIN-G	31
86058	CHORUS tTg-A	31
86058/12	CHORUS tTg-A	31
86060	CHORUS tTg-G	31
86060/12	CHORUS tTg-G	31
86062	CHORUS ASCA-A	30
86062/12	CHORUS ASCA-A	30
86064	CHORUS ASCA-G	30
86064/12	CHORUS ASCA-G	30
86066	CHORUS PR3	30
86066/12	CHORUS PR3	30
86068	CHORUS MPO	30
86068/12	CHORUS MPO	30
86070	CHORUS GBM	30
86070/12	CHORUS GBM	30
86072	CHORUS a-TG	28

CODE	DESCRIPTION	PAGE
86072/12	CHORUS a-TG	28
86074	CHORUS a-TPO	28
86074/12	CHORUS a-TPO	28
86076	CHORUS TG	32
86076/12	CHORUS TG	32
86080	CHORUS LKM-1	29
86080/12	CHORUS LKM-1	29
86084	CHORUS AMA-M2	29
86084/12	CHORUS AMA-M2	29
86086	CHORUS INSULIN	30
86086/12	CHORUS INSULIN	30
86088	CHORUS U1-70 RNP	28
86088/12	CHORUS U1-70 RNP	28
86090	CHORUS INTRINSIC FACTOR	31
86090/12	CHORUS INTRINSIC FACTOR	31
86092	CHORUS snRNP-C	28
86092/12	CHORUS snRNP-C	28
86094	CHORUS ANTI-CCP	27
86094/12	CHORUS ANTI-CCP	27
86098	CHORUS DEAMIDATED GLIADIN-A	30
86098/12	CHORUS DEAMIDATED GLIADIN-A	30
86100	CHORUS DEAMIDATED GLIADIN-G	30
86100/12	CHORUS DEAMIDATED GLIADIN-G	30
86102	CHORUS CH50*	27
86104	CHORUS Tg HIGH SENSITIVITY (NEW)	27
86106	CHORUS PHOSPHOLIPID SCREEN IgG*	29
86108	CHORUS PHOSPHOLIPID SCREEN IgM*	29
86510	CHORUS FSH (NEW)	32
86512	CHORUS LH (NEW)	32
86514	CHORUS PRL (NEW)	32
86516	CHORUS TSH (NEW)	31
86518	CHORUS FT3 (NEW)	31

GENERAL PRODUCT INDEX

CODE	DESCRIPTION	PAGE
86520	CHORUS FT4 (NEW)	31
86550	CHORUS PSA (NEW)	32
86600	CHORUS CALPROTECTIN	26
86600/C03	CHORUS CALPROTECTIN	26
86601	CHORUS CALPROTECTIN K	26
86602	CHORUS COPROCOLLECT	26
86900	CHORUS 25OH VITAMIN D TOTAL	26
88000	CHORUS TOTAL IgE	32
91010	CYTOMEGALOVIRUS IgG	34
91011	CYTOMEGALOVIRUS IgM	34
91013	CMV SCREEN	34
91020	HERPES SIMPLEX 1+2 IgG	34
91021	HERPES SIMPLEX 1+2 IgM	34
91023	HERPES SIMPLEX 1 IgG RECOMBINANT	34
91025	HERPES SIMPLEX 2 IgM	35,36
91026	HERPES SIMPLEX 2 IgG	34,36
91027	HERPES SIMPLEX 1 IgG	34
91028	HERPES SIMPLEX 1 IgM	34
91029	HERPES SIMPLEX 2 IgG RECOMBINANT	34,36
91030	RUBELLA IgG	35
91031	RUBELLA IgM	35
91040	TOXOPLASMA IgG	35
91041	TOXOPLASMA IgM	35
91043	TOXOPLASMA IgA	35
91050	TREPONEMA IgG	36
91051	TREPONEMA IgM	36
91055	EPSTEIN-BARR VCA IgG	37
91056	EPSTEIN-BARR VCA IgM	37
91057	EPSTEIN-BARR EBNA IgG	37
91060	HELICOBACTER PYLORI IgG	37
91062	HELICOBACTER PYLORI IgA	37
91067	TOSCANA VIRUS IgG-IgM	37

CODE	DESCRIPTION	PAGE
91072	MEASLES IgG	35
91073	MEASLES IgM	35
91075	MUMPS IgG	35
91076	MUMPS IgM	36
91078	VARICELLA IgG	36
91079	VARICELLA IgM	36
91092	CYTOMEGALOVIRUS IgG AVIDITY	34
91095	RUBELLA IgG AVIDITY	35
91098	TOXOPLASMA IgG AVIDITY	35
91100	SYPHILIS SCREEN RECOMBINANT	36
91133	CHIKUNGUNYA IgG	37
91134	CHIKUNGUNYA IgM	37
P20550510	BAR-CODE READER	9,11,12
R21430320	BLUETOOTH MINI-CUBE PRINTER	9

Valid as
january 1st 2020



www.diesse.it