

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

## QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS PROMONITOR INFLIXIMAB

Code 86700

Lot No. 861



2023-08



2022-08

### REAGENTS

Label	Description	Lot.	Exp.	NO
D4311	PROMONITOR INFLIXIMAB	861	2023-08	6x6
B4312	CALIBRATORE IFX	861	2024-04	1
C4311	CONTROL SERUM IFX	861	2024-04	1

### CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	345	$\geq 100$
Control Serum (ug/L)	2,1	2,0-3,5

### CALIBRATOR


2	7	4	8	6	1	6	3	2	4	4	2	4	4	5	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature



### CONTROL SERUM

2	7	4	8	6	1	6	1	9	2	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.

DATE 03.10.2022

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.