

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600/C03

Lot No. 878



2023-10



2022-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	878	2023-10	2X6
D9510	CALPROTECTIN	878	2023-10	1X4
B9510	CALIBRATORE CALPROTECTIN	878	2023-11	1
C9510	CONTROL CALPROTECTIN	878	2023-11	1
G9502	SAMPLE DILUENT	878	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	440	≥ 100
Control (mg/kg)	151,3	81,3-225,0

CALIBRATOR

9 0 0 8 7 8 6 2 1 8 2 1 4 8 2 9 0 0

Label checked with barcode reader

YES ☒

Operator Signature



CONTROL SERUM

9 0 0 8 7 8 6 1 5 2 7 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 18.11.2022

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.