

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

## QUALITY CONTROL CERTIFICATE

**KIT NAME** CHORUS PROMONITOR ADALIMUMAB

**Code** 86702

**Lot No.** 887



2023-11



2022-11

### REAGENTS

Label	Description	Lot.	Exp.	NO
D4331	PROMONITOR ADALIMUMAB	887	2023-11	6x6
B4332	CALIBRATORE ADL	887	2023-11	1
C4331	CONTROL SERUM ADL	887	2023-11	1

### CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	307	$\geq 100$
Control Serum (ug/L)	5,3	4,2-6,3

### CALIBRATOR

2	7	6	8	8	7	6	5	2	2	4	1	2	5	8	8	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *tel*



### CONTROL SERUM

2	7	6	8	8	7	6	3	2	4	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *te*



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *te*

DATE *23-11-2022*

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.