

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS FOLATE

Code 86402

Lot No. 902

2024-01

2023-01

REAGENTS

Label	Description	Lot.	Exp	No
D4560	FOLATE	902	2024-01	6X6
B4560	Calibratore FOLATE	902	2025-06	1
C4560	Control FOLATE	902	2025-06	1
G4560	FOLATE Sample Diluent/Negative Control	902	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	477	≥ 100
Control (ng/ml)	5,4	3,9-6,7
FOLATE Sample Diluent/Negative Control (ng/ml)	<2,0	< 2,0

CALIBRATOR

9	9	1	9	0	2	6	5	1	4	2	1	4	8	3	2	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature



CONTROL SERUM

9	9	1	9	0	2	6	1	7	2	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 31/01/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.