

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME SARS-CoV-2 "Neutralizing" Ab - CHORUS

Code	81408		
Lot No.	887	2023-11	2023-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D4240	SARS-CoV-2 "Neutralizing" Ab	887	2023-11	6X6
B4240	CALIBRATORE SARS-CoV-2 "Neutralizing" Ab	887	2024-01	1
C4240	CONTROL SERUM SARS-CoV-2 "Neutralizing" Ab	887	2024-01	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	590	≥ 100
Control Serum (BAU/mL)	201,5	100,0-265,6

CALIBRATOR

1	7	0	8	8	7	6	2	1	5	2	1	9	8	2	6	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



CONTROL SERUM

1	7	0	8	8	7	6	0	7	1	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 09.02.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.