

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600

Lot No. 899



2023-12



2022-12

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	899	2023-12	2 X 6
B9500	CALIBRATORE CALPROTECTIN	899	2024-05	1
C9500	CONTROL CALPROTECTIN	899	2024-05	1
G9502	SAMPLE DILUENT	899	2024-12	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	390	≥ 100
Control (mg/kg)	103,8	68,8-150,0


CALIBRATOR

9 0 0 8 9 9 6 3 1 8 2 3 5 6 3 3 0 0

Label checked with barcode reader

YES ☒

Operator Signature *RB*




CONTROL SERUM

9 0 0 8 9 9 6 1 3 2 3 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature *RB*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....*Blue*

DATE *15.02.2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent reagents should be handled with maximum attention.