

DIESSE DIAGNOSTICA SENESE SpA
 53035 Monteriggioni (SI) – Italia – Via delle Rose 10
 Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione
 20144 Milano (Italy) Via A. Solari 19 sc. 6
 Tel. +39-02-4859121 r.a. Fax +39-02-48008530

PRODUCT NAME: VACU-TEC

Code 10208		
Lot No 3068	2024-01	2023-03
Reference IO-09/88		

ON LINE QUALITY CONTROLS (model CDP not enclosed)

QC Passed

- sodium citrate solution concentration ☒
- sodium citrate volume ☒
- vacuum level (*) ☒
- correct placement of the cap ☒
- no extra body inside the product ☒

(*) Only for VACU-CODE and VACU-TEC tubes ☒

PACKAGING QUALITY CONTROL (model CDP not enclosed)

QC Passed

- control of the number of the product for each pack ☒

QUALITY CONTROLS ON FINISHED PRODUCT (model CDP not enclosed)

QC Passed

- pack integrity ☒
- data printed on the product label ☒

QC Passed **N/A**

Irradiation treatment Certificate n° date

☐ ☒

QC Passed **N/A**

Gamma rays treatment process in compliance with the specifications of the supply Agreement

☐ ☒

Document number: ODV N.A.

*This batch has been manufactured and checked according to quality requirements.
 Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.

DATE 15/03/2023