



DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME **HELICOBACTER PYLORI IgA-CHORUS**

Code **81062**

Lot No. **917**  **2024-03**  **2023-03**

REAGENTS

Label	Description	Lot.	Exp.	NO
D2401	Helicobacter Pylori IgA	917	2024-03	6X6
B2401	CALIBRATORE Helicobacter Pylori IgA	917	2024-03	1
C2401	CONTROL SERUM Helicobacter Pylori IgA	917	2024-03	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	289	≥ 100
Control Serum (INDEX)	4,6	2,5-4,6
Negative Control (INDEX)	0,0	$< 0,7$

CALIBRATOR

0 2 7 9 1 7 6 0 0 0 0 0 2 9 0 0 0 0

Label checked with barcode reader YES ☒ Operator Signature 



CONTROL SERUM

0 2 7 9 1 7 6 2 5 4 6 0 0 0 0 0 0 1

Label checked with barcode reader YES ☒ Operator Signature 



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. 

DATE **24.03.2023**

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent reagents should be handled with maximum attention.