

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600

Lot No. 919-A



2024-03



2023-03

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	919	2024-03	2 X 6
B9500	CALIBRATORE CALPROTECTIN	919-1	2024-08	1
C9500	CONTROL CALPROTECTIN	919-1	2024-08	1
G9502	SAMPLE DILUENT	919	2025-03	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	360	> 100
Control (mg/kg)	142,0	80,0-216,7

CALIBRATOR

9	0	0	9	1	9	6	2	2	0	2	4	4	4	3	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature *Green*




CONTROL SERUM

9	0	0	9	1	9	6	1	6	3	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature *Green*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Green*

DATE 20-06-2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.