

QUALITY CONTROL CERTIFICATE

KIT NAME **HELICOBACTER PYLORI Ag - CHORUS**

Code **81063**

Lot No. 926

2024-04

2023-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D6960	Helicobacter pylori Ag	926	2024-04	2X6
B6960	CALIBRATORE Helicobacter pylori Ag	926	2024-07	2
C6960	CONTROL SERUM Helicobacter pylori Ag	926	2024-07	3
G6960	H. pylori Ag Stool Diluent/Negative Control	C0553500-1	2024-04	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	270	≥ 100
Control + (Index)	4,1	2,6-4,9
Stool Diluent/Negative Control (Index)	0,2	$< 0,5$

CALIBRATOR

9	3	2	9	2	6	6	0	0	0	0	0	2	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES ☒ Operator Signature *[Signature]*



CONTROL SERUM

9	3	2	9	2	6	6	2	6	4	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature *[Signature]*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*

DATE *21-06-2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.