

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

### QUALITY CONTROL CERTIFICATE

KIT NAME TOSCANA VIRUS IgG-CHORUS

|         |       |  |         |  |         |
|---------|-------|--|---------|--|---------|
| Code    | 81067 |  |         |  |         |
| Lot No. | 931   |  | 2024-04 |  | 2023-04 |

### REAGENTS

| Label | Description                     | Lot. | Exp.    | NO  |
|-------|---------------------------------|------|---------|-----|
| D2411 | TOSCANA VIRUS IgG               | 931  | 2024-04 | 6X6 |
| B2412 | CALIBRATORE TOSCANA VIRUS IgG   | 931  | 2024-04 | 1   |
| C2411 | CONTROL SERUM TOSCANA VIRUS IgG | 931  | 2024-04 | 1   |

### CONTROLS


|  | Found results | Expected results |
|--|---------------|------------------|
|--|---------------|------------------|

|                          |     |            |
|--------------------------|-----|------------|
| Calibrator (mOD)         | 690 | $\geq 100$ |
| Control Serum (INDEX)    | 1,6 | 1,2-2,3    |
| Negative Control (INDEX) | 0,0 | $< 0,7$    |

### CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 1 | 3 | 9 | 3 | 1 | 6 | 0 | 0 | 0 | 0 | 0 | 5 | 7 | 0 | 0 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader YES ☒ Operator Signature



### CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 1 | 3 | 9 | 3 | 1 | 6 | 1 | 2 | 2 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader YES ☒ Operator Signature



This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 18/05/2023

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.