

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS FOLATE

Code 86402

Lot No. 926

2024-04

2023-04

REAGENTS

Label	Description	Lot.	Exp	No
D4560	FOLATE	926	2024-04	6X6
B4560	Calibratore FOLATE	926	2026-03	1
C4560	Control FOLATE	926	2026-03	1
G4560	FOLATE Sample Diluent/Negative Control	926	2026-03	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	5,78	> 100
Control (ng/ml)	4,7	2,7-6,5
FOLATE Sample Diluent/Negative Control (ng/ml)	<2,0	< 2,0

CALIBRATOR

9	9	1	9	2	6	6	6	1	4	1	8	5	1	2	6	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature



CONTROL SERUM

9	9	1	9	2	6	6	1	6	2	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 18-05-2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.