

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

**DIESSE**

Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME SYPHILIS SCREEN RECOMBINANT - CHORUS

|         |         |
|---------|---------|
| Code    | 81100   |
| Lot No. | 948     |
|         | 2024-06 |
|         | 2023-07 |

REAGENTS

| Label | Description                               | Lot. | Exp.    | NO  |
|-------|-------------------------------------------|------|---------|-----|
| D2541 | Syphilis Screen Recombinant               | 948  | 2024-06 | 6X6 |
| B2541 | CALIBRATORE Syphilis Screen Recombinant   | 948  | 2024-06 | 1   |
| C2541 | CONTROL SERUM Syphilis Screen Recombinant | 948  | 2024-06 | 1   |

CONTROLS

|  | Found results | Expected results |
|--|---------------|------------------|
|--|---------------|------------------|

|                          |     |            |
|--------------------------|-----|------------|
| Calibrator (mOD)         | 450 | $\geq 100$ |
| Control Serum (INDEX)    | 2,3 | 1,7-3,2    |
| Negative Control (INDEX) | 0,2 | $< 0,7$    |

CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 3 | 2 | 9 | 4 | 8 | 6 | 0 | 0 | 0 | 0 | 0 | 4 | 1 | 0 | 0 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader YES ☒ Operator Signature *RB*



CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 0 | 3 | 2 | 9 | 4 | 8 | 6 | 1 | 7 | 3 | 2 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader YES ☒ Operator Signature *RB*



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE 04/09/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.