

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME Cemp-B - CHORUS

Code 86028

Lot No. 962



2024-08



2023-08

REAGENTS

Label	Description	Lot.	Exp.	NO
D3241	Cemp-B	962	2024-08	6X6
B3241	CALIBRATORE Cemp-B	962	2024-08	1
C3241	CONTROL SERUM Cemp-B	962	2024-08	1


CONTROLS

	Found results	Expected results
Calibrator (mOD)	337	≥ 100
Control Serum (AU/mL)	56,4	39,2-89,3

CALIBRATOR

1	9	6	9	6	2	6	2	1	9	3	2	3	0	4	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



CONTROL SERUM

1	9	6	9	6	2	6	2	3	4	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 05/09/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.