

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME LEGIONELLA PNEUMOPHILA 1-6 IgG -CHORUS

Code 81090

Lot No. 978



2024-07



2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D7381	LEGIONELLA PNEUMOPHILA 1-6 IgG	978	2024-07	6X6
B7381	CALIBRATORE L. PNEUMOPHILA 1-6 IgG	978	2024-07	1
C7381	CONTROL SERUM L. PNEUMOPHILA 1-6 IgG	978	2024-07	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	411	≥ 100
Control Serum (INDEX)	3,7	3,1-5,7
Negative Control (INDEX)	0,0	$< 0,7$

CALIBRATOR


0	6	3	9	7	8	5	0	0	0	0	0	2	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *osm*



CONTROL SERUM


0	6	3	9	7	8	5	3	1	5	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *osm*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *osm*

DATE 02.11.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.