

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME ENZYWELL HERPES SIMPLEX 1 + 2 IgM

Code 91021

LOT.NO 326



2025-01



2023-10

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
PF93000	006190	MICROPLATE	672	2025-04	1
PF91520	006180	CONJUGATE	246	2025-01	1
PF93900/HSV	C4940	NEGATIVE CONTROL	804	2025-04	1
PF92022	005140	POSITIVE CONTROL	225	2025-04	1
PF91821	006150	CUT OFF CONTROL	238	2025-04	1
PF91621	006170	ANTIGEN	223	2026-03	6
PF93611	009660	DILUENT 2	231	2025-03	1
PF93603	1035	WASHING BUFFER 10x	464	2025-04	1
PF93602	1032/RISK0	STOP SOLUTION	218-5	2025-06	1
PF93619	1034	SUBSTRATE-HS	380	2025-03	1

CONTROLS

	FOUND RESULTS	EXPECTED RESULTS
NEGATIVE CONTROL	0.086	< 0,6 Cut Off
CUT OFF CONTROL	0.428	≥ 0.200
POSITIVE CONTROL	1.979	> 1,5 Cut Off
CV %	10%	≤ 15%

*Evaluation of the method according to the package inset in the attached document
This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents

DATE: 03.11.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.