



DIESE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME MYCOPLASMA PNEUMONIAE IgM - CHORUS

Code81035

Lot No.011

2025-05

2024-02

REAGENTS

Label	Description	Lot.	Exp.	NO
D2261	MYCOPLASMA PNEUMONIAE IgM	011	2025-05	6x6
B2261	CALIBRATORE MYCOPLASMA PNEUMONIAE IgM	011	2025-05	1
C2261	CONTROL SERUM MYCOPLASMA PNEUMONIAE IgM	011	2025-05	1


CONTROLS

	Found results	Expected results
Calibrator (mOD)	452	≥ 100
Control Serum (INDEX)	1,5	1,3-2,3
Negative Control (INDEX)	0,1	≤ 0,7

CALIBRATOR

067011700000560000

Label checked with barcode readerYESOperator Signature



CONTROL SERUM

067011713230000001

Label checked with barcode readerYESOperator Signature



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 06/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma.All serum or plasma used has been tested and found non-reactive