

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME ENZYWELL EPSTEIN BARR EBNA IgG

Code 91057

LOT.NO 106

 2025-06

 2024-04

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
PF93071	054800	MICROPLATE	121	2025-09	1
PF93509	005470	CONJUGATE	106	2025-07	1
PF93910	1036	NEGATIVE CONTROL	179	2025-06	1
PF92073	005390	POSITIVE CONTROL	098	2025-10	1
PF91873	005450	CUT OFF CONTROL	098	2025-10	1
PF93619	1034	SUBSTRATE-HS	385	2025-08	1
PF93621	006030	DILUENT 10	166	2025-08	1
PF93603	1035	WASHING BUFFER 10x	473	2025-08	1
PF93602	1032/RISK0	STOP SOLUTION	220	2026-03	1
PF91083	006840	SORBENT G	096	2025-09	1

CONTROLS

	FOUND RESULTS	EXPECTED RESULTS
NEGATIVE CONTROL	0.000	$\leq 0,6$ Cut Off
CUT OFF CONTROL	0.375	≥ 0.200
POSITIVE CONTROL	2.094	$\geq 1,5$ Cut Off
CV %	6%	$\leq 15\%$

*Evaluatin of the method according to the package inset in the attached document
This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents 

DATE: 03.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.

IO 09- 022-V



CODE 91057

CODE

QUALITY CONTROL CERTIFICATE

Annex: Method Evaluation

NAME KIT EPSTEIN BARR EBNA IgG
LOT. 106
DATE OF TEST 30/04/2024

METHOD: ACCORDING TO THE PACKAGE INSERT

1. SOLID PHASE	FOUND		EXPECTED
A) UNIFORMITY BETWEEN PLATES			
1st PLATE	OD 450 nm	0.466	
	CV%	12%	< 20%
2nd PLATE	OD 450 nm	N/A	
	CV%	N/A	< 20%
B) UNIFORMITY WITHIN PLATE			
	OD 450 nm	0.433	
	CV%	10%	< 15%
2. BULK CONTROL date: 29/04/2024			
A) Average OD of Cut Off Control	OD 450 nm	0.365	≥ 0.200
B) Index of Negative and Positive Sera		accepted	According to specification
C) Sensitivity with working calibrators		accepted	According to specification
3. FINAL CONTROL date: 30/04/2024			
A) Blank	OD 450 nm	0.052	≤ 0.100
B) Average OD of Negative Control	OD 450 nm	0.000	≤ 0,6 Cut Off
C) Average OD of Cut Off Control	OD 450 nm	0.375	≥ 0.200
D) Average OD of Positive Control	OD 450 nm	2.094	≥ 1,5 Cut Off
E) HOMOGENEITY OF CUT OFF			
	No.	6	
	OD 450 nm	0.427	
	CV%	6%	≤ 15%
F) PRECISION OF FILLING PROCESS			



SATISFACTORY



UNSATISFACTORY

FINAL RESULT: ACCEPTED

Date: 30/04/2024

Signed:

Released: Q.C.Manager