

**AUTO-DAT
Microplate**



REF 26005

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy





ISTRUZIONI PER L'USO

AUTO-DAT Microplate

Micropiastre da utilizzare per le reazioni di agglutinazione sullo strumento AUTO-DAT (REF 26000).

Solo per uso diagnostico in vitro

1. UTILIZZAZIONE

Micropiastre da utilizzare per le reazioni di agglutinazione sullo strumento AUTO-DAT.

2. COMPOSIZIONE DEL KIT

20 micropiastre da 96 tests a fondo piatto con frazionamento 2x4 in ogni strip da 8 pozzetti.

3. AVVERTENZE PER L'USO

Per bloccare le strip sulla cornice premere al centro delle stesse.



DIESE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy



INSTRUCTIONS FOR USE

AUTO-DAT Microplate

Microplates to be used for agglutination reactions on the AUTO-DAT instrument (REF 26000).

For in vitro diagnostic use only

1. INTENDED USE

Microplates to be used for agglutination reactions on the AUTO-DAT instrument.

2. KIT COMPOSITION

20 pieces of 96 tests microplates, flat bottom, with 2x4 fractioning in each 8-wells strip.













3. WARNINGS

To lock the strips on the frame, press on their centre.



DIESE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy



	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabricante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote