

**DIESE DIAGNOSTICA SENESE SpA**  
 53035 Monteriggioni (SI) – Italia – Via delle Rose 10  
 Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione  
 20144 Milano (Italy) Via A. Solari 19 sc. 6  
 Tel. +39-02-4859121 r.a. Fax +39-02-48008530

**CERTIFICATO DI CONTROLLO QUALITÀ/QUALITY CONTROL CERTIFICATE**

**PRODOTTO/PRODUCT** ESR CONTROL CUBE

**REF:** 10435 (4x9 ml)  
 10436 (2x9 ml)

**REAGENTI/REAGENTS**

	Kit	Level 1	Level 2
Lot:	633	633	633
	2020-08	2020-08	2020-08
	2022-02	2022-02	2022-02

**RISULTATI/RESULTS**

	Level 1		Level 2		Temp.
	Risultato/ Result	Range	Risultato/ Result	Range	
VES-Matic Cube 30	9	1-18	90	65-115	OFF
VES-Matic Cube 80/VES-Matic Cube 200/Ves-Matic Cube Track	7	1-15	69	44-94	OFF
MINI-CUBE/Cube 30 touch	9	1-18	81	56-106	ON

**Data/Date**

05/10/2020

**Testato da/Tested by:**

Katia Tanzini

**Approvato da/Approved by:**

Alessandra Cianchi (Responsabile CQ/QC Head)