



DIESE DIAGNOSTICA SENESE SpA
Via delle Rose 10, 53035 Monteriggioni (SI), Italy
Tel. ++39-0577-587111 fax ++39-0577-318690

Sede Legale e Amministrazione
Via A. Solari 19 20144 Milano Italy
Tel. ++39-02-4859121 r.a. fax ++39-02-48008530

CERTIFICATE OF ANALYSIS

CODE	PRODUCT NAME	LOT	EXPIRY DATE	PREP DATE
10439	SediCheck Cube 30 4 x 3 ml	024	2021-08	2020-08

Level I (NORMAL Level) lot 635	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	5	1-14
MINI-CUBE	Temperature ON		
Cube 30 Touch	Temperature ON		

Level II (ABNORMAL Level) lot 635	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	60	35-85
MINI-CUBE	Temperature ON		
Cube 30 Touch	Temperature ON		

Date: 23-10-2020

Signature (Head Quality Control)