



CHORUS Aspergillus galactomannan Ag



CHORUS Aspergillus galactomannan Ag

For the **qualitative determination** of Aspergillus galactomannan antigen in human **serum** and bronchial alveolar lavage fluid (**BAL**).

Benefits of galactomannan antingen detection

- **Early diagnosis of invasive aspergillosis** - key factor for effective treatment and a possible reduction in mortality
- **Appropriate drugs treatment** - to avoid the rapid worsening of patients respiratory status
- **Same day results** - no need of batching
- **Improved patient care**

CHORUS advantages



Single test
ready to
use device



Random
access
system



Safe and
easy to use



No hidden costs
No waste of
reagents

Main features

- Calibrator and control included in the kit
- Calibration stability: every new lot
- Time to first result: 2 hours
- Sensitivity: 93.4% ($CI_{95\%}$: 86.3 - 96.9)
- Specificity: 94.4% ($CI_{95\%}$: 91.0 - 97.1)



NOME	TEST	CODICE
CHORUS Aspergillus galactomannan Ag	36	81178