

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) - Italia - Via della Rose 10
Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

QUALITY CONTROL CERTIFICATE

KIT NAME: CALPROTECTIN - CHORUS

Code 86600/C03	
Lot No. 852	2023-07 2022-07

REAGENTS

Label	Description	Lot.	Exp	No
D9500	CALPROTECTIN	852	2023-07	2X6
D9510	CALPROTECTIN	852	2023-07	1X4
B9510	CALIBRATORE CALPROTECTIN	852	2023-11	1
C9510	CONTROL CALPROTECTIN	852	2023-11	1
G9502	SAMPLE DILUENT	852	2024-03	1

CONTROLS

Found results	Expected results
---------------	------------------

Calibrator (mOD)	400	≥ 100
Control (mg/kg)	144,7	78,6-200,0

CALIBRATOR

9	0	0	8	5	2	6	2	1	7	2	1	5	1	3	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

Operator Signature *AP*


CONTROL

9	0	0	8	5	2	6	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

Operator Signature *AP*

The test has been carried out in accordance with the package insert

TESTED BY:

Antonella Leo

APPROVED

RESPONSIBLE HEAD, Quality Control Reagents.

DATE

26.07.2022

WARNING: POTENTIAL BIOHAZARDOUS MATERIALI

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.