

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

## QUALITY CONTROL CERTIFICATE

KIT NAME 88000-D1 DERMATOPHAGOIDES PTERONYSSINUS

Code 88000-D1

Lot No. 870



2023-09



2022-11

### REAGENTS

| Label | Description                                           | Lot. | Exp.    | NO  |
|-------|-------------------------------------------------------|------|---------|-----|
| D9791 | 88000-D1 DERMATOPHAGOIDES PTERONYSSINUS               | 870  | 2023-09 | 2x6 |
| B9790 | CALIBRATORE 88000-D1 DERMATOPHAGOIDES PTERONYSSINUS   | 870  | 2023-10 | 1   |
| C9790 | CONTROL SERUM 88000-D1 DERMATOPHAGOIDES PTERONYSSINUS | 870  | 2023-10 | 1   |

### CONTROLS

|  | Found results | Expected results |
|--|---------------|------------------|
|--|---------------|------------------|

|                      |      |            |
|----------------------|------|------------|
| Calibrator (mOD)     | 270  | $\geq 100$ |
| Control Serum (KU/L) | 14,0 | 6,0-22,8   |

### CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 6 | 0 | 2 | 8 | 7 | 0 | 6 | 5 | 2 | 7 | 5 | 7 | 2 | 2 | 7 | 7 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader

YES



Operator Signature

SB



### CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 6 | 0 | 2 | 8 | 7 | 0 | 6 | 3 | 2 | 6 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader

YES



Operator Signature

SB



This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 08-11-2022

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.