

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME ASCA-G - CHORUS

Code	86064		
Lot No.	889	2023-11	2022-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D3491	ASCA-G	889	2023-11	6X6
B3491	CALIBRATORE ASCA-G	889	2023-11	1
C3491	CONTROL SERUM ASCA-G	889	2023-11	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	423	≥ 100
Control Serum (AU/mL)	111,9	55,0-163,2

CALIBRATOR

5	0	8	8	8	9	6	6	2	2	3	2	3	0	7	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES ☒ Operator Signature



CONTROL SERUM

5	0	8	8	8	9	6	2	3	4	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 01.12.2022

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.