

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

### QUALITY CONTROL CERTIFICATE

KIT NAME **CARDIOLIPIN-G - CHORUS**

Code	86046				
Lot No.	899		2023-12		2022-12

### REAGENTS

Label	Description	Lot.	Exp.	NO
D3311	CARDIOLIPIN-G	899	2023-12	6X6
B3311	CALIBRATORE CARDIOLIPIN-G	899	2023-12	1
C3311	CONTROL SERUM CARDIOLIPIN-G	899	2023-12	1

### CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	280	$\geq 100$
Control Serum (GPLU/mL)	42,7	30,0-54,0

### CALIBRATOR

3	3	5	8	9	9	6	2	2	0	4	0	3	6	6	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES ☒ Operator Signature *FB*



### CONTROL SERUM

3	3	5	8	9	9	6	2	0	3	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature *FB*



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Beno*

DATE *15.02.2023*

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent reagents should be handled with maximum attention.