

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

## QUALITY CONTROL CERTIFICATE

KIT NAME Epstein-Barr VCA IgG - CHORUS

Code	81055
Lot No.	926
	2024-07
	2023-04

## REAGENTS

Label	Description	Lot.	Exp.	NO
D2341	Epstein-Barr VCA IgG	926	2024-07	6X6
B2342	CALIBRATORE Epstein-Barr VCA IgG	926	2024-07	1
C2341	CONTROL SERUM Epstein-Barr VCA IgG	926	2024-07	1

## CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	240	$\geq 100$
Control Serum (INDEX)	3,9	2,4-4,4
Negative Control (INDEX)	0,0	$< 0,7$

## CALIBRATOR

0	0	5	9	2	6	7	0	0	0	0	0	3	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature .....



## CONTROL SERUM

0	0	5	9	2	6	7	2	4	4	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature .....



This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. ....

DATE 26.04.2023

## WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.