

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) – Italia – Via delle Rose 10
Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

PRODUCT NAME: VACU-TEC

Code 10600

Lot No 3136



2024-03



2023-05

Reference IO-09/88

ON LINE QUALITY CONTROLS (model CDP not enclosed)

QC Passed

- sodium citrate solution concentration
- sodium citrate volume
- vacuum level (*)
- correct placement of the cap
- no extra body inside the product



(*) Only for VACU-CODE and VACU-TEC tubes



PACKAGING QUALITY CONTROL (model CDP not enclosed)

QC Passed

- control of the number of the product for each pack



QUALITY CONTROLS ON FINISHED PRODUCT (model CDP not enclosed)

QC Passed

- pack integrity
- data printed on the product label



QC Passed

N/A

Irradiation treatment Certificate n° date



QC Passed

N/A

Gamma rays treatment process in compliance with the specifications of the supply Agreement



Document number:

ODV N.A.

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED

RESPONSIBLE HEAD, Quality Control Reagents

DATE 24/05/2023