

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS PROMONITOR ADALIMUMAB

Code 86702

Lot No. 943



2024-06



2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D4331	PROMONITOR ADALIMUMAB	943	2024-06	6x6
B4332	CALIBRATORE ADL	943	2025-01	1
C4331	CONTROL SERUM ADL	943	2025-01	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	446	≥ 100
Control Serum (ug/L)	4,5	3,6-5,3

CALIBRATOR


2 7 6 9 4 3 6 3 2 2 4 5 2 8 8 2 0 0

Label checked with barcode reader

YES



Operator Signature *PS*



CONTROL SERUM


2 7 6 9 4 3 6 2 9 4 3 0 0 0 0 0 0 1

Label checked with barcode reader

YES



Operator Signature *PS*



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. *PS*

DATE 05/07/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.