

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

## QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

|         |         |
|---------|---------|
| Code    | 86600   |
| Lot No. | 948     |
|         | 2024-06 |
|         | 2023-06 |

### REAGENTS

| Label | Description              | Lot. | Exp.    | NO    |
|-------|--------------------------|------|---------|-------|
| D9500 | CALPROTECTIN             | 948  | 2024-06 | 2 X 6 |
| B9500 | CALIBRATORE CALPROTECTIN | 948  | 2025-01 | 1     |
| C9500 | CONTROL CALPROTECTIN     | 948  | 2025-01 | 1     |
| G9502 | SAMPLE DILUENT           | 948  | 2025-06 | 1     |


### CONTROLS

|                  | Found results | Expected results |
|------------------|---------------|------------------|
| Calibrator (mOD) | 440           | $\geq 100$       |
| Control (mg/kg)  | 154,5         | 91,7-237,5       |

### CALIBRATOR

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 9 | 0 | 0 | 9 | 4 | 8 | 6 | 2 | 1 | 6 | 2 | 0 | 4 | 8 | 2 | 8 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader YES ☒ Operator Signature .....



### CONTROL SERUM

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 9 | 0 | 0 | 9 | 4 | 8 | 6 | 1 | 5 | 2 | 7 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader YES ☒ Operator Signature .....



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 20/07/2023

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.