

QUALITY CONTROL CERTIFICATE

KIT NAME LEGIONELLA PNEUMOPHILA 1-6 IgG -CHORUS

Code 81090

Lot No. 952



2024-04



2023-07

REAGENTS

Label	Description	Lot.	Exp.	NO
D7381	LEGIONELLA PNEUMOPHILA 1-6 IgG	952	2024-04	6X6
B7381	CALIBRATORE L. PNEUMOPHILA 1-6 IgG	952	2024-04	1
C7381	CONTROL SERUM L. PNEUMOPHILA 1-6 IgG	952	2024-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	378	> 100
Control Serum (INDEX)	4,4	3,1-5,7
Negative Control (INDEX)	0,1	<0,7

CALIBRATOR


0	6	3	9	5	2	5	0	0	0	0	0	2	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *[Signature]*



CONTROL SERUM

0	6	3	9	5	2	5	3	1	5	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES



Operator Signature *[Signature]*



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. *[Signature]*

DATE 20.07.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.